BROADBAND WHITE LIGHT SOURCE

USER MANUAL

Manual version: 1.0.B

Product code: G2002

Product Version: 1.0

Contents

- 1. Overview 3
- 2. EU Declaration of Conformity 4
- 3. Safety..... 3
 - 3.1 Warning 3
 - 3.2 Use of Equipment..... 3
 - 3.3 Servicing 3
- 4. Unpacking..... 4
 - 4.1 Packing List..... 4
 - 4.2 Damage Inspection..... 5
- 5. Operation..... 5

1. Overview

Ossila's USB powered broadband white light source is a simple and flexible light source for steady-state spectroscopy applications. It can be used in conjunction with a spectrometer to perform transmission and absorption measurements across the visible spectrum.

Key features

- Broadband emission (360 nm – 900 nm approx.)
- LED source for stable output
- USB powered
- Compact design

2. EU Declaration of Conformity

We

Company Name: Ossila BV

Postal Address: Biopartner 3 building, Galileiweg 8

Postcode: 2333 BD Leiden

Country: The Netherlands

Telephone number: +31 (0)71 3322992

Email Address: info@ossila.com

declare that the DoC is issued under our sole responsibility and belongs to the following product:

Product: Broadband White Light Source (G2002A1)

Serial number: G2002A1- xxxx

Object of declaration:

Broadband White Light Source (G2002A1)

The object of declaration described above is in conformity with the relevant Union harmonisation legislation:

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
3. Safety

3.1 Warning

The Ossila Broadband White Light Source is classified as being in Risk Group 1 (Low Risk) for blue light hazard under the IEC/EN 62471 (Photobiological Safety of Lamps and Lamp Systems) Standard. The RG0 (Exempt) exposure limit is exceeded at viewing distances below 50 cm, and at viewing times greater than 25 minutes at 20 cm. **Do not look directly into the source.**

3.2 Use of Equipment

The Ossila Broadband White Light is designed to be used as instructed. It is intended for use under the following conditions:

- Indoors in a laboratory environment (Pollution Degree 2).
- Altitudes up to 2000m.
- Temperatures of 5°C to 40°C; maximum relative humidity of 80% up to 31°C.

3.3 Servicing

If servicing is required, please return the unit to Ossila Ltd. The warranty will be invalidated if:

- Modification or service has taken place by anyone other than an Ossila engineer.
- The Unit has been subjected to chemical damage through improper use.
- The Unit has been operated outside the usage parameters stated in the user documentation associated with the Unit.
- The Unit has been rendered inoperable through accident, misuse, contamination, improper maintenance, modification, or other external causes.

4. Unpacking

4.1 Packing List

The standard items included with the Ossila Broadband White Light are:

- Ossila Broadband White Light.
- USB type-C cable.

Figure 4.1. Ossila Broadband White Light.


Figure 4.2. USB type-A to type-C cable.

4.2 Damage Inspection

Examine the components for evidence of shipping damage. If damage has occurred, please contact Ossila directly for further action.

5. Operation

The Ossila Broadband white light is powered by a USB type-C cable supplying 5V at 80mA. Any PC USB port or plug socket USB adapter can power the device. The source will turn on when the USB cable is plugged into the USB port on the rear of the device. A SMA905 optical fibre can be connected to the output port of the device.

**INNOVATION
AWARD 2017**

IOP Institute of Physics

THE QUEEN'S AWARDS
FOR ENTERPRISE:
INTERNATIONAL TRADE
2018