
# SILICON PHOTOMULTIPLIER USER MANUAL

Manual version: 1.0.0  
Product code: G3002  
Product Version: 1.0  
Software version: 1.0

# Contents

<b>1. EU Declaration of Conformity</b> .....	<b>3</b>
<b>2. Safety</b> .....	<b>6</b>
2.1 Use of Equipment.....	6
2.2 Servicing.....	6
2.3 Health and Safety – Servicing.....	6
<b>3. Requirements</b> .....	<b>7</b>
<b>4. Unpacking</b> .....	<b>7</b>
4.1 Packing List .....	7
4.2 Damage Inspection .....	8
<b>5. Specifications</b> .....	<b>8</b>
<b>6. Operation</b> .....	<b>9</b>
6.1 Command Library .....	9
6.2 Photodetector Console.....	11
<b>7. Troubleshooting</b> .....	<b>17</b>
<b>8. Related Products</b> .....	<b>18</b>

# 1. EU Declaration of Conformity

## We

**Company Name:** Ossila BV

**Postal Address:** Biopartner 3 building, Galileiweg 8

**Postcode:** 2333 BD Leiden

**Country:** The Netherlands

**Telephone number:** +31 (0)71 3322992

**Email Address:** info@ossila.com

**declare that the DoC is issued under our sole responsibility and belongs to the following product:**

**Product:** USB Silicon Photomultiplier (G3002A1)

**Serial number:** G3002A1-xxxx

## Object of declaration:

USB Silicon Photomultiplier (G3002A1)

**The object of declaration described above is in conformity with the relevant Union harmonisation legislation:**

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## 2. Safety

### 2.1 Use of Equipment

The Ossila USB Silicon Photomultiplier is designed to be used as instructed. It is intended for use under the following conditions:

- Indoors in a laboratory environment (Pollution Degree 2)
- Altitudes up to 2000m
- Temperatures of 5°C to 40°C; maximum relative humidity of 80% up to 31°C.

### 2.2 Servicing

If servicing is required, please return the unit to Ossila Ltd. The warranty will be invalidated if:

- Modification or service has been carried out by anyone other than an Ossila engineer.
- The Unit has been subjected to chemical damage through improper use.
- The Unit has been operated outside the usage parameters stated in the user documentation associated with the Unit.
- The Unit has been rendered inoperable through accident, misuse, contamination, improper maintenance, modification, or other external causes.

### 2.3 Health and Safety – Servicing

Servicing should only be performed by an Ossila engineer. Any modification or alteration may damage the equipment, cause injury, or death. It will also void your equipment's warranty.

## 3. Requirements

Table 3.1 details the requirements for the Ossila USB Silicon Photomultiplier, and the minimum computer specifications for the Ossila Photodetector Console software.

Table 3.1. Ossila USB Silicon Photomultiplier and Ossila Photodetector Console requirements.

Power	USB
Operating Systems	Windows 11 (64-bit)
CPU	Dual Core 2 GHz
RAM	4 GB
Available Hard Drive Space	183 MB
Monitor Resolution	1280 x 720
Connectivity	USB-C 2.0

## 4. Unpacking

### 4.1 Packing List

The standard items included with the Ossila USB Silicon Photomultiplier are:

- The Ossila USB Silicon Photomultiplier.
- USB-C cable.
- USB memory stick loaded with software installer, and user manual.

Figure 4.1. Ossila USB Silicon Photomultiplier.

## 4.2 Damage Inspection

Examine the components for evidence of shipping damage. If damage has occurred, please contact Ossila directly for further action. The shipping packaging will come with a shock indicator to show if there has been any mishandling of the package during transportation.

## 5. Specifications

The Ossila USB Silicon Photomultiplier specifications are shown in **Table 5.1**.

**Table 5.1.** Ossila USB Silicon Photomultiplier specifications.

<b>Spectral Range</b>	200 nm – 1000 nm
<b>Irradiance Measurement Range</b>	1 nW/cm <sup>2</sup> – 10 μW/cm <sup>2</sup>
<b>Photodetection Efficiency</b>	Up to 35%
<b>Maximum Counts</b>	16777215
<b>Signal-to-Noise Ratio</b>	>50,000:1
<b>Maximum Data Rate</b>	30 Hz
<b>Detection Area</b>	3 mm x 3 mm
<b>Dimensions</b>	65 mm (L) x 65 mm (W) x 14 mm (D)
<b>Weight</b>	100 g

## 6. Operation

The Ossila USB Silicon Photomultiplier can be controlled either programmatically using serial commands, or through the Ossila Photodetector Console software.

### 6.1 Command Library

The primary interface of the Ossila USB Silicon Photomultiplier is the serial command library. When connected to a PC, the device will appear as a COM port, to which the serial commands can be sent. This section describes the command protocol and lists the available commands.

#### 6.1.1 Command Format

Commands should be sent to the silicon photomultiplier in ASCII format, and responses from the silicon photomultiplier will also be in ASCII format.

All commands sent to the silicon photomultiplier have a start and end delimiter, `<` and `>` respectively. Only commands enclosed by these delimiters will be acknowledged by the device. If any invalid commands are sent between delimiters, the system will return `<Invalid Command>`.

Commands are either action commands with the format `<command>` or setting commands with the format `<command value>`. Setting commands can also be queried using the format `<command?>`.

Setting commands that have a corresponding query command list both commands together. In these cases, the set command value and query command return value have the same format.

Values will be one of the following types depending on the command:

- `float` – Floating-point number
  - l. These can be sent in scientific notation (for example, `123e4` or `123E4`).
- `int` – Integer
  - l. If a floating-point number is sent instead, it will be rounded down to an integer.
- `bool` – Boolean (0 or 1)
  - l. Non-zero integers or floating-point numbers will be interpreted as 1.
- `str` – String
  - l. Only used as returns in some query commands.

## 6.1.2 Action Commands

Command	Returns	Function
<read>	<read <i>int</i> >	Reads the photodetector and returns the current value.
<zero>	<zero <i>int</i> >	Automatically sets the baseline to make the current reading approximately 0.02% of the full scale.

## 6.1.3 Setting Commands

Command	Values>Returns	Function
<tcomp <i>X</i> > <tcomp?>	<i>X</i> - bool - temperature compensation state	1 sets the SiPM to apply temperature compensation to measurements, 0 turns temperature compensation off.
<baseline <i>X</i> > <baseline?>	<i>X</i> - int - baseline value	Sets an offset to measurements. Can be set between -10000 and 10000. These values are arbitrary, <b>not</b> measured counts.
<precision <i>X</i> > <precision?>	<i>X</i> - int - precision level (1, 2, 3)	Sets the precision of the silicon multiplier to: 1 - low 2 - medium 3 - high Higher precision decreases measurement noise and speed.
<temp?>	<temp <i>float</i> >	Returns the temperature of the silicon photomultiplier in °C.
<serial?>	<serial <i>int</i> >	Returns the serial number of the silicon photomultiplier.
<firmware?>	<firmware <i>int.int.int</i> >	Returns the firmware version of the silicon photomultiplier in the format <i>major.minor.patch</i> .

## 6.2 Photodetector Console


Figure 6.1. Ossila Photodetector Console.

Ossila provides free software to allow for simple use of the silicon photomultiplier. A Windows installer is provided on the USB drive supplied with the system and can be downloaded from our website at [www.ossila.com/pages/software-drivers](http://www.ossila.com/pages/software-drivers).

To install the software simply run the installer and follow the on-screen instructions.

### 6.2.1 Control Tool Bar



Figure 6.2. Photodetector Console control tool bar.

The control tool bar has the following sections:

#### (I) Connection Controls

The drop-down box lists the serial numbers and addresses of any connected Ossila USB Silicon Photomultipliers.

To connect to the selected silicon photomultiplier, click the  button. To disconnect click the  button.

To search for connected Ossila USB Silicon Photomultipliers, click the  button. If only a single silicon photomultiplier is detected the software will automatically connect to it.

## (II) Photodetector Controls


Figure 6.3. Silicon photomultiplier controls.

<b>Precision</b>	Set the precision of the silicon photomultiplier.
<b>Baseline Adjustment</b>	Set the baseline adjustment of the silicon photomultiplier.
	Automatically sets the baseline to make the current reading approximately 0.02% of the full scale.
	Toggle temperature compensation on or off.
	Start or stop reading data from the silicon photomultiplier.

## (III) Plot Controls


Figure 6.4. Photodetector Console plot controls.

<b>Data Points</b>	Set the number of data points to display in the plot.
	Zeroes all data to the most recent data point. Click again to restore raw values.
	Toggle displaying a grid on the plot.
	Save the data in the plot to a .csv (comma-separated value) or .txt (text) file.
	Clear all data from the plot.

## Plot Display Controls

By default, the plots will automatically scale the axes to display all the data within it. The view can be controlled manually using the following mouse controls:

- Left or middle click and drag – pan the axes.
- Right click and drag – scale the axes.
- Scroll wheel – scale the axes.

A specific axis can be controlled by using these controls on the axis labels. The axes can be reset by clicking the 'A' button in the bottom-left of the plot, shown in **Figure 6.5** (note, this button will only appear whilst the mouse cursor is over the plots).


Figure 6.5. Button to reset plot axes.

Axis scales and automatic scaling can also be set by right-clicking on the plot and selecting the appropriate options in the context menu.

## (IV) Data Display


Figure 6.6. Photodetector Console data display.

<b>Live Measurement</b>	The most recent data point to be read from the silicon photomultiplier.
<b>Maximum</b>	The maximum value displayed in the plot.
<b>Minimum</b>	The minimum value displayed in the plot.
<b>Cursor</b>	The value of the data point closest to the cursor's position.

All values are in counts and are relative to the zero value if the data has been zeroed.

## 6.2.2 Menu Bar

### (I) File Menu


Figure 6.7. File menu.

**Search for Ossila Photodetectors**

Search the USB ports and network for Ossila USB Silicon Photomultipliers. Discovered units will be displayed in the drop-down box in tool bar.

**Exit**

Close the software.

### (II) View Menu


Figure 6.8. View menu.

**Toggle Theme**

Switch the software between light and dark themes.

**Edit Plot Colours**

Open the edit plot colours settings dialog (see Figure 6.9).


Figure 6.9. Edit plot colours dialog.

### (III) Help Menu


Figure 6.10. Help menu.


<b>User Manual</b>	Open the user manual in your default web browser.
<b>Check for updates</b>	Check whether a newer version of the software is available to download and install. If an update is available a message will be displayed in the status bar with the latest version number. Clicking the message will take you to <a href="#">Ossila's website</a> where you can download the new version.
<b>Check for Updates on Startup</b>	When checked the software will automatically check for software updates each time it starts.
<b>Enable Debug Logging</b>	When checked, the software will write more information to the log file. The log file is written to the temporary directory which can be accessed by typing <code>%temp%</code> into the file explorer's address bar.
<b>View License</b>	Display software license information.
<b>View Third Party License</b>	Display the license information of third-party software that has been used in the software.
<b>About</b>	Display information about the software.

## 7. Troubleshooting

Most of the issues that may arise will be detailed here. However, if you encounter any issues that are not detailed here, then contact us by email at [info@ossila.com](mailto:info@ossila.com). We will respond as soon as possible.

Problem	Possible cause	Action
Cannot connect to the system via USB	The USB cable may not be connected properly.	Ensure the USB cable is firmly plugged in at both ends.
	The USB cable may not be connected to a working USB port.	Try connecting the system to a different USB port on the computer. If using a USB hub, try plugging the unit into a port directly on the PC.
	The USB cable is defective.	Try using a different USB cable
Software does not start	The wrong version of Windows is installed on the computer.	Install the software on a computer with Windows 11.
	The software has not installed properly.	Try reinstalling the software.


## 8. Related Products

### Monochromator

Isolate specific wavelengths between 380 – 1000 nm with <1 nm accuracy.

Product codes: G3001A1

### Tunable LED Light Source

Source variable monochromatic light between 380 and 1000 nm to <1 nm accuracy.

Product code: G3006A1

### LED Light Source

Powerful and consistent 380 – 1000 nm light for your spectroscopy experiments.

Product codes: G3003A1

### Adjustable Slit

A manual variable slit with a width range of 50  $\mu\text{m}$  to 6 mm.


Product code: G3004A1

### Optical Filters

Longpass and shortpass filters in a range of wavelengths.

Product code: G2101 / G2102

### Optical Breadboard Plate

Solid anodized metric aluminium optical breadboard plate.

Product codes: G2007A1