

Ossila

enabling science

MATERIALS CATALOG 2025 Q4




Contents

Chemistry Building Blocks	4-33
• Fluorinated Building Blocks	6-16
• Heterocyclic Building Blocks	17-28
• Non-Heterocyclic Building Blocks	29-33
Monomers	34-47
• Boronate Monomers	36-37
• Organotic Compounds	38
• Dibromo Monomers.....	39-42
• Carbaldehyde Monomers	43-45
• Diamines and Dianhydrides	46-47
Porous Organic Frameworks	48-60
• MOF Ligands	50-51
• COF Ligands	52-59
Low Dimensional Materials	60-69
• Carbon Nanotubes (CNTs)	62
• 2D Materials	63-69
Battery Materials	70-77
• Cathode Active Materials	72
• Anode Active Materials	73-75
• Battery Binders.....	76
Optoelectronic Materials	78-81
• Photonic and Optical Materials	80
• Liquid Crystals	81
Fullerenes and NFAs	82-87
• Fullerene Acceptors	84
• Non-Fullerene Acceptors	85-87
Semiconducting Polymers	88-103
• OPV Polymers	90-92
• OFET and OLED Polymers	93-94
• Interface Polymers	95-96
• Luminosyn™ Polymers	98-103
Perovskite Materials	104-117
• Perovskite Inks.....	106
• Perovskite Precursors	107-108
• Perovskite Interface Materials	109-113
• Self-Assembled Monolayers	114-117
DSSC Dyes and Electrolytes	118-119
OLED Materials	120-159
• Charge Transport Layer Materials	122-127
• Host Materials	128-132
• Dopant Materials	133-138
• TADF Materials	139-147
• High Purity Sublimed Materials	148-159
Organic Conductors	160-185
• Cathode Interlayer Materials	162
• Green Energy Materials	163-165
• Semiconducting Molecules	166-181
• PEDOT:PSS & PEDOT Polymer Blends	182-185
Inorganic Electronic Materials	186-193
• Metal Oxides	188-191
• Advanced Ceramic Materials	190-193


Chemistry Building Blocks

Building blocks for the chemical synthesis of materials with tailored structures and properties. These are the foundational units, or 'component parts', required for creating a wide array of functional materials with extensive applications.



- Fluorinated building blocks
- Heterocyclic building blocks
- Non-heterocyclic building blocks


Popular Products



BDD

High-purity, bromide-capped PCE12 monomer

- Electron deficient
- Benzodithiophene-dione building block
- Fused Thiophene Core



BDTT

High-purity PCE10 and PCE12 monomer

- Capped with trimethyltin
- Electron-rich pendants
- Good solubility



3-Hexylthiophene

High-purity monomer for popular polymer donor P3HT

- A synthesis precursor
- Hexylthiophene building block
- High purity

To help you choose the right material, full technical specifications are available on our website.

Ossila

B371-500mg

BDD

1,3-Bis(5-bromothiophen-2-yl)-5,7-bis(2-ethylhexyl)benzo[1,2-c:4,5-c']dithiophene-4,8-dione

Purity (>98%)
MW 766.73g/mol

CAS# 1415929-78-0

















Ossila

B171-1g

2,5-bis(trimethylsilyl)-thieno[3,2-b]thiophene

















Fluorinated Building Blocks

	2-Bromo-6-fluoronaphthalene SKU: B1331 CAS Number: 324-41-4 Options: • Quantity: 1 g, 5 g		2-(Difluoromethyl)benzimidazole SKU: B1571 CAS Number: 705-09-9 Options: • Quantity: 5 g, 10 g
	5,5'-dibromo-3,3'-difluoro-2,2'-bithiophene SKU: B1381 CAS Number: 1619967-08-6 Options: • Quantity: 500 mg, 1 g, 5 g		6-Fluoro-4-chromanone SKU: B1581 CAS Number: 66892-34-0 Options: • Quantity: 10 g, 25 g
	5,5'-Bis(trifluoromethyl)-2,2'-bipyridine SKU: B1451 CAS Number: 142946-80-3 Options: • Quantity: 1 g, 2 g, 5 g		2-Fluoro-9H-carbazole SKU: B1591 CAS Number: 391-53-7 Options: • Quantity: 1 g, 2 g
	4-(Trifluoromethyl)-2-pyridone SKU: B1481 CAS Number: 50650-59-4 Options: • Quantity: 1 g, 5 g, 10 g		7-Hydroxy-4-(trifluoromethyl)coumarin SKU: B1601 CAS Number: 575-03-1 Options: • Quantity: 1 g, 5 g
	4-Fluoro-1H-indazole SKU: B1501 CAS Number: 341-23-1 Options: • Quantity: 1 g, 5 g		4,7-Dibromo-5,6-difluoro-2,1,3-benzothiadiazole SKU: B161 CAS Number: 1295502-53-2 Options: • Quantity: 1 g, 5 g
	PCE11 Monomer SKU: B151 CAS Number: 1504626-07-6 Options: • Quantity: 1 g, 2 g		7-(Diethylamino)-4-(trifluoromethyl)coumarin SKU: B1611 CAS Number: 41934-47-8 Options: • Quantity: 250 mg, 500 mg
	3-(4-Methyl-1H-imidazol-1-yl)-5-(trifluoromethyl)aniline SKU: B1511 CAS Number: 641571-11-1 Options: • Quantity: 5 g, 10 g		2-Furoyltrifluoroacetone SKU: B1621 CAS Number: 326-90-9 Options: • Quantity: 25 g, 50 g
	3,3-Difluoroazetidine hydrochloride SKU: B1561 CAS Number: 288315-03-7 Options: • Quantity: 5 g		2-(Trifluoromethyl)thioxanthen-9-one SKU: B1631 CAS Number: 1693-28-3 Options: • Quantity: 10 g

	5-(Trifluoromethyl)-2-pyridone SKU: B1641 CAS Number: 33252-63-0 Options: • Quantity: 25 g, 50 g		6-Fluoroindole SKU: B1721 CAS Number: 399-51-9 Options: • Quantity: 10 g, 25 g, 50 g
	3-(Trifluoromethyl)-2-pyridone SKU: B1651 CAS Number: 22245-83-6 Options: • Quantity: 10 g, 25 g, 50 g		6-Fluoroindole-2-carboxylic acid SKU: B1731 CAS Number: 3093-97-8 Options: • Quantity: 1 g, 2 g
	1-[2-(Trifluoromethyl)phenyl]imidazole SKU: B1661 CAS Number: 25371-96-4 Options: • Quantity: 1 g, 2 g		5-(Trifluoromethoxy)isatin SKU: B1741 CAS Number: 169037-23-4 Options: • Quantity: 5 g, 10 g
	4-(Trifluoromethyl)-1H-imidazole SKU: B1671 CAS Number: 3468-69-8 Options: • Quantity: 5 g, 10 g		5-Fluoroisatin SKU: B1751 CAS Number: 56341-41-4 Options: • Quantity: 10 g, 25 g
	4-Fluoroindole SKU: B1681 CAS Number: 387-43-9 Options: • Quantity: 5 g, 10 g		5-Fluoroisatin SKU: B1761 CAS Number: 443-69-6 Options: • Quantity: 25 g, 50 g
	3-(4-Fluorophenyl)-1-isopropylindole SKU: B1691 CAS Number: 93957-49-4 Options: • Quantity: 25 g, 50 g		Benzyl (3-Fluoro-4-morpholinophenyl)carbamate SKU: B1771 CAS Number: 168828-81-7 Options: • Quantity: 25 g, 50 g
	7-Fluoroindole SKU: B1701 CAS Number: 387-44-0 Options: • Quantity: 5 g, 10 g		3-Fluoro-4-morpholinoaniline SKU: B1781 CAS Number: 93246-53-8 Options: • Quantity: 25 g, 50 g
	5-Fluoroindole SKU: B1711 CAS Number: 399-52-0 Options: • Quantity: 5 g, 10 g		7-Chloro-1-cyclopropyl-6-fluoro-1,4-dihydro-4-oxo-1,8-naphthyridine-3-carboxylic acid SKU: B1791 CAS Number: 100361-18-0 Options: • Quantity: 25 g, 50 g

















	5-Bromo-3-fluorosalicylaldehyde SKU: B1831 CAS Number: 251300-28-4 Options: • Quantity: 2 g, 5 g, 10 g		2-(2,4-Difluorophenyl)-5-methylpyridine SKU: B1911 CAS Number: 583052-21-5 Options: • Quantity: 1 g, 2 g, 5 g
	Ethyl (4-fluorobenzoyl)acetate SKU: B1841 CAS Number: 1999-00-4 Options: • Quantity: 5 g, 10 g, 25 g		2-(2,4-Difluorophenyl)-5-fluoropyridine SKU: B1921 CAS Number: 1426047-01-9 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g
	2,4,6-Trifluoronitrobenzene SKU: B1851 CAS Number: 315-14-0 Options: • Quantity: 5 g, 10 g, 25 g		2,2-Bis(3-amino-4-hydroxyphenyl)hexafluoropropane (6FAP) SKU: B1971 CAS Number: 83558-87-6 Options: • Quantity: 25 g, 50 g, 100 g, 200 g
	(5S)-(-)-5-(Trifluoromethyl)-2-pyrrolidin-one SKU: B1861 CAS Number: 1287211-10-2 Options: • Quantity: 1 g, 2 g, 5 g		2,2'-Bis(trifluoromethyl)benzidine (TFMB) SKU: B2001 CAS Number: 341-58-2 Options: • Quantity: 25 g, 50 g, 100 g
	8-Fluoro-4-hydroxy-2-methylquinoline SKU: B1871 CAS Number: 5288-22-2 Options: • Quantity: 1 g, 2 g, 5 g		4,4'-(Hexafluoroisopropylidene)diphthalic anhydride (6FDA) SKU: B2061 CAS Number: 1107-00-2 Options: • Quantity: 10 g, 25 g, 50 g
	6-Chloro-8-fluoroquinoline SKU: B1881 CAS Number: 52200-53-0 Options: • Quantity: 1 g, 5 g, 10 g		1,4-Bis(4-amino-2-trifluoromethylphenoxy)benzene (6FAPB) SKU: B2071 CAS Number: 94525-05-0 Options: • Quantity: 10 g, 25 g, 50 g
	4-Bromo-8-fluoroquinoline SKU: B1891 CAS Number: 927800-38-2 Options: • Quantity: 500 mg, 1 g, 2 g		2,2-Bis[4-(4-aminophenoxy)phenyl]hexafluoroethane (4-BDAF) SKU: B2081 CAS Number: 69563-88-8 Options: • Quantity: 25 g, 50 g, 100 g
	4-Chloro-7-trifluoromethylquinoline-3-carboxylic acid ethyl ester SKU: B1901 CAS Number: 1168-42-3 Options: • Quantity: 500 mg, 1 g, 2 g		2,2'-Bis(trifluoromethyl)-4,4'-diaminodiphenyl ether (6FODA) SKU: B2091 CAS Number: 344-48-9 Options: • Quantity: 10 g, 25 g, 50 g, 100 g

















	2,2-Bis(3-(3-aminobenzoylamino)-4-hydroxyphenyl)hexafluoropropane (m-6FD-AP) SKU: B2161 CAS Number: 220426-92-6 Options: • Quantity: 5 g, 10 g, 25 g		2,4-Difluorobenzonitrile SKU: B2251 CAS Number: 3939-09-1 Options: • Quantity: 25 g, 50 g
	p-6FDAP SKU: B2171 CAS Number: 129197-38-2 Options: • Quantity: 5 g, 10 g, 25 g		2-Amino-4-fluorophenol SKU: B2261 CAS Number: 399-97-3 Options: • Quantity: 5 g, 10 g
	AB-TFMB SKU: B2181 CAS Number: 1449757-11-2 Options: • Quantity: 2 g, 5 g, 10 g		4-Fluoro-3-nitrobenzoic acid SKU: B2271 CAS Number: 453-71-4 Options: • Quantity: 25 g, 50 g
	6-Fluoro-3-(4-piperidinyl)-1,2-benzoxazole hydrochloride SKU: B2201 CAS Number: 84163-13-3 Options: • Quantity: 10 g, 25 g, 50 g		3,3-Difluoropyrrolidine hydrochloride SKU: B2281 CAS Number: 163457-23-6 Options: • Quantity: 2 g, 5 g, 10 g
	1-(tert-Butoxycarbonyl)-3-fluoro-4-piperidone SKU: B2211 CAS Number: 211108-50-8 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g		2-Fluoroadenine SKU: B2291 CAS Number: 700-49-2 Options: • Quantity: 1 g, 2 g, 5 g
	5-Bromo-2-fluoropyridine SKU: B2221 CAS Number: 766-11-0 Options: • Quantity: 25 g, 50 g		4-Bromo-2,6-difluoroaniline SKU: B2301 CAS Number: 67567-26-4 Options: • Quantity: 5 g, 10 g, 25 g
	4,4'-Bis(trifluoromethyl)-2,2'-bipyridine SKU: B2231 CAS Number: 142946-79-0 Options: • Quantity: 500 mg, 1 g, 2 g		(R)-(-)-3-Fluoropyrrolidine hydrochloride SKU: B2311 CAS Number: 136725-55-8 Options: • Quantity: 1 g, 2 g, 5 g
	2,4-Difluoro-(1H-1,2,4-triazolyl)acetophenone SKU: B2241 CAS Number: 86404-63-9 Options: • Quantity: 10 g, 25 g, 50 g		5-Fluorosalicylic acid SKU: B2321 CAS Number: 345-16-4 Options: • Quantity: 10 g, 25 g, 50 g

















	Methyl pentafluorobenzoate SKU: B2331 CAS Number: 36629-42-2 Options: • Quantity: 10 g, 25 g, 50 g		2-Bromo-4-fluoropyridine SKU: B2411 CAS Number: 357927-50-5 Options: • Quantity: 5 g, 10 g
	(S)-(+)-3-Fluoropyrrolidine hydrochloride SKU: B2341 CAS Number: 136725-53-6 Options: • Quantity: 1 g, 2 g, 5 g		6-Fluoronicotinic acid SKU: B2421 CAS Number: 403-45-2 Options: • Quantity: 1 g, 2 g, 5 g
	N-Boc-trans-4-fluoro-L-proline SKU: B2351 CAS Number: 203866-14-2 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g		3,5-Difluorobenzonitrile SKU: B2431 CAS Number: 64248-63-1 Options: • Quantity: 10 g, 25 g, 50 g
	Bis(4-fluorophenyl)methanol SKU: B2361 CAS Number: 365-24-2 Options: • Quantity: 50 g, 100 g, 250 g		1-Bis(4-fluorophenyl)methyl piperazine SKU: B2441 CAS Number: 27469-60-9 Options: • Quantity: 10 g, 25 g
	2-Fluoro-5-nitrobenzonitrile SKU: B2371 CAS Number: 17417-09-3 Options: • Quantity: 10 g, 25 g, 50 g		2-Fluoro-5-nitrobenzoic acid SKU: B2451 CAS Number: 7304-32-7 Options: • Quantity: 10 g, 25 g, 50 g
	2-Bromo-6-fluorobenzaldehyde SKU: B2381 CAS Number: 360575-28-6 Options: • Quantity: 5 g, 10 g, 25 g		4-Fluorobenzenesulfonamide SKU: B2461 CAS Number: 402-46-0 Options: • Quantity: 10 g, 25 g, 50 g
	2-Chloro-5-fluoropyrimidine SKU: B2391 CAS Number: 62802-42-0 Options: • Quantity: 5 g, 10 g		3,4,5-Trifluoroaniline SKU: B2481 CAS Number: 163733-96-8 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	3-Fluoro-2-hydroxybenzaldehyde SKU: B2401 CAS Number: 394-50-3 Options: • Quantity: 5 g, 10 g		4-Fluororesorcinol SKU: B2491 CAS Number: 103068-41-3 Options: • Quantity: 1 g, 2 g






	2-Bromo-3'-fluoroacetophenone SKU: B2501 CAS Number: 53631-18-8 Options: • Quantity: 5 g, 10 g		α,α-Difluorophenylacetic acid SKU: B2581 CAS Number: 360-03-2 Options: • Quantity: 2 g, 5 g, 10 g
	3,4,5-Trifluorobenzoic acid SKU: B2511 CAS Number: 121602-93-5 Options: • Quantity: 5 g, 10 g, 25 g		3-Fluoro-4-hydroxybenzaldehyde SKU: B2591 CAS Number: 405-05-0 Options: • Quantity: 5 g, 10 g
	N-Boc-4,4-difluoro-L-proline SKU: B2521 CAS Number: 203866-15-3 Options: • Quantity: 500 mg, 1 g, 2 g		2-Amino-3-fluoropyridine SKU: B2601 CAS Number: 21717-95-3 Options: • Quantity: 2 g, 5 g, 10 g
	2-Fluoro-4-(trifluoromethyl)aniline SKU: B2511 CAS Number: 69409-98-9 Options: • Quantity: 5 g, 10 g		3-Fluoro-2-pyridinecarbonitrile SKU: B2611 CAS Number: 97509-75-6 Options: • Quantity: 5 g, 10 g, 25 g
	2,6-Dichloro-5-fluoro-3-pyridinecarbonitrile SKU: B2541 CAS Number: 82671-02-1 Options: • Quantity: 25 g, 50 g		2-Bromo-4-fluorobenzonitrile SKU: B2621 CAS Number: 36282-26-5 Options: • Quantity: 10 g, 25 g, 50 g
	4-Fluoro-β-nitrostyrene SKU: B2531 CAS Number: 706-08-1 Options: • Quantity: 5 g, 10 g		3'-Fluoro-4'-methoxyacetophenone SKU: B2631 CAS Number: 455-91-4 Options: • Quantity: 25 g
	2,2-Difluorocyclopropanecarboxylic acid SKU: B2561 CAS Number: 107873-03-0 Options: • Quantity: 1 g, 2 g, 5 g		3-Fluoro-p-anisidine SKU: B2641 CAS Number: 366-99-4 Options: • Quantity: 10 g, 25 g, 50 g
	3,3-Difluoropiperidine hydrochloride SKU: B2571 CAS Number: 496807-97-7 Options: • Quantity: 1 g, 2 g, 5 g		2-Amino-6-fluorobenzonitrile SKU: B2651 CAS Number: 77326-36-4 Options: • Quantity: 25 g, 50 g

	3-Fluoro-2-methylbenzoic acid SKU: B2661 CAS Number: 699-90-1 Options: • Quantity: 5 g, 10 g		3-Chloro-2-fluorobenzoic acid SKU: B2741 CAS Number: 161957-55-7 Options: • Quantity: 10 g, 25 g, 50 g
	2-Fluoro-4-nitroaniline SKU: B2671 CAS Number: 369-35-7 Options: • Quantity: 10 g, 25 g, 50 g		2-Fluoro-4-hydroxybenzoic acid SKU: B2751 CAS Number: 65145-13-3 Options: • Quantity: 10 g, 25 g
	4-Chloro-3-fluoroaniline SKU: B2681 CAS Number: 367-22-6 Options: • Quantity: 5 g, 10 g		2,6-Difluorocinnamic aldehyde SKU: B2761 CAS Number: 117338-43-9 Options: • Quantity: 250 mg, 500 mg
	2,4,5-Trifluorobenzonitrile SKU: B2721 CAS Number: 98349-22-5 Options: • Quantity: 10 g, 25 g, 50 g		5-Cyano-2-fluoropyridine SKU: B2771 CAS Number: 3939-12-6 Options: • Quantity: 1 g, 2 g, 5 g
	4-Fluoro-2-methylphenol SKU: B2701 CAS Number: 452-72-2 Options: • Quantity: 10 g, 25 g, 50 g		3,5-Difluoropyridine-2-carbonitrile SKU: B2781 CAS Number: 298709-29-2 Options: • Quantity: 5 g, 10 g
	2-Amino-5-fluorobenzonitrile SKU: B2711 CAS Number: 61272-77-3 Options: • Quantity: 5 g, 10 g		3-Fluoro-4-(trifluoromethyl)benzoic acid SKU: B2791 CAS Number: 115754-21-7 Options: • Quantity: 10 g, 25 g, 50 g
	3-Fluorobenzoic acid hydrazide SKU: B2721 CAS Number: 499-55-8 Options: • Quantity: 10 g, 25 g		4-Fluoro-2-methylbenzonitrile SKU: B2801 CAS Number: 147754-12-9 Options: • Quantity: 10 g, 25 g, 50 g
	2-Fluoro-3-methoxybenzaldehyde SKU: B2731 CAS Number: 103438-88-6 Options: • Quantity: 5 g, 10 g		4-Methoxy-3-(trifluoromethyl)aniline SKU: B2811 CAS Number: 393-15-7 Options: • Quantity: 1 g, 2 g, 5 g

















	3-Fluoro-4-methoxybenzoic acid SKU: B2821 CAS Number: 403-20-3 Options: • Quantity: 5 g, 10 g, 25 g		6-Fluoropyridine-2-carbonitrile SKU: B2901 CAS Number: 3939-15-9 Options: • Quantity: 1 g, 2 g, 5 g
	2-Fluoro-4-(trifluoromethyl)benzoic acid SKU: B2831 CAS Number: 115029-24-8 Options: • Quantity: 5 g, 10 g		3-Chloro-5-fluoroaniline SKU: B2911 CAS Number: 4863-91-6 Options: • Quantity: 2 g, 5 g, 10 g
	1-Bromo-2-chloro-4-fluorobenzene SKU: B2841 CAS Number: 110407-59-5 Options: • Quantity: 10 g, 25 g, 50 g		5-Bromo-2-fluoropyrimidine SKU: B2921 CAS Number: 62802-38-4 Options: • Quantity: 1 g, 2 g, 5 g
	2-Bromo-5-fluorobenzonitrile SKU: B2851 CAS Number: 57381-39-2 Options: • Quantity: 10 g, 25 g, 50 g		4-Hydroxy-2-(trifluoromethyl)benzoic acid SKU: B2931 CAS Number: 320-32-1 Options: • Quantity: 1 g, 2 g
	2-Fluoropyridine-3-carbonitrile SKU: B2861 CAS Number: 3939-13-7 Options: • Quantity: 10 g, 25 g, 50 g		Methyl 6-fluoropyridine-2-carboxylate SKU: B2941 CAS Number: 455-71-0 Options: • Quantity: 1 g, 2 g
	2,3,5,6-Tetrafluoro-4-hydroxybenzoic acid SKU: B2871 CAS Number: 652-34-6 Options: • Quantity: 1 g, 2 g, 5 g		1-(3-Fluorophenyl)-1H-pyrazole-4-carbaldehyde SKU: B2951 CAS Number: 936940-82-8 Options: • Quantity: 1 g, 2 g
	2,4,6-Trifluorophenylacetic acid SKU: B2881 CAS Number: 209991-63-9 Options: • Quantity: 250 mg, 500 mg, 1 g		4-Amino-2-(trifluoromethyl)pyridine SKU: B2961 CAS Number: 147149-98-2 Options: • Quantity: 1 g, 2 g, 5 g
	4-Fluoro-4'-methoxybiphenylene SKU: B2891 CAS Number: 345-89-1 Options: • Quantity: 5 g, 10 g, 25 g		4-Bromo-6-fluoroquinoline SKU: B2971 CAS Number: 661463-17-8 Options: • Quantity: 500 mg, 1 g, 2 g

















	6-Fluoro-4-hydroxy-2-(trifluoromethyl)quinoline SKU: B2981 CAS Number: 31009-34-4 Options: • Quantity: 500 mg, 1 g, 2 g		trans-3,4-Difluorocinnamic acid SKU: B3281 CAS Number: 112897-97-9 Options: • Quantity: 10 g, 25 g, 50 g
	5-Difluoromethoxy-2-mercaptobenzimidazole SKU: B2991 CAS Number: 97963-62-7 Options: • Quantity: 25 g, 50 g		2-Fluoro-5-methylbenzoic acid SKU: B3291 CAS Number: 321-12-0 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	4-Bromo-3-fluoriodobenzene SKU: B3221 CAS Number: 136434-77-0 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		5-Chloro-2-fluorobenzoic acid SKU: B3301 CAS Number: 394-30-9 Options: • Quantity: 2 g, 5 g, 10 g, 25 g
	3-Bromo-5-fluorophenol SKU: B3231 CAS Number: 433939-27-6 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		4-Fluoro-3-methylbenzoic acid SKU: B3311 CAS Number: 403-15-6 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	4-Bromo-3-fluorobenzonitrile SKU: B3241 CAS Number: 133059-44-6 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		4'-Fluoro-3'-nitroacetophenone SKU: B3321 CAS Number: 400-93-1 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	2-Chloro-5-fluoronitrobenzene SKU: B3251 CAS Number: 345-17-5 Options: • Quantity: 2 g, 5 g, 10 g, 25 g		2-Fluoro-6-methylbenzoic acid SKU: B3331 CAS Number: 90259-27-1 Options: • Quantity: 2 g, 5 g, 10 g
	3-Fluoro-4-iodopyridine SKU: B3261 CAS Number: 22282-75-3 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g		5-Chloro-2-fluorobenzonitrile SKU: B3341 CAS Number: 57381-34-7 Options: • Quantity: 2 g, 5 g, 10 g, 25 g
	4-Fluoro-2-iodobenzoic acid SKU: B3271 CAS Number: 56096-89-0 Options: • Quantity: 1 g, 2 g, 5 g, 10 g		3-Fluoro-5-(trifluoromethyl)benzoic acid SKU: B3351 CAS Number: 161622-05-5 Options: • Quantity: 2 g, 5 g, 10 g, 25 g

















	3'-Bromo-4'-fluoroacetophenone SKU: B3361 CAS Number: 1007-15-4 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		2,4,5-Trifluorophenylacetic acid SKU: B3441 CAS Number: 209995-38-0 Options: • Quantity: 2 g, 5 g, 10 g, 25 g
	5-Bromo-4-fluoro-2-methylaniline SKU: B3371 CAS Number: 627871-16-3 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		Pentafluorobenzoic acid SKU: B3451 CAS Number: 602-94-8 Options: • Quantity: 25 g, 50 g, 100 g
	5-Fluoro-2-methylbenzoic acid SKU: B3381 CAS Number: 33184-16-6 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		2-Fluorobiphenyl SKU: B3461 CAS Number: 321-60-8 Options: • Quantity: 10 g, 25 g, 50 g, 100 g
	2-Amino-5-chloro-3-fluoropyridine SKU: B3391 CAS Number: 246847-98-3 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		Potassium nonafluoro-1-butanedisulfonate SKU: B3471 CAS Number: 29420-49-3 Options: • Quantity: 25 g, 50 g, 100 g
	7-Fluoro-4-hydroxycoumarin SKU: B3401 CAS Number: 2145-27-9 Options: • Quantity: 100 mg, 250 mg, 500 mg		2-Ethylhexyl 4,6-dibromo-3-fluorothiophene-2-carboxylate SKU: B361 CAS Number: 1237479-38-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	5-Bromo-3-fluoropyridine-2-carboxylic acid SKU: B3411 CAS Number: 669066-91-5 Options: • Quantity: 1 g, 2 g, 5 g		3,3'-Difluoro-5,5'-bis(trimethylstannyl)-2,2'-bithiophene SKU: B391 CAS Number: 161967-09-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	5-Fluorobenzoxazole-2-thiol SKU: B3421 CAS Number: 13451-78-0 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g		1,3,5-Tris(3-fluoro-4-formylphenyl)benzene SKU: B3921 CAS Number: 2088007-04-7 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4-Fluoro-1,2-phenylenediamine SKU: B3431 CAS Number: 367-31-7 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		DFBTA-2ThBr-HD SKU: B421 CAS Number: 1887135-96-7 Options: • Quantity: 500 mg, 1 g

















	BDTT-F-Sn SKU: B521 CAS Number: 1514905-25-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	2FIC SKU: B551 CAS Number: 2083617-82-5 Options: • Quantity: 500 mg, 1 g, 5 g
	MeHDOQx-DBrDF SKU: B771 CAS Number: 2414359-15-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	HDOQx-DBrDF SKU: B781 CAS Number: 2269476-12-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	FIC SKU: B801 CAS Number: 507484-54-0 Options: • Quantity: 500 mg, 1 g, 5 g

















Heterocyclic Building Blocks

	4,7-Dibromo-2,1,3-benzothiadiazole SKU: B101 CAS Number: 15155-41-6 Options: • Quantity: 5 g, 10 g		5,5'-Dibromo-3,3'-bis(2-(2-(2-methoxyethoxy)ethoxy)ethoxy)-2,2'-bithiophene SKU: B1121 CAS Number: 1974310-58-1 Options: • Quantity: 500 mg, 1 g, 5 g
	3,4-(2,2-Dimethylpropylenedioxy)thiophene SKU: B1041 CAS Number: 255901-50-9 Options: • Quantity: 1 g, 5 g, 10 g		3,6-Di-tert-butylcarbazole SKU: B1131 CAS Number: 37500-95-1 Options: • Quantity: 10 g, 25 g, 100 g, 250 g
	2,6-Dibromodithieno[3,2-b:2',3'-d]thiophene SKU: B1051 CAS Number: 67061-69-2 Options: • Quantity: 500 mg, 1 g, 2 g		3,6-Diphenylcarbazole SKU: B1141 CAS Number: 56525-79-2 Options: • Quantity: 10 g, 25 g, 100 g, 250 g
	2,2':5',2''-Terthiophene SKU: B1061 CAS Number: 1081-34-1 Options: • Quantity: 10 g, 25 g, 50 g		Indolo[3,2-b]carbazole SKU: B1151 CAS Number: 6336-32-9 Options: • Quantity: 5 g, 10 g, 25 g
	5,5''-Dibromo-2,2':5',2''-terthiophene SKU: B1071 CAS Number: 98057-08-0 Options: • Quantity: 1 g, 5 g		10,15-Dihydro-5H-diindolo[3,2-a:3',2'-c]-arbazole SKU: B1161 CAS Number: 109005-10-9 Options: • Quantity: 1 g, 5 g
	3,3',5,5'-Tetrabromo-2,2'-bithiophene SKU: B1081 CAS Number: 125143-53-5 Options: • Quantity: 10 g, 25 g, 50 g, 100 g		2-Bromo-9H-carbazole SKU: B1171 CAS Number: 3652-90-2 Options: • Quantity: 25 g, 50 g, 100 g, 250 g
	2,3,5,6-Tetrabromothieno[3,2-b]thiophene SKU: B1091 CAS Number: 124638-53-5 Options: • Quantity: 1 g, 5 g, 10 g		3-Bromo-9H-carbazole SKU: B1181 CAS Number: 1592-95-6 Options: • Quantity: 25 g, 50 g, 100 g, 250 g
	6,9-bis(5-bromo-4-(2-butyloctyl)thiophen-2-yl)dithieno[3,2-f:2',3'-h]quinoxaline SKU: B1111 CAS Number: 2754291-36-4 Options: • Quantity: 100 mg, 250 mg		NTCDA - 1,4,5,8-Naphthalenetetracarboxylic dianhydride SKU: B1191 CAS Number: 81-30-1 Options: • Quantity: 100 g, 250 g

















	Bromothiophene Benzothiadiazole SKU: B121 CAS Number: 288071-87-4 Options: • Quantity: 1 g, 5 g		5,5'-dibromo-3,3'-difluoro-2,2'-bithiophene SKU: B1381 CAS Number: 1619967-08-6 Options: • Quantity: 500 mg, 1 g, 5 g
	2-Bromodibenzothiophene SKU: B1221 CAS Number: 22439-61-8 Options: • Quantity: 25 g, 50 g, 100 g		2,2'-Bipyridine-4,4'-dicarboxylic acid SKU: B1391 CAS Number: 6813-38-3 Options: • Quantity: 5 g, 10 g, 25 g
	2-Decylthiophene SKU: B1241 CAS Number: 24769-39-9 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		Thienothiophene, Thieno[3,2-b]thiophene SKU: B141 CAS Number: 251-41-2 Options: • Quantity: 1 g, 5 g, 10 g
	2-Dodecylthiophene SKU: B1261 CAS Number: 4861-61-4 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		2,2'-Bipyridine-5,5'-dicarboxylic acid SKU: B1411 CAS Number: 1802-30-8 Options: • Quantity: 5 g, 10 g, 25 g
	Tetrabromothiophene SKU: B1281 CAS Number: 3958-03-0 Options: • Quantity: 50 g, 100 g, 250 g		Indolo[2,3-a]carbazole SKU: B1421 CAS Number: 60511-85-5 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	9-(heptadecan-9-yl)-2,7-bis(4,4,5,5-tetra-tert-butylphenyl)-9H-fluorene SKU: B131 CAS Number: 958261-51-3 Options: • Quantity: 1 g, 2 g, 5 g		3,3'-Bicarbazole SKU: B1431 CAS Number: 1984-49-2 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	2,5-Dibromo-3,4-dinitrothiophene SKU: B1351 CAS Number: 52431-30-8 Options: • Quantity: 5 g, 10 g, 25 g		5,5'-Bis(trifluoromethyl)-2,2'-bipyridine SKU: B1451 CAS Number: 142946-80-3 Options: • Quantity: 1 g, 2 g, 5 g
	4,4'-Dibromo-2,2'-bipyridine SKU: B1361 CAS Number: 18511-71-2 Options: • Quantity: 5 g, 10 g, 25 g		5,5'-Dimethyl-2,2'-bipyridine SKU: B1471 CAS Number: 1762-34-1 Options: • Quantity: 10 g, 25 g















	4-(Trifluoromethyl)-2-pyridone SKU: B1481 CAS Number: 50650-59-4 Options: • Quantity: 1 g, 5 g, 10 g		5,5'-Dibromo-4,4'-dihexyl-2,2'-bithiophene SKU: B1551 CAS Number: 214493-03-5 Options: • Quantity: 1 g, 2 g
	4,4'-Dimethyl-2,2'-bipyridyl SKU: B1491 CAS Number: 1134-35-6 Options: • Quantity: 10 g, 25 g		3,3-Difluoroazetidine hydrochloride SKU: B1561 CAS Number: 288315-03-7 Options: • Quantity: 5 g
	4-Fluoro-1H-indazole SKU: B1501 CAS Number: 341-23-1 Options: • Quantity: 1 g, 5 g		2-(Difluoromethyl)benzimidazole SKU: B1571 CAS Number: 705-09-9 Options: • Quantity: 5 g, 10 g
	PCE11 Monomer SKU: B1511 CAS Number: 1504626-07-6 Options: • Quantity: 1 g, 2 g		6-Fluoro-4-chromanone SKU: B1581 CAS Number: 66892-34-0 Options: • Quantity: 10 g, 25 g
	3-(4-Methyl-1H-imidazol-1-yl)-5-(trifluoromethyl)aniline SKU: B1511 CAS Number: 641571-11-1 Options: • Quantity: 5 g, 10 g		2-Fluoro-9H-carbazole SKU: B1591 CAS Number: 391-53-7 Options: • Quantity: 1 g, 2 g
	2-Bromo-5-decylthiophene SKU: B1521 CAS Number: 514188-72-8 Options: • Quantity: 1 g, 5 g, 10 g		7-Hydroxy-4-(trifluoromethyl)coumarin SKU: B1601 CAS Number: 575-03-1 Options: • Quantity: 1 g, 5 g
	2-Bromo-5-dodecylthiophene SKU: B1531 CAS Number: 153561-74-1 Options: • Quantity: 1 g, 5 g, 10 g		7-(Diethylamino)-4-(trifluoromethyl)coumarin SKU: B1611 CAS Number: 41934-47-8 Options: • Quantity: 250 mg, 500 mg
	5,5'-Dibromo-2,2'-bipyridine SKU: B1541 CAS Number: 15862-18-7 Options: • Quantity: 1 g, 5 g, 10 g		2-Furoyltrifluoroacetone SKU: B1621 CAS Number: 326-90-9 Options: • Quantity: 25 g, 50 g

















	2-(Trifluoromethyl)thioxanthen-9-one SKU: B1631 CAS Number: 1693-28-3 Options: • Quantity: 10 g		2,5-bis(trimethylstannyl)-thieno[3,2-b]thiophene SKU: B171 CAS Number: 469912-82-1 Options: • Quantity: 1 g, 5 g
	5-(Trifluoromethyl)-2-pyridone SKU: B1641 CAS Number: 33252-63-0 Options: • Quantity: 25 g, 50 g		5-Fluoroindole SKU: B1711 CAS Number: 399-52-0 Options: • Quantity: 5 g, 10 g
	3-(Trifluoromethyl)-2-pyridone SKU: B1651 CAS Number: 22245-83-6 Options: • Quantity: 10 g, 25 g, 50 g		6-Fluoroindole SKU: B1721 CAS Number: 399-51-9 Options: • Quantity: 10 g, 25 g, 50 g
	1-[2-(Trifluoromethyl)phenyl]imidazole SKU: B1641 CAS Number: 25371-96-4 Options: • Quantity: 1 g, 2 g		6-Fluoroindole-2-carboxylic acid SKU: B1741 CAS Number: 3093-97-8 Options: • Quantity: 1 g, 2 g
	4-(Trifluoromethyl)-1H-imidazole SKU: B1671 CAS Number: 3468-69-8 Options: • Quantity: 5 g, 10 g		5-(Trifluoromethoxy)isatin SKU: B1741 CAS Number: 169037-23-4 Options: • Quantity: 5 g, 10 g
	4-Fluoroindole SKU: B1681 CAS Number: 387-43-9 Options: • Quantity: 5 g, 10 g		5-Fluoroisatin SKU: B1751 CAS Number: 56341-41-4 Options: • Quantity: 10 g, 25 g
	3-(4-Fluorophenyl)-1-isopropylindole SKU: B1681 CAS Number: 93957-49-4 Options: • Quantity: 25 g, 50 g		5-Fluoroisatin SKU: B1761 CAS Number: 443-69-6 Options: • Quantity: 25 g, 50 g
	7-Fluoroindole SKU: B1701 CAS Number: 387-44-0 Options: • Quantity: 5 g, 10 g		Benzyl (3-Fluoro-4-morpholinophenyl)carbamate SKU: B1771 CAS Number: 168828-81-7 Options: • Quantity: 25 g, 50 g

















	3-Fluoro-4-morpholinoaniline SKU: B1781 CAS Number: 93246-53-8 Options: • Quantity: 25 g, 50 g		4-Bromo-8-fluoroquinoline SKU: B1891 CAS Number: 927800-38-2 Options: • Quantity: 500 mg, 1 g, 2 g
	7-Chloro-1-cyclopropyl-6-fluoro-1,4-dihydro-4-oxo-1,8-naphthyridine-3-carboxylic acid SKU: B1791 CAS Number: 100361-18-0 Options: • Quantity: 25 g, 50 g		4-Chloro-7-trifluoromethylquinoline-3-carboxylic acid ethyl ester SKU: B1901 CAS Number: 1168-42-3 Options: • Quantity: 500 mg, 1 g, 2 g
	2,4,6-Trimethyl-N,N-bis(4-(4,4,5,5-tetra-tert-butylphenyl)phenyl)aniline SKU: B1801 CAS Number: 1282616-14-1 Options: • Quantity: 1 g, 2 g, 5 g, 10 g		2-(2,4-Difluorophenyl)-5-methylpyridine SKU: B1911 CAS Number: 583052-21-5 Options: • Quantity: 1 g, 2 g, 5 g
	2,1,3-Benzothiadiazole-4,7-bis(boronic acid pinacol ester) SKU: B181 CAS Number: 934365-16-9 Options: • Quantity: 1 g, 5 g		2-(2,4-Difluorophenyl)-5-fluoropyridine SKU: B1921 CAS Number: 1426047-01-9 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g
	4-Butyl-N,N-bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolane-4-phenyl)aniline SKU: B1821 CAS Number: 444289-55-8 Options: • Quantity: 1 g, 2 g, 5 g, 10 g		3,3',4,4'-Benzophenonetetracarboxylic dianhydride (BTDA) SKU: B1931 CAS Number: 2421-28-5 Options: • Quantity: 100 g, 200 g, 500 g
	(5S)-(-)-5-(Trifluoromethyl)-2-pyrrolidinone SKU: B1881 CAS Number: 1287211-10-2 Options: • Quantity: 1 g, 2 g, 5 g		4,4'-Oxydiphthalic anhydride (ODPA) SKU: B1961 CAS Number: 1823-59-2 Options: • Quantity: 10 g, 25 g, 50 g
	8-Fluoro-4-hydroxy-2-methylquinoline SKU: B1871 CAS Number: 5288-22-2 Options: • Quantity: 1 g, 2 g, 5 g		4,4'-(4,4'-isopropylidenediphenoxy)bis(p-ththalic anhydride) (BPADA) SKU: B1991 CAS Number: 38103-06-9 Options: • Quantity: 50 g, 100 g, 200 g
	6-Chloro-8-fluoroquinoline SKU: B1881 CAS Number: 52200-53-0 Options: • Quantity: 1 g, 5 g, 10 g		1,2,4,5-cyclohexanetetracarboxylic dianhydride (HPMDA) SKU: B2041 CAS Number: 2754-41-8 Options: • Quantity: 5 g, 10 g, 25 g

















	1,2,3,4-Butanetetrahydrophthalic dianhydride (BDA) SKU: B2101 CAS Number: 4534-73-0 Options: • Quantity: 25 g, 50 g, 100 g		5-Bromo-2-fluoropyridine SKU: B2221 CAS Number: 766-11-0 Options: • Quantity: 25 g, 50 g
	2,6-dibromo-4,4-bis(2-ethylhexyl)-4H-cyclopenta[1,2-b:5,4-b']dithiophene SKU: B211 CAS Number: 365547-21-3 Options: • Quantity: 1 g, 5 g		4,4'-Bis(trifluoromethyl)-2,2'-bipyridine SKU: B2231 CAS Number: 142946-79-0 Options: • Quantity: 500 mg, 1 g, 2 g
	Bicyclo(2.2.2)oct-7-ene-2,3,5,6-tetracarboxylic acid dianhydride (COeDA) SKU: B2111 CAS Number: 1719-83-1 Options: • Quantity: 25 g, 50 g, 100 g		3,3-Difluoropyrrolidine hydrochloride SKU: B2281 CAS Number: 163457-23-6 Options: • Quantity: 2 g, 5 g, 10 g
	2,3,3',4'-Biphenyltetracarboxylic dianhydride (a-BPDA) SKU: B2131 CAS Number: 36978-41-3 Options: • Quantity: 10 g, 25 g, 50 g		2-Fluoroadenine SKU: B2291 CAS Number: 700-49-2 Options: • Quantity: 1 g, 2 g, 5 g
	Dicyclohexyl-3,4,3',4'-tetracarboxylic dianhydride (HBPA) SKU: B2141 CAS Number: 122640-83-9 Options: • Quantity: 10 g, 25 g, 50 g		(R)-(-)-3-Fluoropyrrolidine hydrochloride SKU: B2311 CAS Number: 136725-55-8 Options: • Quantity: 1 g, 2 g, 5 g
	4-Phenylethynylphthalic anhydride (PEPA) SKU: B2151 CAS Number: 119389-05-8 Options: • Quantity: 10 g, 25 g, 50 g		(S)-(+)-3-Fluoropyrrolidine hydrochloride SKU: B2341 CAS Number: 136725-53-6 Options: • Quantity: 1 g, 2 g, 5 g
	6-Fluoro-3-(4-piperidinyl)-1,2-benzisoxazole hydrochloride SKU: B2201 CAS Number: 84163-13-3 Options: • Quantity: 10 g, 25 g, 50 g		N-Boc-trans-4-fluoro-L-proline SKU: B2351 CAS Number: 203866-14-2 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g
	5,5'-bis(trimethylstannyl)-2,2'-bithiophene SKU: B221 CAS Number: 143367-56-0 Options: • Quantity: 1 g, 5 g		2-Chloro-5-fluoropyrimidine SKU: B2391 CAS Number: 62802-42-0 Options: • Quantity: 5 g, 10 g


	3,6-bis(5-bromothiophen-2-yl)-2,5-bis(2-H,5H)-dione SKU: B241 CAS Number: 1260685-63-9 Options: • Quantity: 1 g, 2 g		3,3-Difluoropiperidine hydrochloride SKU: B2571 CAS Number: 496807-97-7 Options: • Quantity: 1 g, 2 g, 5 g
	2-Bromo-4-fluoropyridine SKU: B2411 CAS Number: 357927-50-5 Options: • Quantity: 5 g, 10 g		2-Amino-3-fluoropyridine SKU: B2601 CAS Number: 21717-95-3 Options: • Quantity: 2 g, 5 g, 10 g
	6-Fluoronicotinic acid SKU: B2421 CAS Number: 403-45-2 Options: • Quantity: 1 g, 2 g, 5 g		Benzo[1,2-b:4,5-b']dithiophene-4,8-dione SKU: B261 CAS Number: 32281-36-0 Options: • Quantity: 5 g, 10 g, 25 g
	1-Bis(4-fluorophenyl)methyl piperazine SKU: B2441 CAS Number: 27469-60-9 Options: • Quantity: 10 g, 25 g		3-Fluoro-2-pyridinecarbonitrile SKU: B2611 CAS Number: 97509-75-6 Options: • Quantity: 5 g, 10 g, 25 g
	2-Fluoro-5-nitrobenzoic acid SKU: B2451 CAS Number: 7304-32-7 Options: • Quantity: 10 g, 25 g, 50 g		4H-Cyclopenta[1,2-b:5,4-b']dithiophene SKU: B271 CAS Number: 389-58-2 Options: • Quantity: 1 g, 5 g
	2,7-Dibromocarbazole SKU: B2521 CAS Number: 136630-39-2 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		5-Cyano-2-fluoropyridine SKU: B2771 CAS Number: 3939-12-6 Options: • Quantity: 1 g, 2 g, 5 g
	N-Boc-4,4-difluoro-L-proline SKU: B2521 CAS Number: 203866-15-3 Options: • Quantity: 500 mg, 1 g, 2 g		3,5-Difluoropyridine-2-carbonitrile SKU: B2781 CAS Number: 298709-29-2 Options: • Quantity: 5 g, 10 g
	2,6-Dichloro-5-fluoro-3-pyridinecarbonitrile SKU: B2541 CAS Number: 82671-02-1 Options: • Quantity: 25 g, 50 g		Thieno[3,2-b]thiophene-2-carbonitrile SKU: B281 CAS Number: 40985-58-8 Options: • Quantity: 1 g, 2 g, 5 g

	2-Fluoropyridine-3-carbonitrile SKU: B2861 CAS Number: 3939-13-7 Options: • Quantity: 10 g, 25 g, 50 g		3-Hexylthiophene SKU: B301 CAS Number: 1693-86-3 Options: • Quantity: 25 g
	6-Fluoropyridine-2-carbonitrile SKU: B2901 CAS Number: 3939-15-9 Options: • Quantity: 1 g, 2 g, 5 g		2,2'-bithiophene-5-carbonitrile SKU: B311 CAS Number: 16278-99-2 Options: • Quantity: 5 g
	5-Bromo-2-fluoropyrimidine SKU: B2921 CAS Number: 62802-38-4 Options: • Quantity: 1 g, 2 g, 5 g		2,2'-bithiophene SKU: B321 CAS Number: 492-97-7 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	1-(3-Fluorophenyl)-1H-pyrazole-4-carbaldehyde SKU: B2951 CAS Number: 936940-82-8 Options: • Quantity: 1 g, 2 g		2,6-Dibromonaphthalene-1,4,5,8-tetracarboxylic Dianhydride SKU: B3211 CAS Number: 83204-68-6 Options: • Quantity: 1 g, 5 g, 10 g
	4-Amino-2-(trifluoromethyl)pyridine SKU: B2961 CAS Number: 147149-98-2 Options: • Quantity: 1 g, 2 g, 5 g		3-Fluoro-4-iodopyridine SKU: B3261 CAS Number: 22282-75-3 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g
	4-Bromo-6-fluoroquinoline SKU: B2971 CAS Number: 661463-17-8 Options: • Quantity: 500 mg, 1 g, 2 g		2-Amino-5-chloro-3-fluoropyridine SKU: B3391 CAS Number: 246847-98-3 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	6-Fluoro-4-hydroxy-2-(trifluoromethyl)quinoline SKU: B2981 CAS Number: 31009-34-4 Options: • Quantity: 500 mg, 1 g, 2 g		7-Fluoro-4-hydroxycoumarin SKU: B3401 CAS Number: 2145-27-9 Options: • Quantity: 100 mg, 250 mg, 500 mg
	5-Difluoromethoxy-2-mercaptobenzimidazole SKU: B2991 CAS Number: 97963-62-7 Options: • Quantity: 25 g, 50 g		5-Bromo-3-fluoropyridine-2-carboxylic acid SKU: B3411 CAS Number: 669066-91-5 Options: • Quantity: 1 g, 2 g, 5 g


	5-Fluorobenzoxazole-2-thiol SKU: B3421 CAS Number: 13451-78-0 Options: • Quantity: 500 mg, 1 g, 2 g, 5 g		DFBTA-2ThBr-HD SKU: B421 CAS Number: 1887135-96-7 Options: • Quantity: 500 mg, 1 g
	BDTT SKU: B351 CAS Number: 1352642-37-5 Options: • Quantity: 1 g, 5 g		NDI-2EH-2ThBr SKU: B4231 CAS Number: 1320282-67-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	BDD SKU: B371 CAS Number: 1415929-78-0 Options: • Quantity: 500 mg, 1 g, 5 g, 10 g, 25 g		BDT-Th(TPS)-Sn SKU: B431 CAS Number: 2035466-88-5 Options: • Quantity: 500 mg, 1 g
	Benzodithiophene-2ThBr-4EH SKU: B381 CAS Number: 1439937-07-1 Options: • Quantity: 250 mg, 500 mg, 1 g		Quinoxaline-2PhO-2Br SKU: B451 CAS Number: 565228-31-1 Options: • Quantity: 250 mg, 500 mg, 1 g
	TPD-C8 SKU: B401 CAS Number: 566939-58-0 Options: • Quantity: 1 g		DPPDPy2Br SKU: B461 CAS Number: 1455028-34-8 Options: • Quantity: 100 mg, 250 mg, 1 g
	2,4,6-Tris(4-hydrazinylphenyl)-1,3,5-triazine SKU: B4031 CAS Number: 804548-56-9 Options: • Quantity: 100 mg, 250 mg, 500 mg		3,3',6,6'-Tetrabromo-9,9'-bicarbazole SKU: B4661 CAS Number: 18628-03-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	NDI-2OD-2Br SKU: B411 CAS Number: 1100243-35-3 Options: • Quantity: 250 mg, 500 mg, 1 g		DPP-EH-2ThBr SKU: B471 CAS Number: 1000623-95-9 Options: • Quantity: 1 g, 5 g
	9-(4-Ethynylphenyl)carbazole SKU: B4151 CAS Number: 262861-81-4 Options: • Quantity: 250 mg, 500 mg, 1 g		DPP-HD-2ThBr SKU: B481 CAS Number: 1000623-98-2 Options: • Quantity: 1 g, 5 g

	DPP-BO-2ThBr SKU: B491 CAS Number: 1224709-68-5 Options: • Quantity: 1 g		DTS-EH-Sn SKU: B641 CAS Number: 1089687-06-8 Options: • Quantity: 1 g, 5 g
	NDI-2HD-2Br SKU: B501 CAS Number: 1459168-68-3 Options: • Quantity: 250 mg, 500 mg, 1 g		D18-2Br SKU: B651 CAS Number: 2433725-51-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BDTT-Cl-Sn SKU: B531 CAS Number: 2239295-69-1 Options: • Quantity: 250 mg, 500 mg, 1 g		PTRZ-2Cl SKU: B671 CAS Number: 1700-02-3 Options: • Quantity: 10 g, 25 g, 100 g
	IDTB6-2CHO SKU: B531 CAS Number: 1884694-93-2 Options: • Quantity: 500 mg, 1 g		BPPTRZ-2Cl SKU: B681 CAS Number: 1472062-94-4 Options: • Quantity: 5 g, 25 g, 100 g
	IDTB6-2CHO SKU: B571 CAS Number: 1878125-76-8 Options: • Quantity: 500 mg, 1 g		ICz-Br SKU: B691 CAS Number: 1174032-81-5 Options: • Quantity: 1 g, 5 g, 10 g
	EHThCHO-Br SKU: B591 CAS Number: 1356861-54-5 Options: • Quantity: 1 g, 5 g		CDT-EH-Sn SKU: B731 CAS Number: 920504-00-3 Options: • Quantity: 1 g, 5 g
	EHOTHCHO-Br SKU: B601 CAS Number: 2055812-54-7 Options: • Quantity: 1 g, 5 g		BBT-2C6ThBr SKU: B741 Options: • Quantity: 250 mg, 500 mg, 1 g
	BDTOEHSn SKU: B621 CAS Number: 1160823-78-8 Options: • Quantity: 1 g, 5 g		Ethyl 2,5-dibromothiophene-3-carboxylate SKU: B751 CAS Number: 289470-44-6 Options: • Quantity: 1 g, 5 g, 10 g, 25 g


	2,5-Dibromo-3-chlorothiophene SKU: B761 CAS Number: 32431-91-7 Options: • Quantity: 1 g, 5 g, 10 g, 25 g		3TCN-HD-2Br SKU: B881 CAS Number: 2767223-03-8 Options: • Quantity: 500 mg, 1 g
	HDOQx-DBrDF SKU: B781 CAS Number: 2269476-12-0 Options: • Quantity: 250 mg, 500 mg, 1 g		BT-C12-2Br SKU: B891 CAS Number: 753470-95-0 Options: • Quantity: 1 g, 5 g, 10 g
	BT-C14-2Br SKU: B791 CAS Number: 888491-16-5 Options: • Quantity: 1 g, 5 g, 10 g, 25 g		BT-HD-2Br SKU: B901 Options: • Quantity: 500 mg, 1 g
	Br-BT-CHO SKU: B811 CAS Number: 1071224-34-4 Options: • Quantity: 1 g, 5 g, 10 g		5,5'-Dibromo-2,2'-bithiophene SKU: B921 CAS Number: 4805-22-5 Options: • Quantity: 5 g, 10 g, 25 g
	IDTB6-2Sn SKU: B841 CAS Number: 1252555-61-5 Options: • Quantity: 500 mg, 1 g, 5 g		DTT, Dithieno[3,2-b:2',3'-d]thiophene SKU: B931 CAS Number: 3593-75-7 Options: • Quantity: 1 g, 5 g
	IDTB6-2Br SKU: B851 CAS Number: 1049034-71-0 Options: • Quantity: 250 mg, 500 mg, 1 g		3,3'-Dibromo-2,2'-bithiophene SKU: B941 CAS Number: 51751-44-1 Options: • Quantity: 5 g, 10 g, 25 g
	IDTTB6-2Sn SKU: B861 CAS Number: 1420071-65-3 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		3,6-Dibromothieno[3,2-b]thiophene (TT36) SKU: B951 CAS Number: 392662-65-6 Options: • Quantity: 5 g, 10 g, 25 g
	IDTTB6-2Br SKU: B871 CAS Number: 2227361-27-3 Options: • Quantity: 250 mg, 500 mg, 1 g		2,5-Dihydro-3,6-di-2-thienyl-pyrrolo[3,4-c]pyrrole-1,4-dione SKU: B961 CAS Number: 850583-75-4 Options: • Quantity: 5 g, 10 g, 25 g



6,6'-Dibromoisindigo
 SKU: B971 | CAS Number: 1147124-21-7
 Options:
 • Quantity: 1 g, 5 g, 10 g


















2,8-Dibromodibenzothiophene
 SKU: B981 | CAS Number: 31574-87-5
 Options:
 • Quantity: 10 g, 25 g, 50 g



































3,7-Dibromodibenzothiophene 5,5-dioxide
 SKU: B991 | CAS Number: 83834-12-2
 Options:
 • Quantity: 250 mg, 500 mg, 1 g

Non-Heterocyclic Building Blocks

	2,7-Dibromo-9,9'-spirobifluorene SKU: B1001 CAS Number: 171408-84-7 Options: • Quantity: 10 g, 25 g, 50 g, 100 g		2,3-Dibromonaphthalene SKU: B1251 CAS Number: 13214-70-5 Options: • Quantity: 1 g, 5 g, 10 g, 25 g
	2,2',7,7'-Tetrabromo-9,9'-spirobifluorene SKU: B1011 CAS Number: 128055-74-3 Options: • Quantity: 10 g, 25 g, 50 g, 100 g		1,8-Dibromonaphthalene SKU: B1271 CAS Number: 17135-74-9 Options: • Quantity: 10 g, 25 g, 50 g
	9,9'-Spirobifluorene SKU: B1021 CAS Number: 159-66-0 Options: • Quantity: 10 g, 25 g, 50 g, 100 g		9,10-Dibromoanthracene SKU: B1291 CAS Number: 523-27-3 Options: • Quantity: 50 g, 100 g, 250 g
	2,7-Dibromophenanthrene-9,10-dione SKU: B1011 CAS Number: 84405-44-7 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		2,7-Dibromonaphthalene SKU: B1301 CAS Number: 58556-75-5 Options: • Quantity: 10 g, 25 g, 50 g
	N,N-Bis(4-bromophenyl)-2,4,6-trimethyl-aniline SKU: B1101 CAS Number: 663943-27-9 Options: • Quantity: 1 g, 5 g, 10 g		9-Bromoanthracene SKU: B1311 CAS Number: 1564-64-3 Options: • Quantity: 10 g, 25 g, 50 g, 100 g
	1,4-Dibromonaphthalene SKU: B1201 CAS Number: 83-53-4 Options: • Quantity: 50 g, 100 g, 250 g		2,6-Dibromonaphthalene SKU: B1321 CAS Number: 13720-06-4 Options: • Quantity: 5 g, 10 g, 25 g
	2-Bromonaphthalene SKU: B1211 CAS Number: 580-13-2 Options: • Quantity: 50 g, 100 g, 250 g		2-Bromo-6-fluoronaphthalene SKU: B1331 CAS Number: 324-41-4 Options: • Quantity: 1 g, 5 g
	1,5-Dibromonaphthalene SKU: B1231 CAS Number: 7351-74-8 Options: • Quantity: 1 g, 5 g, 10 g		2-Bromoanthracene SKU: B1341 CAS Number: 7321-27-9 Options: • Quantity: 10 g, 25 g, 50 g

	4,4'-Dibromo-2,2'-dinitrophenyl SKU: B1371 CAS Number: 91371-12-9 Options: • Quantity: 5 g, 10 g, 25 g		4,4'-(1,3-Phenylenedioxy)dianiline (TPE-R) SKU: B2011 CAS Number: 2479-46-1 Options: • Quantity: 50 g, 100 g, 200 g
	2,6-Dihydroxyanthraquinone SKU: B1401 CAS Number: 84-60-6 Options: • Quantity: 10 g, 25 g, 50 g		2,2'-Dimethyl[1,1'-biphenyl]-4,4'-diamine SKU: B2021 CAS Number: 84-67-3 Options: • Quantity: 50 g, 100 g, 200 g
	2,6-Dibromoanthraquinone SKU: B1461 CAS Number: 633-70-5 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		3,4'-Oxydianiline (3,4'-ODA) SKU: B2031 CAS Number: 2657-87-6 Options: • Quantity: 25 g, 50 g, 100 g
	2,7-Dibromo-9,9-bis[3,3'-(N,N-dimethylamino)propyl]fluorene SKU: B191 CAS Number: 673474-73-2 Options: • Quantity: 1 g, 5 g		4,4'-Bis(4-aminophenoxy)biphenyl (BAPB) SKU: B2051 CAS Number: 13080-85-8 Options: • Quantity: 25 g, 50 g, 100 g
	2,2-Bis[4-(4-aminophenoxy)phenyl]propane (BAPP) SKU: B1941 CAS Number: 13080-86-9 Options: • Quantity: 50 g, 100 g, 200 g, 500 g		4,4'-Bis(3-aminophenoxy)biphenyl (43BAPOBP) SKU: B2121 CAS Number: 105112-76-3 Options: • Quantity: 10 g, 25 g, 50 g
	1,4-Bis(4-aminophenoxy)benzene (TPE-Q) SKU: B1951 CAS Number: 3491-12-1 Options: • Quantity: 25 g, 50 g, 100 g		4-Bromo-N-(4-bromophenyl)-N-(4-sec-butylphenyl)aniline SKU: B2191 CAS Number: 287976-94-7 Options: • Quantity: 1 g, 2 g, 5 g
	9,9-Bis(4-aminophenyl)fluorene (FDA) SKU: B1981 CAS Number: 15499-84-0 Options: • Quantity: 50 g, 100 g, 200 g		4-Bromo-2,6-difluoroaniline SKU: B2301 CAS Number: 67567-26-4 Options: • Quantity: 5 g, 10 g, 25 g
	2,2'-(9,9-diethyl-9H-fluorene-2,7-diyl)bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolane) SKU: B201 CAS Number: 196207-58-6 Options: • Quantity: 5 g		5-Fluorosalicylic acid SKU: B2321 CAS Number: 345-16-4 Options: • Quantity: 10 g, 25 g, 50 g

	3,4,5-Trifluorobenzoic acid SKU: B2511 CAS Number: 121602-93-5 Options: • Quantity: 5 g, 10 g, 25 g		3-Fluorobenzoic hydrazide SKU: B2721 CAS Number: 499-55-8 Options: • Quantity: 10 g, 25 g
	2-Fluoro-4-(trifluoromethyl)aniline SKU: B2531 CAS Number: 69409-98-9 Options: • Quantity: 5 g, 10 g		2-Bromo-5-fluorobenzonitrile SKU: B2851 CAS Number: 57381-39-2 Options: • Quantity: 10 g, 25 g, 50 g
	4-Fluoro-β-nitrostyrene SKU: B2551 CAS Number: 706-08-1 Options: • Quantity: 5 g, 10 g		4-Bromo-3-fluorobenzonitrile SKU: B3241 CAS Number: 133059-44-6 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	2,2-Difluorocyclopropanecarboxylic acid SKU: B2581 CAS Number: 107873-03-0 Options: • Quantity: 1 g, 2 g, 5 g		2-Chloro-5-fluoronitrobenzene SKU: B3251 CAS Number: 345-17-5 Options: • Quantity: 2 g, 5 g, 10 g, 25 g
	α,α-Difluorophenylacetic acid SKU: B2581 CAS Number: 360-03-2 Options: • Quantity: 2 g, 5 g, 10 g		4-Fluoro-2-iodobenzoic acid SKU: B3271 CAS Number: 56096-89-0 Options: • Quantity: 1 g, 2 g, 5 g, 10 g
	2-Bromo-4-fluorobenzonitrile SKU: B2621 CAS Number: 36282-26-5 Options: • Quantity: 10 g, 25 g, 50 g		trans-3,4-Difluorocinnamic acid SKU: B3281 CAS Number: 112897-97-9 Options: • Quantity: 10 g, 25 g, 50 g
	3'-Fluoro-4'-methoxyacetophenone SKU: B2631 CAS Number: 455-91-4 Options: • Quantity: 25 g		2-Fluoro-5-methylbenzoic acid SKU: B3291 CAS Number: 321-12-0 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	2-Amino-5-fluorobenzonitrile SKU: B2711 CAS Number: 61272-77-3 Options: • Quantity: 5 g, 10 g		5-Chloro-2-fluorobenzoic acid SKU: B3301 CAS Number: 394-30-9 Options: • Quantity: 2 g, 5 g, 10 g, 25 g


	9,9-Dioctyl-2,7-dibromofluorene SKU: B331 CAS Number: 198964-46-4 Options: • Quantity: 10 g		5-Fluoro-2-methylbenzoic acid SKU: B3381 CAS Number: 33184-16-6 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	4-Fluoro-3-methylbenzoic acid SKU: B3311 CAS Number: 403-15-6 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		4-Fluoro-1,2-phenylenediamine SKU: B3431 CAS Number: 367-31-7 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	4'-Fluoro-3'-nitroacetophenone SKU: B3311 CAS Number: 400-93-1 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		2,4,5-Trifluorophenylacetic acid SKU: B3441 CAS Number: 209995-38-0 Options: • Quantity: 2 g, 5 g, 10 g, 25 g
	2-Fluoro-6-methylbenzoic acid SKU: B3331 CAS Number: 90259-27-1 Options: • Quantity: 2 g, 5 g, 10 g		Pentafluorobenzoic acid SKU: B3451 CAS Number: 602-94-8 Options: • Quantity: 25 g, 50 g, 100 g
	5-Chloro-2-fluorobenzonitrile SKU: B3341 CAS Number: 57381-34-7 Options: • Quantity: 2 g, 5 g, 10 g, 25 g		2-Fluorobiphenyl SKU: B3461 CAS Number: 321-60-8 Options: • Quantity: 10 g, 25 g, 50 g, 100 g
	3-Fluoro-5-(trifluoromethyl)benzoic acid SKU: B3351 CAS Number: 161622-05-5 Options: • Quantity: 2 g, 5 g, 10 g, 25 g		Potassium nonafluoro-1-butanedisulfonate SKU: B3471 CAS Number: 29420-49-3 Options: • Quantity: 25 g, 50 g, 100 g
	3'-Bromo-4'-fluoroacetophenone SKU: B3361 CAS Number: 1007-15-4 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		1,3,5-Tris(4-hydroxyphenyl)benzene SKU: B4121 CAS Number: 15797-52-1 Options: • Quantity: 250 mg, 500 mg, 1 g
	5-Bromo-4-fluoro-2-methylaniline SKU: B3371 CAS Number: 627871-16-3 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		1,2-Bis(4-ethynylphenyl)-1,2-diphenylethane SKU: B4131 CAS Number: 1240785-42-5 Options: • Quantity: 250 mg, 500 mg, 1 g

	2,2'-(1,2-Phenylene)bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolane) SKU: B4141 CAS Number: 269410-07-3 Options: • Quantity: 250 mg, 500 mg, 1 g		2HIC, 3-(Dicyanomethylidene)indan-1-one SKU: B541 CAS Number: 1080-74-6 Options: • Quantity: 5 g, 10 g, 25 g
	[1,1'3',1''-Terphenyl]-4,4'',5',5''-tricarbaldehyde SKU: B4161 CAS Number: 1977585-07-1 Options: • Quantity: 250 mg, 500 mg, 1 g		2CIIC SKU: B581 CAS Number: 2197167-50-1 Options: • Quantity: 500 mg, 1 g, 5 g
	2,5-Bis(allyloxy)terephthalohydrazide SKU: B4171 CAS Number: 2227151-69-9 Options: • Quantity: 250 mg, 500 mg, 1 g		2-Bromotriphenylene SKU: B661 CAS Number: 19111-87-6 Options: • Quantity: 5 g, 25 g, 100 g
	Tetrakis(4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl)methane SKU: B4181 CAS Number: 875772-13-7 Options: • Quantity: 250 mg, 500 mg, 1 g		BrIC SKU: B581 CAS Number: 507484-47-1 Options: • Quantity: 500 mg, 1 g, 5 g
	4,4',4''-([1,1'-Biphenyl]-4,4'-diylbis(azacetyl))tetrabenzaldehyde SKU: B4191 CAS Number: 865448-72-2 Options: • Quantity: 250 mg, 500 mg, 1 g		CIIC SKU: B831 CAS Number: 2364382-35-2 / 507484-48-2 Options: • Quantity: 500 mg, 1 g, 5 g
	[1,1'-Biphenyl]-2,2'-dicarbonitrile SKU: B4201 CAS Number: 4341-02-0 Options: • Quantity: 250 mg, 500 mg, 1 g		1,4-Di-tert-butyl-2,5-dimethoxybenzene SKU: M2391 CAS Number: 7323-63-9 Options: • Quantity: 25 g, 100 g, 250 g, 500 g, 1 kg
	(4-(1,2,2-Triphenylvinyl)phenyl)methanol SKU: B4211 CAS Number: 1015082-83-3 Options: • Quantity: 500 mg, 1 g		
	2,3,6,7,10,11-Hexamethoxytriphenylene SKU: B4221 CAS Number: 808-57-1 Options: • Quantity: 1 g, 5 g		


Monomers

Monomers act as the fundamental building blocks for the synthesis of diverse polymeric materials. Their reactive sites enable the incorporation of side chains or additives, which optimize film morphology, interfacial properties, and overall device performance.

- Boronates
- Organotins
- Diamines
- Dibromos
- Dianhydrides
- Carbaldehydes




Popular Products



Benzodithiophene-2THBr-4EH

Bromide-capped monomer for polymers including PBDD4T and PBDD4T-2F


- For the synthesis of low band-gap polymer semiconductors



BDTOEHsSn

Monomer for low bandgap polymer such as PTB7

- Benzodithiophene (BDT) derivative
- Ethylhexyl side chains



D18-2BR

Bromide-capped monomer for low bandgap D18 polymer

- Rigid electron accepting dithienobenzothiadiazole (DTBT) core

To help you choose the right material, full technical specifications are available on our website.

Ossila

B621-250mg

Purity
MW 7

BDTOEHSn

4,8-Bis[(2-ethylhexyl)oxy]-2,6-bis[
stannyl]benzo[1,2-b:4,5-b']dith




CAS# 1160823-78-8

Solpro Business Park Windsor













Fatal if swallowed, in contact with skin or on clothing. Very toxic to aquatic life. Avoid release to the environment. POISON CENTER/ doctor. Immediately call a person to fresh air and water. POISON CENTER/ doctor.

Boronates
















	9-(heptadecan-9-yl)-2,7-bis(4,4,5,5-tetra-4-bazole) SKU: B131 CAS Number: 958261-51-3 Options: • Quantity: 1 g, 2 g, 5 g		5'-(4,4,5,5-Tetramethyl-1,3,2-dioxaborol-an-2-yl)-[1,1':3'',1''-terphenyl]-4,4''-diamine SKU: B3661 CAS Number: 2754262-59-2 Options: • Quantity: 500 mg, 1 g, 5 g
	2,4,6-Trimethyl-N,N-bis(4-(4,4,5,5-tetra-4-bazole) SKU: B1801 CAS Number: 1282616-14-1 Options: • Quantity: 1 g, 2 g, 5 g, 10 g		Tris(4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl)amine SKU: B3681 CAS Number: 267221-90-9 Options: • Quantity: 1 g, 5 g, 10 g
	2,1,3-Benzothiadiazole-4,7-bis(boronic acid pinacol ester) SKU: B1801 CAS Number: 934365-16-9 Options: • Quantity: 1 g, 5 g		((1,3,5-Triazine-2,4,6-tryl)tris(benzene-4,1-diyl))triboronic acid SKU: B3731 CAS Number: 910231-21-9 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4-Butyl-N,N-bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolane-4-phenyl)aniline SKU: B1821 CAS Number: 444289-55-8 Options: • Quantity: 1 g, 2 g, 5 g, 10 g		3,5-Diformylphenylboronic acid SKU: B3751 CAS Number: 480424-62-2 Options: • Quantity: 500 mg, 1 g, 5 g
	2,2'-(9,9-dioctyl-9H-fluorene-2,7-diyl)bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolane) SKU: B201 CAS Number: 196207-58-6 Options: • Quantity: 5 g		1,3,5-Tris(4-phenylboronic acid)benzene SKU: B3931 CAS Number: 900795-73-5 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	2,5-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)terephthalaldehyde SKU: B3631 CAS Number: 2305615-53-4 Options: • Quantity: 500 mg, 1 g, 5 g		2,4,6-Tris(4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl)-1,3,5-triazine SKU: B3981 CAS Number: 1447947-87-6 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4,4'-(4-(4,4,5,5-Tetramethyl-1,3,2-dioxahyde SKU: B3641 CAS Number: 2009169-65-5 Options: • Quantity: 500 mg, 1 g, 5 g		4,4,5,5-Tetramethyl-2-[4-(triphenylvinyl)-phenyl]-1,3,2-dioxaborolane SKU: B4061 CAS Number: 1260865-91-5 Options: • Quantity: 500 mg, 1 g
	5'-(4,4,5,5-Tetramethyl-1,3,2-dioxaborol-aldehyde SKU: B3651 CAS Number: 1033752-94-1 Options: • Quantity: 500 mg, 1 g, 5 g		2,2'-(1,2-Phenylene)bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolane) SKU: B4141 CAS Number: 269410-07-3 Options: • Quantity: 250 mg, 500 mg, 1 g

	Tetrakis(4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl)methane SKU: B4181 CAS Number: 875772-13-7 Options: Quantity: 250 mg, 500 mg, 1 g
	FL3N-2BAPE SKU: B4241 CAS Number: 953390-94-8 Options: Quantity: 500 mg, 1 g, 5 g
	1,3,6,8-Tetrakis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)pyrene SKU: B4631 CAS Number: 1398053-00-3 Options: Quantity: 500 mg, 1 g

















Organotin Compounds

















	2,5-bis(trimethylstannyl)-thieno[3,2-b]thiophene SKU: B171 CAS Number: 469912-82-1 Options: • Quantity: 1 g, 5 g		DTS-EH-Sn SKU: B641 CAS Number: 1089687-06-8 Options: • Quantity: 1 g, 5 g
	5,5'-bis(trimethylstannyl)-2,2'-bithiophene SKU: B221 CAS Number: 143367-56-0 Options: • Quantity: 1 g, 5 g		CDT-EH-Sn SKU: B731 CAS Number: 920504-00-3 Options: • Quantity: 1 g, 5 g
	BDTT SKU: B351 CAS Number: 1352642-37-5 Options: • Quantity: 1 g, 5 g		IDTB6-2Sn SKU: B841 CAS Number: 1252555-61-5 Options: • Quantity: 500 mg, 1 g, 5 g
	3,3'-Difluoro-5,5'-bis(trimethylstannyl)-2,2'-bithiophene SKU: B391 CAS Number: 1619967-09-7 Options: • Quantity: 250 mg, 500 mg, 1 g		IDTTB6-2Sn SKU: B861 CAS Number: 1420071-65-3 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	BDT-Th(TPS)-Sn SKU: B431 CAS Number: 2035466-88-5 Options: • Quantity: 500 mg, 1 g		
	BDTT-F-Sn SKU: B521 CAS Number: 1514905-25-9 Options: • Quantity: 250 mg, 500 mg, 1 g		
	BDTT-Cl-Sn SKU: B531 CAS Number: 2239295-69-1 Options: • Quantity: 250 mg, 500 mg, 1 g		
	BDTOEHSn SKU: B621 CAS Number: 1160823-78-8 Options: • Quantity: 1 g, 5 g		

Dibromo Monomers

















	2,7-Dibromo-9,9'-spirobifluorene SKU: B1001 CAS Number: 171408-84-7 Options: Quantity: 10 g, 25 g, 50 g, 100 g		1,4-Dibromonaphthalene SKU: B1201 CAS Number: 83-53-4 Options: Quantity: 50 g, 100 g, 250 g
	4,7-Dibromo-2,1,3-benzothiadiazole SKU: B101 CAS Number: 15155-41-6 Options: Quantity: 5 g, 10 g		Bromothiophene Benzothiadiazole SKU: B121 CAS Number: 288071-87-4 Options: Quantity: 1 g, 5 g
	2,7-Dibromophenanthrene-9,10-dione SKU: B1031 CAS Number: 84405-44-7 Options: Quantity: 5 g, 10 g, 25 g, 50 g		1,5-Dibromonaphthalene SKU: B1231 CAS Number: 7351-74-8 Options: Quantity: 1 g, 5 g, 10 g
	2,6-Dibromodithieno[3,2-b:2',3'-d]thiophene SKU: B1051 CAS Number: 67061-69-2 Options: Quantity: 500 mg, 1 g, 2 g		9,10-Dibromoanthracene SKU: B1291 CAS Number: 523-27-3 Options: Quantity: 50 g, 100 g, 250 g
	5,5''-Dibromo-2,2':5,2''-terthiophene SKU: B1071 CAS Number: 98057-08-0 Options: Quantity: 1 g, 5 g		2,7-Dibromonaphthalene SKU: B1301 CAS Number: 58556-75-5 Options: Quantity: 10 g, 25 g, 50 g
	N,N-Bis(4-bromophenyl)-2,4,6-trimethyl-aniline SKU: B1101 CAS Number: 663943-27-9 Options: Quantity: 1 g, 5 g, 10 g		2,6-Dibromonaphthalene SKU: B1321 CAS Number: 13720-06-4 Options: Quantity: 5 g, 10 g, 25 g
	6,9-bis(5-bromo-4-(2-butyloctyl)thiophen-2-yl)dithieno[3,2-f:2',3'-h]quinoxaline SKU: B1111 CAS Number: 2754291-36-4 Options: Quantity: 100 mg, 250 mg		2,5-Dibromo-3,4-dinitrothiophene SKU: B1351 CAS Number: 52431-30-8 Options: Quantity: 5 g, 10 g, 25 g
	5,5'-Dibromo-3,3'-bis(2-(2-methoxyethoxy)ethoxy)-2,2'-bithiophene SKU: B1121 CAS Number: 1974310-58-1 Options: Quantity: 500 mg, 1 g, 5 g		4,4'-Dibromo-2,2'-bipyridine SKU: B1361 CAS Number: 18511-71-2 Options: Quantity: 5 g, 10 g, 25 g

















	4,4'-Dibromo-2,2'-dinitrophenyl SKU: B1371 CAS Number: 91371-12-9 Options: • Quantity: 5 g, 10 g, 25 g		2,6-dibromo-4,4-bis(2-ethylhexyl)-4H-clopenta[1,2-b:5,4-b']dithiophene SKU: B211 CAS Number: 365547-21-3 Options: • Quantity: 1 g, 5 g
	5,5'-dibromo-3,3'-difluoro-2,2'-bithiophene SKU: B1381 CAS Number: 1619967-08-6 Options: • Quantity: 500 mg, 1 g, 5 g		4-Bromo-N-(4-bromophenyl)-N-(4-sec-butylphenyl)aniline SKU: B2191 CAS Number: 287976-94-7 Options: • Quantity: 1 g, 2 g, 5 g
	2,6-Dibromoanthraquinone SKU: B1461 CAS Number: 633-70-5 Options: • Quantity: 5 g, 10 g, 25 g, 50 g		2,6-Dibromonaphthalene-1,4,5,8-tetracarboxylic Dianhydride SKU: B3211 CAS Number: 83204-68-6 Options: • Quantity: 1 g, 5 g, 10 g
	PCE11 Monomer SKU: B1511 CAS Number: 1504626-07-6 Options: • Quantity: 1 g, 2 g		9,9-Dioctyl-2,7-dibromofluorene SKU: B3311 CAS Number: 198964-46-4 Options: • Quality: 10 g
	5,5'-Dibromo-2,2'-bipyridine SKU: B1541 CAS Number: 15862-18-7 Options: • Quantity: 1 g, 5 g, 10 g		2-Ethylhexyl 4,6-dibromo-3-fluorothiophene-3-carboxylate SKU: B3611 CAS Number: 1237479-38-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	5,5'-Dibromo-4,4'-dihexyl-2,2'-bithiophene SKU: B1551 CAS Number: 214493-03-5 Options: • Quantity: 1 g, 2 g		3,6-Dibromo-9-(2-ethylhexyl)-9H-carbazole SKU: B3671 CAS Number: 173063-52-0 Options: • Quantity: 1 g, 5 g, 10 g
	4,7-Dibromo-5,6-difluoro-2,1,3-benzothiadiazole SKU: B161 CAS Number: 1295502-53-2 Options: • Quantity: 1 g, 5 g		BDD SKU: B371 CAS Number: 1415929-78-0 Options: • Quantity: 500 mg, 1 g, 5 g, 10 g, 25 g
	2,7-Dibromo-9,9-bis[3,3'-(N,N-dimethylamino)propyl]fluorene SKU: B191 CAS Number: 673474-73-2 Options: • Quantity: 1 g, 5 g		2,3-Bis(4-bromophenyl)fumaric nitrile SKU: B3761 CAS Number: 82193-93-9 Options: • Quantity: 500 mg, 1 g, 5 g
















	Benzodithiophene-2THBr-4EH SKU: B381 CAS Number: 1439937-07-1 Options: • Quantity: 250 mg, 500 mg, 1 g		DPPDPy2Br SKU: B461 CAS Number: 1455028-34-8 Options: • Quantity: 100 mg, 250 mg, 1 g
	TPD-C8 SKU: B401 CAS Number: 566939-58-0 Options: • Quantity: 1 g		3,3',6,6'-Tetrabromo-9,9'-bicarbazole SKU: B4661 CAS Number: 18628-03-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	NDI-2OD-2Br SKU: B411 CAS Number: 1100243-35-3 Options: • Quantity: 250 mg, 500 mg, 1 g		3,7-Dibromo-10H-phenothiazine SKU: B4671 CAS Number: 21667-32-3 Options: • Quantity: 5 g, 10 g, 25 g
	DFBTA-2ThBr-HD SKU: B411 CAS Number: 1887135-96-7 Options: • Quantity: 500 mg, 1 g		DPP-EH-2ThBr SKU: B411 CAS Number: 1000623-95-9 Options: • Quantity: 1 g, 5 g
	NDI-2EH-2THBr SKU: B4231 CAS Number: 1320282-67-4 Options: • Quantity: 250 mg, 500 mg, 1 g		DPP-HD-2ThBr SKU: B481 CAS Number: 1000623-98-2 Options: • Quantity: 1 g, 5 g
	g344-BTBr SKU: B4251 CAS Number: 1404493-47-5 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g, 10 g		DPP-BO-2ThBr SKU: B491 CAS Number: 1224709-68-5 Options: • Quantity: 1 g
	2FT-2Sn SKU: B4252 CAS Number: 870718-97-1 Options: • Quantity: 250 mg, 500 mg, 1 g		NDI-2HD-2Br SKU: B501 CAS Number: 1459168-68-3 Options: • Quantity: 250 mg, 500 mg, 1 g
	Quinoxaline-2PhO-2Br SKU: B451 CAS Number: 565228-31-1 Options: • Quantity: 250 mg, 500 mg, 1 g		D18-2Br SKU: B651 CAS Number: 2433725-51-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

	BBT-2C6ThBr SKU: B741 Options: • Quantity: 250 mg, 500 mg, 1 g		3TCN-HD-2Br SKU: B881 CAS Number: 2767223-03-8 Options: • Quantity: 500 mg, 1 g
	Ethyl 2,5-dibromothiophene-3-carboxylate SKU: B751 CAS Number: 289470-44-6 Options: • Quantity: 1 g, 5 g, 10 g, 25 g		BT-C12-2Br SKU: B891 CAS Number: 753470-95-0 Options: • Quantity: 1 g, 5 g, 10 g
	2,5-Dibromo-3-chlorothiophene SKU: B761 CAS Number: 32431-91-7 Options: • Quantity: 1 g, 5 g, 10 g, 25 g		BT-HD-2Br SKU: B901 Options: • Quantity: 500 mg, 1 g
	MeHDOQx-DBrDF SKU: B751 CAS Number: 2414359-15-0 Options: • Quantity: 250 mg, 500 mg, 1 g		5,5'-Dibromo-2,2'-bithiophene SKU: B921 CAS Number: 4805-22-5 Options: • Quantity: 5 g, 10 g, 25 g
	HDOQx-DBrDF SKU: B781 CAS Number: 2269476-12-0 Options: • Quantity: 250 mg, 500 mg, 1 g		3,3'-Dibromo-2,2'-bithiophene SKU: B941 CAS Number: 51751-44-1 Options: • Quantity: 5 g, 10 g, 25 g
	BT-C14-2Br SKU: B791 CAS Number: 888491-16-5 Options: • Quantity: 1 g, 5 g, 10 g, 25 g		3,6-Dibromothieno[3,2-b]thiophene (TT36) SKU: B951 CAS Number: 392662-65-6 Options: • Quantity: 5 g, 10 g, 25 g
	IDTB6-2Br SKU: B851 CAS Number: 1049034-71-0 Options: • Quantity: 250 mg, 500 mg, 1 g		2,8-Dibromodibenzothiophene SKU: B981 CAS Number: 31574-87-5 Options: • Quantity: 10 g, 25 g, 50 g
	IDTTB6-2Br SKU: B871 CAS Number: 2227361-27-3 Options: • Quantity: 250 mg, 500 mg, 1 g		3,7-Dibromodibenzothiophene 5,5-dioxide SKU: B991 CAS Number: 83834-12-2 Options: • Quantity: 250 mg, 500 mg, 1 g

















Carbaldehyde Monomers















	4,4',4''-(Ethene-1,1,2,2-tetrayl)tetrabenzaldehyde SKU: B3481 CAS Number: 2170451-48-4 Options: • Quantity: 500 mg, 1 g, 5 g		2,5-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)terephthalaldehyde SKU: B3631 CAS Number: 2305615-53-4 Options: • Quantity: 500 mg, 1 g, 5 g
	[1,1'-Biphenyl]-3,4',5-tricarbaldehyde SKU: B3491 CAS Number: 187281-19-2 Options: • Quantity: 500 mg, 1 g, 5 g		4,4'-(4,4,5,5-Tetramethyl-1,3,2-dioxahydro-2,5-dithienyl)dibenzaldehyde SKU: B3641 CAS Number: 2009169-65-5 Options: • Quantity: 500 mg, 1 g, 5 g
	Hexa(4-formylphenyl)benzene SKU: B3521 CAS Number: 1862220-96-9 Options: • Quantity: 500 mg, 1 g, 5 g		5'-(4,4,5,5-Tetramethyl-1,3,2-dioxaborol-2-yl)-4,4'-biphenyl-2-carbaldehyde SKU: B3651 CAS Number: 1033752-94-1 Options: • Quantity: 500 mg, 1 g, 5 g
	[3,3-Bipyridine]-6,6-dicarboxaldehyde SKU: B3541 CAS Number: 1264748-06-2 Options: • Quantity: 500 mg, 1 g, 5 g		Tetrakis(4-formylphenyl)silane SKU: B3691 CAS Number: 1055999-34-2 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4,4'-(Ethyne-1,2-diyl)dibenzaldehyde SKU: B3571 CAS Number: 84907-55-1 Options: • Quantity: 500 mg, 1 g, 5 g		[1,1'-Biphenyl]-3,3',5,5'-tetracarbaldehyde SKU: B3701 CAS Number: 150443-85-9 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4,4'-(2,2-Diphenylethene-1,1-diyl)dibenzaldehyde SKU: B3581 CAS Number: 1601465-06-8 Options: • Quantity: 500 mg, 1 g, 5 g		3,5-Diformylphenylboronic acid SKU: B3751 CAS Number: 480424-62-2 Options: • Quantity: 500 mg, 1 g, 5 g
	5,5'-(Ethyne-1,2-diyl)diisophthalaldehyde SKU: B3611 CAS Number: 2237965-02-3 Options: • Quantity: 500 mg, 1 g, 5 g		1,2,4,5-Tetrakis(4-formylphenyl)benzene SKU: B3771 CAS Number: 883835-33-4 Options: • Quantity: 1 g, 5 g
	1,3,5-Tris(4-formyl-3-hydroxyphenyl)benzene SKU: B3621 CAS Number: 2244985-34-8 Options: • Quantity: 500 mg, 1 g, 5 g		4',4''',4''''-Nitrilotris([1,1'-biphenyl]-4-carbaldehyde) SKU: B3801 CAS Number: 872689-79-7 Options: • Quantity: 500 mg, 1 g, 5 g

	5-(4-(Diphenylamino)phenyl)thiophene-2-carbaldehyde SKU: B3811 CAS Number: 291279-14-6 Options: • Quantity: 1 g, 5 g		1,3,5-Tris([4-phenylboronic acid])benzene SKU: B3931 CAS Number: 900795-73-5 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4,4''-p-Terphenyldicarboxaldehyde SKU: B3821 CAS Number: 62940-38-9 Options: • Quantity: 1 g, 5 g		1,3,6,8-Tetrakis(4-formylphenyl)pyrene SKU: B3941 CAS Number: 1415238-25-3 Options: • Quantity: 250 mg, 500 mg, 1 g
	4,4',4'''-(Benzene-1,3,5-triyltris(ethyne-2,1-diyl))tribenzaldehyde SKU: B3831 CAS Number: 2289758-98-9 Options: • Quantity: 500 mg, 1 g, 5 g		2',5'-Dimethoxy-[1,1':4',1''-terphenyl]-4,4''-dicarbaldehyde SKU: B3951 CAS Number: 111759-27-4 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	9,10-Bis(4-formylphenyl)anthracene SKU: B3851 CAS Number: 324750-99-4 Options: • Quantity: 250 mg, 500 mg, 1 g		4,4',4'''-(pyridine-2,4,6-triyl)tribenzaldehyde SKU: B3961 CAS Number: 2230887-28-8 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4-(1,2,2-Triphenylvinyl)benzaldehyde SKU: B3881 CAS Number: 1289218-74-1 Options: • Quantity: 1 g, 5 g		2,5-Bis(octyloxy)terephthalaldehyde SKU: B4041 CAS Number: 123440-34-6 Options: • Quantity: 500 mg, 1 g
	4,4'-(Benzo[c][1,2,5]thiadiazole-4,7-diyl)di-ibenzaldehyde SKU: B3891 CAS Number: 914651-17-5 Options: • Quantity: 500 mg, 1 g		1,1':3',1''-Terphenyl-4,4''-dicarbaldehyde SKU: B4051 CAS Number: 171820-02-3 Options: • Quantity: 500 mg, 1 g
	5'-(4-Formylphenyl)-2',4',6'-trimethyl-[1,—1':3',1''-terphenyl]-4,4''-dicarbaldehyde SKU: B3911 CAS Number: 2195343-70-3 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		[1,1':4',1''-terphenyl]-3,3',5,5''-tetracarbaldehyde SKU: B4081 CAS Number: 2380275-62-5 Options: • Quantity: 250 mg, 500 mg, 1 g
	1,3,5-Tris(3-fluoro-4-formylphenyl)benzene SKU: B3921 CAS Number: 2088007-04-7 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		4',4''',4''''-Nitrilotris([1,1'-biphenyl]-3,5-dicarbaldehyde) SKU: B4091 CAS Number: 2030105-85-0 Options: • Quantity: 250 mg, 500 mg, 1 g

	[1,1':3',1''-Terphenyl]-4,4'',5'-tricarbaldehyde SKU: B4161 CAS Number: 1977585-07-1 Options: • Quantity: 250 mg, 500 mg, 1 g		4',4''',4''''-(1,3,5-Triazine-2,4,6-triyl)tris([1,—1'-biphenyl]-4-carbaldehyde) SKU: B4611 CAS Number: 1221509-80-3 Options: • Quantity: 250 mg, 500 mg, 1 g
	4,4',4'',4'''-([1,1'-Biphenyl]-4,4'-diylbis(azasetriyl))tetrabenzaldehyde SKU: B4191 CAS Number: 865448-72-2 Options: • Quantity: 250 mg, 500 mg, 1 g		IDTB6-2CHO SKU: B561 CAS Number: 1884694-93-2 Options: • Quantity: 500 mg, 1 g
	2,9-Bis[p-(formyl)phenyl]-1,10-phenanthroline SKU: B4521 CAS Number: 120085-99-6 Options: • Quantity: 500 mg, 1 g		IDTTB6-2CHO SKU: B571 CAS Number: 1878125-76-8 Options: • Quantity: 500 mg, 1 g
	5,5',5''-(Benzene-1,3,5-triyl)tripicolinaldehyde SKU: B4531 CAS Number: 1482413-54-6 Options: • Quantity: 250 mg, 500 mg, 1 g		EHTbCHO-Br SKU: B591 CAS Number: 1356861-54-5 Options: • Quantity: 1 g, 5 g
	Tris[4-(5-formyl-2-thienyl)phenyl]amine SKU: B4551 CAS Number: 883236-47-3 Options: • Quantity: 500 mg, 1 g		EHOTbCHO-Br SKU: B601 CAS Number: 2055812-54-7 Options: • Quantity: 1 g, 5 g
	[1,1':3',1''-Terphenyl]-3,3'',5,5''-tetracarbaldehyde SKU: B4571 CAS Number: 201734-76-1 Options: • Quantity: 500 mg, 1 g		TTP-TTI-2CHO SKU: B611 CAS Number: 2304444-53-7 Options: • Quantity: 100 mg, 250 mg
	4',4''',4''''-(Ethene-1,1,2,2-tetrayl)tetrakis([1,1'-biphenyl]-4-carbaldehyde) SKU: B4581 CAS Number: 1624970-54-2 Options: • Quantity: 500 mg, 1 g		Br-BT-CHO SKU: B811 CAS Number: 1071224-34-4 Options: • Quantity: 1 g, 5 g, 10 g
	5''-(4'-Formyl-[1,1'-biphenyl]-4-yl)[1,1':—carbaldehyde SKU: B4591 CAS Number: 805246-78-0 Options: • Quantity: 500 mg, 1 g		


Diamines and Dianhydrides

	NTCDA - 1,4,5,8-Naphthalenetetracarboxylic dianhydride SKU: B1191 CAS Number: 81-30-1 Options: Quantity: 100 g, 250 g		2,2'-Bis(trifluoromethyl)benzidine (TFMB) SKU: B2001 CAS Number: 341-58-2 Options: Quantity: 25 g, 50 g, 100 g
	3,3',4,4'-Benzophenonetetracarboxylic dianhydride (BTDA) SKU: B1931 CAS Number: 2421-28-5 Options: Quantity: 100 g, 200 g, 500 g		4,4'-(1,3-Phenylenedioxy)dianiline (TPE-R) SKU: B2011 CAS Number: 2479-46-1 Options: Quantity: 50 g, 100 g, 200 g
	2,2-Bis(4-(4-aminophenoxy)phenyl)propane (BAPP) SKU: B1941 CAS Number: 13080-86-9 Options: Quantity: 50 g, 100 g, 200 g, 500 g		2,2'-Dimethyl[1,1'-biphenyl]-4,4'-diamine SKU: B2021 CAS Number: 84-67-3 Options: Quantity: 50 g, 100 g, 200 g
	1,4-Bis(4-aminophenoxy)benzene (TPE-Q) SKU: B1951 CAS Number: 3491-12-1 Options: Quantity: 25 g, 50 g, 100 g		3,4'-Oxydianiline (3,4'-ODA) SKU: B2031 CAS Number: 2657-87-6 Options: Quantity: 25 g, 50 g, 100 g
	4,4'-Oxydiphthalic anhydride (ODPA) SKU: B1961 CAS Number: 1823-59-2 Options: Quantity: 10 g, 25 g, 50 g		1,2,4,5-cyclohexanetetracarboxylic dianhydride (HPMDA) SKU: B2041 CAS Number: 2754-41-8 Options: Quantity: 5 g, 10 g, 25 g
	2,2-Bis(3-amino-4-hydroxyphenyl)hexafluoropropane (6FAP) SKU: B1971 CAS Number: 83558-87-6 Options: Quantity: 25 g, 50 g, 100 g, 200 g		4,4'-Bis(4-aminophenoxy)biphenyl (BAPB) SKU: B2051 CAS Number: 13080-85-8 Options: Quantity: 25 g, 50 g, 100 g
	9,9-Bis(4-aminophenyl)fluorene (FDA) SKU: B1981 CAS Number: 15499-84-0 Options: Quantity: 50 g, 100 g, 200 g		4,4'-(Hexafluoroisopropylidene)diphthalic anhydride (6FDA) SKU: B2061 CAS Number: 1107-00-2 Options: Quantity: 10 g, 25 g, 50 g
	4,4'-(4,4'-isopropylidenediphenoxy)bis(p-phthalic anhydride) (BPADA) SKU: B2071 CAS Number: 38103-06-9 Options: Quantity: 50 g, 100 g, 200 g		1,4-Bis(4-amino-2-trifluoromethylphenoxy)benzene (6FAPB) SKU: B2071 CAS Number: 94525-05-0 Options: Quantity: 10 g, 25 g, 50 g

	2,2-Bis[4-(4-aminophenoxy)phenyl]hexafluoropropane (4-BDAF) SKU: B2081 CAS Number: 69563-88-8 Options: Quantity: 25 g, 50 g, 100 g		p-6FDAP SKU: B2171 CAS Number: 129197-38-2 Options: Quantity: 5 g, 10 g, 25 g
	2,2'-Bis(trifluoromethyl)-4,4'-diaminodiphenyl ether (6FODA) SKU: B2091 CAS Number: 344-48-9 Options: Quantity: 10 g, 25 g, 50 g, 100 g		AB-TFMB SKU: B2181 CAS Number: 1449757-11-2 Options: Quantity: 2 g, 5 g, 10 g
	1,2,3,4-Butanetetracarboxylic dianhydride (BDA) SKU: B2101 CAS Number: 4534-73-0 Options: Quantity: 25 g, 50 g, 100 g		4-Fluoro-1,2-phenylenediamine SKU: B3431 CAS Number: 367-31-7 Options: Quantity: 5 g, 10 g, 25 g, 50 g
	Bicyclo(2.2.2)oct-7-ene-2,3,5,6-tetracarboxylic acid dianhydride (COeDA) SKU: B2111 CAS Number: 1719-83-1 Options: Quantity: 25 g, 50 g, 100 g		5'-(4,4,5,5-Tetramethyl-1,3,2-dioxaborolan-2-yl)-[1,1':3',1''-terphenyl]-4,4''-diamine SKU: B3661 CAS Number: 2754262-59-2 Options: Quantity: 500 mg, 1 g, 5 g
	4,4'-Bis(3-aminophenoxy)biphenyl (43BAPOBP) SKU: B2121 CAS Number: 105112-76-3 Options: Quantity: 10 g, 25 g, 50 g		4,4'-(Benzo[c][1,2,5]thiadiazole-4,7-diyl)dianiline SKU: B3901 CAS Number: 1203707-77-0 Options: Quantity: 500 mg, 1 g
	2,3,3',4'-Biphenyltetracarboxylic dianhydride (a-BPDA) SKU: B2131 CAS Number: 36978-41-3 Options: Quantity: 10 g, 25 g, 50 g		[1,1':4',1'':4'',1''':4''''-Quaterphenyl]-4,4''''-diamine SKU: B4601 CAS Number: 53693-67-7 Options: Quantity: 500 mg, 1 g
	Dicyclohexyl-3,4,3',4'-tetracarboxylic dianhydride (HBPDA) SKU: B2141 CAS Number: 122640-83-9 Options: Quantity: 10 g, 25 g, 50 g		
	2,2-Bis(3-(3-aminobenzoylamino)-4-hydroxyphenyl)hexafluoropropane (m-6FDAP) SKU: B2161 CAS Number: 220426-92-6 Options: Quantity: 5 g, 10 g, 25 g		


Porous Organic Frameworks

Porous organic frameworks (POFs), including metal-organic frameworks (MOFs) and covalent organic frameworks (COFs), feature tunable pores and high stability, making them valuable in catalysis, gas storage, separation, and sensing. MOFs use ligands functionalized with carboxylic acids, hydroxyl groups, nitrogen-containing heterocycles, and more to control their structure and properties. COFs, by contrast, are constructed from boronic acids, aldehyde-amines, triazines, and other π -conjugated building blocks to form robust, metal-free 2D or 3D frameworks.



- MOF ligands
- COF ligands


Popular Products



1,4-Di(4-pyridyl)benzene

Rigid, linear COF ligand and pillar molecule


- DPyBz
- 4-pyridine end units



H2BTDB

Rigid, linear ligand for photoluminescent properties

- Benzothiadiazole core



4,4',4'',4'''-(Pyrene-1,3,6,8-tetrayl)tetraaniline

Amine-functionalized bridging ligand

- Pyrene core
- Four aniline side arms

To help you choose the right material, full technical specifications are available on our website.



Ossila

B4561-500ml

















5'-Bis(4-aminophenyl)-2,2',1-terphenyl]-4,4''












purity (>97%)
MW 442.55 g/mol

Solpro Business Park

Harmful
Irritation
If swallowed
unwell
EYES: do not touch
Contact with skin may cause
POISONOUS
Persistent

MOF Ligands

















	2,2'-Bipyridine-4,4'-dicarboxylic acid SKU: B1391 CAS Number: 6813-38-3 Options: Quantity: 5 g, 10 g, 25 g		Diethyl 4,4'-(ethyne-1,2-diyl)dibenzoate SKU: B3051 CAS Number: 83536-13-4 Options: Quantity: 1 g, 5 g
	2,2'-Bipyridine-5,5'-dicarboxylic acid SKU: B1411 CAS Number: 1802-30-8 Options: Quantity: 5 g, 10 g, 25 g		1,3,5-Tris(4'-carboxy[1,1'-biphenyl]-4-yl)-benzene SKU: B3061 CAS Number: 911818-75-2 Options: Quantity: 1 g, 5 g
	Bicyclo(2.2.2)oct-7-ene-2,3,5,6-tetracarboxylic acid dianhydride (COeDA) SKU: B2111 CAS Number: 1719-83-1 Options: Quantity: 25 g, 50 g, 100 g		4,4',4'',4'''-(9H-Carbazole-1,3,6,8-tetrayl)-tetrabenzoic acid SKU: B3071 CAS Number: 2097131-99-0 Options: Quantity: 500 mg, 1 g, 5 g
	[1,1':4,1'':4'':1''']-Quaterphenyl]-4,4''-dicarboxylic acid SKU: B3001 CAS Number: 143613-17-6 Options: Quantity: 500 mg, 1 g, 5 g		[1,1':4,1'':4'':1''']-Terphenyl]-2',3,3'',5,5',5'''-hexacarboxylic acid SKU: B3081 CAS Number: 1542274-12-3 Options: Quantity: 500 mg, 1 g, 5 g
	4,4'-(Anthracene-9,10-diylbis(ethyne-2,1-diyl))dibenzoic acid SKU: B3011 CAS Number: 1562777-29-0 Options: Quantity: 500 mg, 1 g, 5 g		Tetra(4-(4-pyridyl)phenyl)methane SKU: B3091 CAS Number: 1319736-15-6 Options: Quantity: 500 mg, 1 g, 5 g
	3,3'',5,5'''-tetrakis(4-carboxyphenyl)-p-terphenyl SKU: B3021 CAS Number: 1816997-25-7 Options: Quantity: 500 mg, 1 g, 5 g		1,1':4,1'':4'':1''']-Quaterphenyl]-3,3'',5,5'''-tetracarboxylic acid SKU: B3101 CAS Number: 921619-91-2 Options: Quantity: 500 mg, 1 g, 5 g
	2,4,6-Tris-(p-carboxyphenyl)pyridine SKU: B3031 CAS Number: 107063-53-6 Options: Quantity: 500 mg, 1 g, 5 g		1,4-Di(4-pyridyl)benzene SKU: B3111 CAS Number: 113682-56-7 Options: Quantity: 1 g, 5 g
	2'-amino-1,1':4,1''-terphenyl-4,4''-dicarboxylic acid SKU: B3041 CAS Number: 1312703-28-8 Options: Quantity: 500 mg, 1 g, 5 g		[p-Terphenyl]-4,4''-dicarboxylic acid SKU: B3121 CAS Number: 13653-84-4 Options: Quantity: 1 g, 5 g, 10 g

















	2',5'-Dimethoxy[1,1':4',1''-terphenyl]-4,4''-dicarboxylic acid SKU: B3131 CAS Number: 1392416-19-1 Options: • Quantity: 500 mg, 1 g, 5 g		Tris(4-ethynylphenyl)amine SKU: B3531 CAS Number: 189178-09-4 Options: • Quantity: 500 mg, 1 g, 5 g
	4'-(pyridin-4-yl)-[1,1'-biphenyl]-4-carboxylic acid SKU: B3141 CAS Number: 1393711-96-0 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		[3,3-Bipyridine]-6,6-dicarboxaldehyde SKU: B3541 CAS Number: 1264748-06-2 Options: • Quantity: 500 mg, 1 g, 5 g
	2-Amino-[1,1'-biphenyl]-4,4''-dicarboxylic acid SKU: B3151 CAS Number: 1240557-01-0 Options: • Quantity: 1 g, 5 g, 10 g		Benzene-1,3,5-triyltriboronic acid SKU: B3551 CAS Number: 89641-21-4 Options: • Quantity: 500 mg, 1 g, 5 g
	4,4'-(Ethyne-1,2-diyl)dibenzoic acid SKU: B3161 CAS Number: 16819-43-5 Options: • Quantity: 250 mg, 500 mg, 1 g		
	5,5'-(Pyridine-3,5-diyl)diisophthalic acid SKU: B3171 CAS Number: 1433029-60-7 Options: • Quantity: 250 mg, 500 mg, 1 g		
	4,4'-(Benzo[c][1,2,5]thiadiazole-4,7-diyl)di-ibenzonic acid SKU: B3181 CAS Number: 1581774-76-6 Options: • Quantity: 1 g, 5 g		
	2,6-Di(pyridin-4-yl)naphthalene SKU: B3191 CAS Number: 950520-39-5 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		
	4',4''',4''''-nitritotris([(1,1'-biphenyl)-3,5-di-carboxylic acid]) SKU: B3201 CAS Number: 1347748-59-7 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		

















COF Ligands

	3,3',4,4'-Benzophenonetetracarboxylic dianhydride (BTDA) SKU: B1931 CAS Number: 2421-28-5 Options: Quantity: 100 g, 200 g, 500 g		4,4'-(1,3-Phenylenedioxy)dianiline (TPE-R) SKU: B2011 CAS Number: 2479-46-1 Options: Quantity: 50 g, 100 g, 200 g
	2,2-Bis[4-(4-aminophenoxy)phenyl]propane (BAPP) SKU: B1941 CAS Number: 13080-86-9 Options: Quantity: 50 g, 100 g, 200 g, 500 g		2,2'-Dimethyl[1,1'-biphenyl]-4,4'-diamine SKU: B2021 CAS Number: 84-67-3 Options: Quantity: 50 g, 100 g, 200 g
	1,4-Bis(4-aminophenoxy)benzene (TPE-Q) SKU: B1951 CAS Number: 3491-12-1 Options: Quantity: 25 g, 50 g, 100 g		3,4'-Oxydianiline (3,4'-ODA) SKU: B2031 CAS Number: 2657-87-6 Options: Quantity: 25 g, 50 g, 100 g
	4,4'-Oxydiphthalic anhydride (ODPA) SKU: B1961 CAS Number: 1823-59-2 Options: Quantity: 10 g, 25 g, 50 g		1,2,4,5-cyclohexanetetracarboxylic dianhydride (HPMDA) SKU: B2041 CAS Number: 2754-41-8 Options: Quantity: 5 g, 10 g, 25 g
	2,2-Bis(4-amino-4-hydroxyphenyl)hexafluoropropane (6FAP) SKU: B1971 CAS Number: 83558-87-6 Options: Quantity: 25 g, 50 g, 100 g, 200 g		4,4'-Bis(4-aminophenoxy)biphenyl (BAPB) SKU: B2051 CAS Number: 13080-85-8 Options: Quantity: 25 g, 50 g, 100 g
	9,9-Bis(4-aminophenyl)fluorene (FDA) SKU: B1981 CAS Number: 15499-84-0 Options: Quantity: 50 g, 100 g, 200 g		4,4'-(Hexafluoroisopropylidene)diphthalic anhydride (6FDA) SKU: B2061 CAS Number: 1107-00-2 Options: Quantity: 10 g, 25 g, 50 g
	4,4'-(4,4'-isopropylidenediphenoxy)bis(phthalic anhydride) (BPADA) SKU: B1991 CAS Number: 38103-06-9 Options: Quantity: 50 g, 100 g, 200 g		1,4-Bis(4-amino-2-trifluoromethylphenoxy)benzene (6FAPB) SKU: B2071 CAS Number: 94525-05-0 Options: Quantity: 10 g, 25 g, 50 g
	2,2'-Bis(trifluoromethyl)benzidine (TFMB) SKU: B2001 CAS Number: 341-58-2 Options: Quantity: 25 g, 50 g, 100 g		2,2-Bis[4-(4-aminophenoxy)phenyl]hexafluoroisopropane (4-BDAF) SKU: B2081 CAS Number: 69563-88-8 Options: Quantity: 25 g, 50 g, 100 g















	2,2'-Bis(trifluoromethyl)-4,4'-diaminodiphenyl ether (6FDA) SKU: B2091 CAS Number: 344-48-9 Options: Quantity: 10 g, 25 g, 50 g, 100 g		[1,1'-Biphenyl]-3,4',5-tricarbaldehyde SKU: B3491 CAS Number: 187281-19-2 Options: Quantity: 500 mg, 1 g, 5 g
	1,2,3,4-Butanetetracarboxylic dianhydride (BDA) SKU: B2101 CAS Number: 4534-73-0 Options: Quantity: 25 g, 50 g, 100 g		2,5-Diethoxyterephthalohydrazide SKU: B3501 CAS Number: 1136292-71-1 Options: Quantity: 500 mg, 1 g, 5 g
	Bicyclo(2.2.2)oct-7-ene-2,3,5,6-tetracarboxylic acid dianhydride (CoEDA) SKU: B2111 CAS Number: 1719-83-1 Options: Quantity: 25 g, 50 g, 100 g		5-(4,4,5,5-Tetramethyl-1,3,2-dioxaborolan-2-yl)isophthalaldehyde SKU: B3511 CAS Number: 945865-80-5 Options: Quantity: 500 mg, 1 g, 5 g
	4,4'-Bis(3-aminophenoxy)biphenyl (43BAPOBP) SKU: B2121 CAS Number: 105112-76-3 Options: Quantity: 10 g, 25 g, 50 g		Hexa(4-formylphenyl)benzene SKU: B3521 CAS Number: 1862220-96-9 Options: Quantity: 500 mg, 1 g, 5 g
	2,2-Bis(3-(3-aminobenzoylamino)-4-hydroxyphenyl)hexafluoropropane (m-6FD-AP) SKU: B2161 CAS Number: 220426-92-6 Options: Quantity: 5 g, 10 g, 25 g		Tris(4-ethynylphenyl)amine SKU: B3531 CAS Number: 189178-09-4 Options: Quantity: 500 mg, 1 g, 5 g
	p-6FDAP SKU: B2171 CAS Number: 129197-38-2 Options: Quantity: 5 g, 10 g, 25 g		[3,3'-Bipyridine]-6,6'-dicarboxaldehyde SKU: B3541 CAS Number: 1264748-06-2 Options: Quantity: 500 mg, 1 g, 5 g
	AB-TFMB SKU: B2181 CAS Number: 1449757-11-2 Options: Quantity: 2 g, 5 g, 10 g		Benzene-1,3,5-triyltriboronic acid SKU: B3551 CAS Number: 89641-21-4 Options: Quantity: 500 mg, 1 g, 5 g
	4,4',4'',4'''-(Ethene-1,1,2,2-tetrayl)tetrabenzaldehyde SKU: B3481 CAS Number: 2170451-48-4 Options: Quantity: 500 mg, 1 g, 5 g		2,2'-(5'-(4-(Cyanomethyl)phenyl)-[1,1':3'-'',1''-terphenyl]-4,4''-diyl)diacetoneitrile SKU: B3561 CAS Number: 2311888-44-3 Options: Quantity: 500 mg, 1 g, 5 g

	4,4'-(Ethyne-1,2-diyl)dibenzaldehyde SKU: B3571 CAS Number: 84907-55-1 Options: Quantity: 500 mg, 1 g, 5 g		5'-(4,4,5,5-Tetramethyl-1,3,2-dioxaborol-2-yl)-4,4'-dibenzaldehyde SKU: B3651 CAS Number: 1033752-94-1 Options: Quantity: 500 mg, 1 g, 5 g
	4,4'-(2,2-Diphenylethene-1,1-diyl)dibenzaldehyde SKU: B3581 CAS Number: 1601465-06-8 Options: Quantity: 500 mg, 1 g, 5 g		5'-(4,4,5,5-Tetramethyl-1,3,2-dioxaborolan-2-yl)-[1,1':3',1''-terphenyl]-4,4''-diamine SKU: B3661 CAS Number: 2754262-59-2 Options: Quantity: 500 mg, 1 g, 5 g
	1,3,5-Tris(4-formylphenyl)benzene SKU: B3591 CAS Number: 118688-53-2 Options: Quantity: 1 g, 5 g, 10 g		Tris(4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl)amine SKU: B3681 CAS Number: 267221-90-9 Options: Quantity: 1 g, 5 g, 10 g
	3,3'-Dihydroxy-[1,1'-biphenyl]-4,4'-dicarbaldehyde SKU: B3601 CAS Number: 14969-32-5 Options: Quantity: 500 mg, 1 g, 5 g		Tetrakis(4-formylphenyl)silane SKU: B3691 CAS Number: 1055999-34-2 Options: Quantity: 250 mg, 500 mg, 1 g, 5 g
	5,5'-(Ethyne-1,2-diyl)diisophthalaldehyde SKU: B3611 CAS Number: 2237965-02-3 Options: Quantity: 500 mg, 1 g, 5 g		[1,1'-Biphenyl]-3,3',5,5'-tetracarbaldehyde SKU: B3701 CAS Number: 150443-85-9 Options: Quantity: 250 mg, 500 mg, 1 g, 5 g
	1,3,5-Tris(4-formyl-3-hydroxyphenyl)benzene SKU: B3621 CAS Number: 2244985-34-8 Options: Quantity: 500 mg, 1 g, 5 g		4,4',4'',4'''-(Pyrene-1,3,6,8-tetrayl)tetrakis(ethyne-2,1-diyl)tetraaniline SKU: B3711 CAS Number: 1404196-75-3 Options: Quantity: 250 mg, 500 mg, 1 g, 5 g
	2,5-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)terephthalaldehyde SKU: B3631 CAS Number: 2305615-53-4 Options: Quantity: 500 mg, 1 g, 5 g		((1,3,5-Triazine-2,4,6-triyl)tris(benzene-4,1-diyl))triboronic acid SKU: B3731 CAS Number: 910231-21-9 Options: Quantity: 250 mg, 500 mg, 1 g, 5 g
	4,4'-(4-(4,4,5,5-Tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl)-4,4'-dibenzaldehyde SKU: B3641 CAS Number: 2009169-65-5 Options: Quantity: 500 mg, 1 g, 5 g		3,5-Diformylphenylboronic acid SKU: B3751 CAS Number: 480424-62-2 Options: Quantity: 500 mg, 1 g, 5 g

	2,3-Bis(4-bromophenyl)fumaronitrile SKU: B3761 CAS Number: 82193-93-9 Options: • Quantity: 500 mg, 1 g, 5 g		4,4''-Diethynyl-1,1'-biphenyl SKU: B3841 CAS Number: 38215-38-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	1,2,4,5-Tetrakis(4-formylphenyl)benzene SKU: B3771 CAS Number: 883835-33-4 Options: • Quantity: 1 g, 5 g		9,10-Bis(4-formylphenyl)anthracene SKU: B3851 CAS Number: 324750-99-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	2,4,6-Tris(4-ethynylphenyl)-1,3,5-triazine SKU: B3781 CAS Number: 425629-22-7 Options: • Quantity: 500 mg, 1 g, 5 g		1-Ethynylpyrene SKU: B3861 CAS Number: 34993-56-1 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	1,3,5-Tris(4-ethynylphenyl)benzene SKU: B3791 CAS Number: 71866-86-9 Options: • Quantity: 1 g, 5 g		Benzene-1,3,5-tricarbohydrazide SKU: B3871 CAS Number: 36997-31-6 Options: • Quantity: 500 mg
	4',4''',4''''-Nitrilotris([1,1'-biphenyl]-4-carbaldehyde) SKU: B3801 CAS Number: 872689-79-7 Options: • Quantity: 500 mg, 1 g, 5 g		4-(1,2,2-Triphenylvinyl)benzaldehyde SKU: B3881 CAS Number: 1289218-74-1 Options: • Quantity: 1 g, 5 g
	5-(4-(Diphenylamino)phenyl)thiophene-2-carbaldehyde SKU: B3811 CAS Number: 291279-14-6 Options: • Quantity: 1 g, 5 g		4,4'-(Benzo[c][1,2,5]thiadiazole-4,7-diyl)d-ibenzaldehyde SKU: B3891 CAS Number: 914651-17-5 Options: • Quantity: 500 mg, 1 g
	4,4''-p-Terphenyldicarboxaldehyde SKU: B3821 CAS Number: 62940-38-9 Options: • Quantity: 1 g, 5 g		4,4'-(Benzo[c][1,2,5]thiadiazole-4,7-diyl)d-aniline SKU: B3901 CAS Number: 1203707-77-0 Options: • Quantity: 500 mg, 1 g
	4,4',4''-(Benzene-1,3,5-triyltris(ethyne-2,1-diyl))tribenzaldehyde SKU: B3831 CAS Number: 2289758-98-9 Options: • Quantity: 500 mg, 1 g, 5 g		5'-(4-Formylphenyl)-2',4',6'-trimethyl-[1,1':3',1''-terphenyl]-4,4''-dicarbaldehyde SKU: B3911 CAS Number: 2195343-70-3 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g

	1,3,5-Tris(3-fluoro-4-formylphenyl)benzene SKU: B3921 CAS Number: 2088007-04-7 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		[1,1':4',1''-Terphenyl]-4,4''-diacetonitrile SKU: B4001 CAS Number: 1000559-50-1 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	1,3,5-Tris[(4-phenylboronic acid)]benzene SKU: B3931 CAS Number: 900795-73-5 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		4,4',4'',4'''-(Pyrene-1,3,6,8-tetrayl)tetraaniline SKU: B4011 CAS Number: 1610471-69-6 Options: • Quantity: 250 mg, 500 mg, 1 g
	1,3,6,8-Tetrakis(4-formylphenyl)pyrene SKU: B3941 CAS Number: 1415238-25-3 Options: • Quantity: 250 mg, 500 mg, 1 g		1,2,4,5-Tetra(4-aminophenyl)ethynylene SKU: B4021 CAS Number: 2375652-88-1 Options: • Quantity: 250 mg, 500 mg, 1 g
	2',5'-Dimethoxy-[1,1':4',1''-terphenyl]-4,4''-dicarbaldehyde SKU: B3951 CAS Number: 111759-27-4 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		2,4,6-Tris(4-hydrazinylphenyl)-1,3,5-triazine SKU: B4031 CAS Number: 804548-56-9 Options: • Quantity: 100 mg, 250 mg, 500 mg
	4,4',4''-(pyridine-2,4,6-triyl)tribenzaldehyde SKU: B3961 CAS Number: 2230887-26-8 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		2,5-Bis(octyloxy)terephthalaldehyde SKU: B4041 CAS Number: 123440-34-6 Options: • Quantity: 500 mg, 1 g
	1,3,6,8-Tetraethynyl-9H-carbazole SKU: B3971 CAS Number: 2415951-75-4 Options: • Quantity: 250 mg		1,1':3',1''-Terphenyl-4,4''-dicarbaldehyde SKU: B4051 CAS Number: 171820-02-3 Options: • Quantity: 500 mg, 1 g
	2,4,6-Tris[4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl]-1,3,5-triazine SKU: B3981 CAS Number: 1447947-87-6 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		2,5-Dimethoxyterephthalohydrazide SKU: B4071 CAS Number: 114503-42-3 Options: • Quantity: 500 mg, 1 g
	4,4'-Biphenyldiacetonitrile SKU: B3991 CAS Number: 7255-83-6 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		[1,1':4',1''-Terphenyl]-3,3'',5,5''-tetracarbaldehyde SKU: B4081 CAS Number: 2380275-62-5 Options: • Quantity: 250 mg, 500 mg, 1 g

	4',4'''',4'''''-Nitrilotris([1,1'-biphenyl]-3,5-dicarbaldehyde) SKU: B4091 CAS Number: 2030105-85-0 Options: • Quantity: 250 mg, 500 mg, 1 g		2,5-Bis(allyloxy)terephthalohydrazide SKU: B4171 CAS Number: 2227151-69-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	2,5-Dipropoxyterephthalohydrazide SKU: B4101 CAS Number: 2143456-42-0 Options: • Quantity: 500 mg, 1 g		Tetrakis(4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl)methane SKU: B4181 CAS Number: 875772-13-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	2,5-Dibutoxyterephthalohydrazide SKU: B4111 CAS Number: 1136292-72-2 Options: • Quantity: 250 mg, 500 mg, 1 g		4,4',4'',4'''-([1,1'-Biphenyl]-4,4'-diylbis(azanetriyl))tetrabenzaldehyde SKU: B4191 CAS Number: 865448-72-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	1,3,5-Tris(4-hydroxyphenyl)benzene SKU: B4121 CAS Number: 15797-52-1 Options: • Quantity: 250 mg, 500 mg, 1 g		[1,1'-Biphenyl]-2,2'-dicarbonitrile SKU: B4201 CAS Number: 4341-02-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	1,2-Bis(4-ethynylphenyl)-1,2-diphenylethane SKU: B4131 CAS Number: 1240785-42-5 Options: • Quantity: 250 mg, 500 mg, 1 g		(4-(1,2,2-Triphenylvinyl)phenyl)methanol SKU: B4211 CAS Number: 1015082-83-3 Options: • Quantity: 500 mg, 1 g
	2,2'-(1,2-Phenylene)bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolane) SKU: B4141 CAS Number: 269410-07-3 Options: • Quantity: 250 mg, 500 mg, 1 g		2,3,6,7,10,11-Hexamethoxytriphenylene SKU: B4221 CAS Number: 808-57-1 Options: • Quantity: 1 g, 5 g
	9-(4-Ethynylphenyl)carbazole SKU: B4151 CAS Number: 262861-81-4 Options: • Quantity: 250 mg, 500 mg, 1 g		1,1,2,2-Tetrakis(4-ethynylphenyl)ethene SKU: B4511 CAS Number: 4863-90-5 Options: • Quantity: 250 mg, 500 mg, 1 g
	[1,1'3',1''-Terphenyl]-4,4'',5'-tricarbaldehyde SKU: B4161 CAS Number: 1977585-07-1 Options: • Quantity: 250 mg, 500 mg, 1 g		1,3,5-Tris(4-aminophenyl)benzene SKU: B4511 CAS Number: 118727-34-7 Options: • Quantity: 1 g, 5 g, 10 g

	2,9-Bis[p-(formyl)phenyl]-1,10-phenanthroline SKU: B4521 CAS Number: 120085-99-6 Options: Quantity: 500 mg, 1 g		[1,1':4',1'':4'',1''':4''']-Quaterphenyl]-4,4''-diamine SKU: B4601 CAS Number: 53693-67-7 Options: Quantity: 500 mg, 1 g
	5,5',5''-(Benzene-1,3,5-triyl)tripicolinaldehyde SKU: B4531 CAS Number: 1482413-54-6 Options: Quantity: 250 mg, 500 mg, 1 g		4',4''',4''''-(1,3,5-Triazine-2,4,6-triyl)tris([1,1'-biphenyl]-4-carbaldehyde) SKU: B4611 CAS Number: 1221509-80-3 Options: Quantity: 250 mg, 500 mg, 1 g
	4,4',4''-(Benzene-1,3,5-triyltris(ethyne-2,1-diyl))trianiline SKU: B4541 CAS Number: 326002-91-9 Options: Quantity: 500 mg, 1 g		1,3,5-Tris[4-amino(1,1-biphenyl-4-yl)]benzene SKU: B4621 CAS Number: 1400987-00-9 Options: Quantity: 500 mg, 1 g
	Tris[4-(5-formyl-2-thienyl)phenyl]amine SKU: B4551 CAS Number: 883236-47-3 Options: Quantity: 500 mg, 1 g		1,3,6,8-Tetrakis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)pyrene SKU: B4631 CAS Number: 1398053-00-3 Options: Quantity: 500 mg, 1 g
	4',5'-Bis(4-aminophenyl)-[1,1':2',1'':terphenyl]-4,4''-diamine SKU: B4561 CAS Number: 2458125-05-6 Options: Quantity: 500 mg, 1 g		4',4''',4''''-(1,3,5-triazine-2,4,6-triyl)tris([1,1'-biphenyl]-4-amine) SKU: B4641 CAS Number: 2130745-76-3 Options: Quantity: 500 mg, 1 g
	[1,1':3',1'':Terphenyl]-3,3'',5,5''-tetracarbaldehyde SKU: B4571 CAS Number: 201734-76-1 Options: Quantity: 500 mg, 1 g		1,3,5-Tris(4-cyanophenyl)benzene SKU: B4651 CAS Number: 382137-78-2 Options: Quantity: 500 mg, 1 g
	4',4''',4''''-(Ethene-1,1,2,2-tetrayl)tetrakis([1,1'-biphenyl]-4-carbaldehyde) SKU: B4581 CAS Number: 1624970-54-2 Options: Quantity: 500 mg, 1 g		
	5''-(4'-Formyl-[1,1'-biphenyl]-4-yl)-[1,1':-carbaldehyde SKU: B4591 CAS Number: 805246-78-0 Options: Quantity: 500 mg, 1 g		

Low Dimensional Materials

Low dimensional materials exhibit reduced dimensions compared to their bulk counterparts. These materials include two-dimensional (2D) systems like graphene and transition metal dichalcogenides (TMDs), and one-dimensional (1D) systems like carbon nanotubes and nanowires. The low dimensionality results from reduced thickness in one or more dimensions. This characteristic unlocks unique properties for the material, enabling applications in energy storage and next-generation electronics.



Popular Products



Black Phosphorus

Available in crystal or powder form


- Versatile
- High Purity
- Low Cost



Graphene Materials

Available as powder, dispersions, and more

- Large Surface Area
- Great Semiconductor




Single-Walled Carbon Nanotubes




Useful properties for versatile applications

- High Purity
- Range of dimensions available

To help you choose the right material, full technical specifications are available on our website.



Carbon Nanotubes (CNTs)

	Multi-Walled Carbon Nanotubes (MWCNT) SKU: M2008 CAS Number: 308068-56-6 Options*: • Purity: 96%, 99% • Quantity: 5 g, 10 g, 25 g, 50 g		OH Functionalized Double-Walled Carbon Nanotubes SKU: M2018 CAS Number: 308068-56-6 Options*: • Purity: 60%, 65% • Quantity: 250 mg, 500 mg, 1 g
	COOH Functionalized Multi-Walled Carbon Nanotubes SKU: M2010 CAS Number: 308068-56-6 Options*: • Purity: 96%, 99.9% • Quantity: 500 mg, 1 g, 5 g, 10 g		Graphene Carbon Nanotubes Composite SKU: M2373 CAS Number: 308068-56-6 Options: • Quantity: 5 g, 10 g
	OH Functionalized Multi-Walled Carbon Nanotubes SKU: M2011 CAS Number: 308068-56-6 Options*: • Purity: 96%, 99.9% • Quantity: 500 mg, 1 g, 5 g, 10 g		Graphitized Multi-Walled Carbon Nanotubes (G-MWCNT) SKU: M2374 CAS Number: 308068-56-6 Options*: • Quantity: 1 g, 5 g, 10 g
	Single-Walled Carbon Nanotubes (SWCNT) SKU: M2013 CAS Number: 308068-56-6 Options*: • Purity: > 65%, > 95%, > 96%, ≥ 95% • Quantity: 250 mg, 500 mg, 1 g		Ni-MWCNTs (Nickel Coated Multi-walled Carbon Nanotubes) SKU: M2412 CAS Number: 308068-56-6 / 7440-02-0 Options: • Quantity: 1 g, 2 g, 5 g • Size: 4 - 16 nm, 18 - 28 nm, 28 - 48 nm
	COOH Functionalized Single-Walled Carbon Nanotubes SKU: M2014 CAS Number: 308068-56-6 Options*: • Purity: > 92% • Quantity: 250 mg, 500 mg, 1 g		Aligned Multi-Walled Carbon Nanotubes SKU: M2438 CAS Number: 308068-56-6 Options: • Length: 9 - 14 μm, 25 - 95 μm • Quantity: 500 mg, 1 g, 2 g, 5 g
	OH Functionalized Single-Walled Carbon Nanotubes SKU: M2015 CAS Number: 308068-56-6 Options*: • Purity: 95% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	Double-Walled Carbon Nanotubes (DWCNT) SKU: M2016 CAS Number: 308068-56-6 Options*: • Purity: 65%, 91% • Quantity: 500 mg, 1 g, 5 g		
	COOH Functionalized Double-Walled Carbon Nanotubes SKU: M2017 CAS Number: 308068-56-6 Options*: • Purity: > 60%, > 65% • Quantity: 250 mg, 500 mg, 1 g		

* Some options have been omitted, please see our website to browse our full range

2D Materials

	Graphene Oxide Powder SKU: M0880 Options*: - Flake size: 4.5 μm Flake Size, ~ 2 μm Flake Size - Weight/volume: 1 g, 2 g, 250 mg, 500 mg		Molybdenum Ditelluride (MoTe₂) Powder and Crystal SKU: M2109 CAS Number: 12058-20-7 Options*: Weight/size: 500 mg, 1 g, Medium ($\geq 25 \text{ mm}^2$)
	Reduced Graphene Oxide Powders SKU: M0921 CAS Number: 1034343-98-0 Options*: - Amount: 1 g, 5 g - Flake size: 1 - 5 μm		Tungsten Disulfide (WS₂) Crystal SKU: M2110 CAS Number: 12138-09-9 Options: Size: Large ($\geq 100 \text{ mm}^2$), Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Black Phosphorus Powder SKU: M2106 CAS Number: 7723-14-0 Options: - Packaging: Ampoule; Purity $\geq 99.999\%$, Vacuum sealed; Purity $\geq 99.995\%$ - Quantity: 250 mg, 500 mg, 1 g		Tungsten Disulfide (WS₂) Powder SKU: M2110 CAS Number: 12138-09-9 Options: Quantity: 500 mg, 1 g
	Black Phosphorus Crystal SKU: M2106 CAS Number: 7723-14-0 Options*: Purity: $\geq 99.995\%$, $\geq 99.999\%$		Tungsten Diselenide (WSe₂) Powder and Crystal SKU: M2111 CAS Number: 12067-46-8 Options*: Weight/size: 500 mg, 1 g, Large ($\geq 100 \text{ mm}^2$), Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Molybdenum Disulfide (MoS₂) Crystal SKU: M2107 CAS Number: 1317-33-5 Options: Crystal size: Large ($20 \times 20 \text{ mm}^2$), Medium ($15 \times 15 \text{ mm}^2$), Small ($10 \times 10 \text{ mm}^2$)		Tungsten Ditelluride (WTe₂) Powder and Crystal SKU: M2112 CAS Number: 12067-76-4 Options*: Weight/size: 500 mg, 1 g
	Molybdenum Disulfide Powder SKU: M2107 CAS Number: 1317-33-5 Options: Quantity: 500 mg, 1 g		Tin(II) Sulfide (SnS) Crystal SKU: M2113 CAS Number: 1314-95-0 Options: Size: Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Molybdenum Diselenide (MoSe₂) Crystal SKU: M2108 CAS Number: 12058-18-3 Options: Size: Large ($\geq 100 \text{ mm}^2$), Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)		Tin(IV) Sulfide (SnS₂) Powder and Crystal SKU: M2114 CAS Number: 1315-01-1 Options*: Weight/size: 500 mg, 1 g, Large ($\geq 100 \text{ mm}^2$), Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Molybdenum Diselenide (MoSe₂) Powder SKU: M2108 CAS Number: 12058-18-3 Options: Quantity: 500 mg, 1 g		Tin(IV) Selenide (SnSe₂) Powder and Crystal SKU: M2115 CAS Number: 20770-09-6 Options*: Weight/size: 500 mg, 1 g, Large ($\geq 100 \text{ mm}^2$), Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)

* Some options have been omitted, please see our website to browse our full range

	Tin(II) Selenide (SnSe) Crystal SKU: M2116 CAS Number: 1315-06-6 Options: • Size: Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)		Rhenium Disulfide (ReS2) Powder and Crystal SKU: M2139 CAS Number: 12038-63-0 Options*: • Weight/size: 500 mg, 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Hexagonal Boron Nitride (h-BN) Crystals SKU: M2133 CAS Number: 10043-11-5 Options: • Package/substrate: Pack of 10 crystals, Pack of 5 crystals		Rhenium Diselenide (ReSe2) Powder and Crystal SKU: M2140 CAS Number: 12038-64-1 Options*: • Weight/size: 500 mg, 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Hexagonal Boron Nitride (h-BN) Powder SKU: M2133 CAS Number: 10043-11-5 Options: • Particle size: 1 – 2 μm , 100 – 200 nm, 30 – 40 μm , 40 – 50 μm , 790 nm, 800 nm • Quantity: 10 g, 25 g, 50 g, 100 g, 250 g, 500 g		Molybdenum Tungsten Disulfide (MoWS2) Powder SKU: M2141 CAS Number: 109657-36-5 Options: • Quantity: 500 mg, 1 g
	Hexagonal Boron Nitride Nanopowders SKU: M2133 CAS Number: 10043-11-5 Options: • Particle size: 100 – 200 nm, 790 nm, 800 nm • Quantity: 10 g, 25 g, 50 g, 100 g, 250 g, 500 g		Molybdenum Tungsten Diselenide (MoWSe2) Powder SKU: M2142 CAS Number: 126414-37-7 Options: • Quantity: 500 mg, 1 g
	Indium(II) Selenide (InSe) Powder and Crystal SKU: M2135 CAS Number: 1312-42-1 Options*: • Weight/size: 500 mg, 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)		Molybdenum Sulfide Selenide (MoSSe) Powder SKU: M2143 CAS Number: 132004-88-7 Options: • Quantity: 500 mg, 1 g
	Indium(III) Selenide (In2Se3) Powder and Crystal SKU: M2136 CAS Number: 12056-07-4 Options*: • Weight/size: 500 mg, 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)		Tantalum Disulfide (TaS2) Powder and Crystal SKU: M2144 CAS Number: 12143-72-5 Options*: • Weight/size: 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Germanium Selenide (GeSe) Crystal SKU: M2137 CAS Number: 12065-10-0 Options: • Title: Default Title		Tantalum Diselenide (TaSe2) Powder and Crystal SKU: M2145 CAS Number: 12039-55-3 Options*: • Weight/size: 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Gallium Selenide (GaSe) Powder and Crystals SKU: M2138 CAS Number: 12024-11-2 Options: • Form: Gallium Selenide Crystal, Gallium Selenide Powder • Size/weight: 1 g, 500 mg, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)		Titanium(IV) Sulfide (TiS2) Powder and Crystal SKU: M2147 CAS Number: 12039-13-3 Options*: • Weight/size: 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)

* Some options have been omitted, please see our website to browse our full range

	Titanium(IV) Selenide (TiSe₂) Powder and Crystal SKU: M2151 CAS Number: 12067-45-7 Options*: • Weight/size: 1 g, Small ($\geq 10 \text{ mm}^2$)		Molybdenum Diselenide Monolayer Film SKU: M2169 CAS Number: 12058-18-3 Options: • Quantity: 1 Piece - SiO ₂ /Si
	Hafnium Disulfide (HfS₂) Powder and Crystal SKU: M2152 CAS Number: 18855-94-2 Options*: • Weight/size: 1 g, Small ($\geq 10 \text{ mm}^2$)		Molybdenum Diselenide Few-Layer Film SKU: M2170 CAS Number: 12058-18-3 Options: • Quantity: 1 Piece - SiO ₂ /Si
	Hafnium Diselenide (HfSe₂) Powder and Crystal SKU: M2153 CAS Number: 12162-21-9 Options*: • Weight/size: 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)		Tungsten Disulfide (WS₂) Monolayer Film SKU: M2171 CAS Number: 12138-09-9 Options: • Quantity: 1 Piece - Sapphire
	Vanadium Diselenide (VSe₂) Powder and Crystal SKU: M2154 CAS Number: 12299-51-3 Options*: • Weight/size: 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)		Tungsten Disulfide (WS₂) Few-Layer Film SKU: M2172 CAS Number: 12138-09-9 Options: • Quantity: 1 Piece - Sapphire
	Graphene Film SKU: M2155 CAS Number: 1034343-98-0 Options: • Film type: Few-Layer Film, Monolayer Film • Quantity: 1, 2 pieces, 4 pieces • Substrate type: Cu Foil (5 cm x 10 cm), PET (1 cm x 1 cm), Quartz (1 cm x 1 cm), Sapphire (1 cm x 1 cm), SiO ₂ /Si (1 cm x 1 cm)		Zirconium Disulfide (ZrS₂) Powder and Crystal SKU: M2201 CAS Number: 12039-15-5 Options*: • Weight/size: 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Hexagonal Boron Nitride (h-BN) Film SKU: M2161 CAS Number: 10043-11-5 Options*: • Package/substrate: PET (1 cm x 1 cm) - 1 Each, PET (1 cm x 1 cm) - 2 Each, Quartz (1 cm x 1 cm) - 2 Each, Quartz (1 cm x 1 cm) - 1 Each, Quartz (1 cm x 1 cm) - 2 Each, Sapphire (1 cm x 1 cm) - 1 Each, Sapphire (1 cm x 1 cm) - 2 Each, SiO ₂ /Si (1 cm x 1 cm) - 1 Each, SiO ₂ /Si (1 cm x 1 cm) - 2 Each		Zirconium Diselenide (ZrSe₂) Powder and Crystal SKU: M2202 CAS Number: 12166-47-1 Options*: • Weight/size: 500 mg, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Molybdenum Disulfide Monolayer Film SKU: M2165 CAS Number: 1317-33-5 Options: • Quantity: 1 Piece - Isolated Triangular Flakes on Sapphire, 1 Piece - Monolayer Film on Sapphire, 1 Piece - Monolayer Film on SiO ₂ /Si		Zirconium Trisulfide (ZrS₃) Powder and Crystal SKU: M2203 CAS Number: 12166-31-3 Options*: • Weight/size: 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)
	Molybdenum Disulfide Few-Layer Film SKU: M2167 CAS Number: 1317-33-5 Options: • Quantity: 1 Piece - Sapphire, 1 Piece - SiO ₂ /Si		Zirconium Triselenide (ZrSe₃) Powders and Crystals SKU: M2204 CAS Number: 12166-53-9 Options*: • Weight/size: 1 g, Medium ($\geq 25 \text{ mm}^2$), Small ($\geq 10 \text{ mm}^2$)

* Some options have been omitted, please see our website to browse our full range

	Bismuth Sulfide (Bi₂S₃) Powders and Crystal SKU: M2205 CAS Number: 1345-07-9 Options*: • Weight/size: 1 g, 5 g, 10 g, Medium, Small		Tungsten Diselenide (WSe₂) Few-layer Film SKU: M2224 CAS Number: 12067-46-8 Options: • Size: 1 * 1 cm • Substrate: Sapphire - Double Side Polished
	Bismuth Selenide (Bi₂Se₃) Powders and Crystal SKU: M2206 CAS Number: 2068-69-8 Options*: • Weight/size: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)		Gallium Sulfide (GaS) Powder and Crystals SKU: M2299 CAS Number: 12024-10-1 Options: • Form: Gallium Sulfide Crystal (α-GaS), Gallium Sulfide Powder • Size/weight: 1 g, 500 mg, Medium (≥25 mm ²), Small (≥10 mm ²)
	Iron Phosphorus Trisulfide (FePS₃) Powder and Crystals SKU: M2207 Options*: • Weight/size: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)		Niobium Disulfide (NbS₂) Powder and Crystal SKU: M2300 CAS Number: 12034-77-4 Options*: • Weight/size: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)
	Iron Phosphorus Triselenide (FePSe₃) Powder and Crystal SKU: M2208 Options*: • Weight/size: 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)		Niobium Diselenide (NbSe₂) Powder and Crystals SKU: M2302 CAS Number: 12034-77-4 Options*: • Weight/size: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)
	Manganese Phosphorus Trisulfide (MnPS₃) Powder and Crystal SKU: M2209 CAS Number: 20642-09-5 Options*: • Weight/size: 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)		Germanium Sulfide (GeS) Powder and Crystals SKU: M2303 CAS Number: 12025-32-0 Options*: • Weight/size: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)
	Manganese Phosphorus Triselenide (MnPSe₃) Powder and Crystal SKU: M2210 Options*: • Weight/size: 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)		Germanium Disulfide (GeS₂) Powder and Crystals SKU: M2305 CAS Number: 12025-34-2 Options*: • Weight/size: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)
	Nickel Phosphorus Trisulfide (NiPS₃) Powder and Crystal SKU: M2211 Options*: • Weight/size: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)		Germanium Diselenide (GeSe₂) Powder and Crystals SKU: M2307 CAS Number: 12065-11-1 Options*: • Weight/size: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)
	Tungsten Diselenide (WSe₂) Monolayer Film SKU: M2223 CAS Number: 12067-46-8 Options: • Size: 1 * 1 cm • Substrate: Sapphire - Double Side Polished		Chromium(III) Iodide (CrI₃) Powder and Crystals SKU: M2323 CAS Number: 13569-75-0 Options*: • Quantity: 500 mg, 1 g, Medium (≥25 mm ²), Small (≥10 mm ²)

* Some options have been omitted, please see our website to browse our full range

	Silver Indium Phosphorus Sulfide (AgInP2—S6) Powder SKU: M2325 Options: • Quantity: 500 mg, 1 g		Graphite Nanopowder SKU: M2393 CAS Number: 7782-42-5 Options*: • Quantity: 5 g, 10 g, 25 g
	Silver Indium Phosphorus Selenide (AgIn—P2Se6) Powder SKU: M2326 Options: • Quantity: 500 mg, 1 g		Titanium Dioxide (TiO2) Nanoparticles SKU: M2394 CAS Number: 1317-70-0 / 1317-80-2 / 13463-67-7 Options*: • Form: Powder, Solution in Water • Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml
	Silver Tantalum Sulfide (AgTaS3) Powder SKU: M2328 Options: • Quantity: 500 mg, 1 g		Titanium Dioxide Anatase Nanoparticles SKU: M2394 CAS Number: 1317-70-0 Options*: • Form: Powder, Solution in Water • Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml
	Graphene Sheet SKU: M2371 CAS Number: 1034343-98-0 Options: • Quantity: 10 Pieces, 5 Pieces • Size: 10 × 10 cm, 5 × 5 cm		Rutile Titanium Dioxide (TiO2) Nanoparticles SKU: M2394 CAS Number: 1317-80-2 Options*: • Quantity/volume: 50 g, 100 g, 250 g
	Boron-doped Graphene Powder SKU: M2372 CAS Number: 1034343-98-0 / 7440-42-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Titanium Dioxide (TiO2) P25 Nanoparticles SKU: M2394 CAS Number: 13463-67-7 Options*: • Quantity/volume: 50 g, 100 g, 250 g
	Graphene Carbon Nanotubes Composite SKU: M2373 CAS Number: 308068-56-6 Options: • Quantity: 5 g, 10 g		Yttrium Oxide (Y₂O₃) Nanopowder SKU: M2397 CAS Number: 1314-36-9 Options: • Quantity: 50 g, 100 g, 250 g
	Graphene Nanoplatelet Powder SKU: M2373 CAS Number: 1034343-98-0 Options: • Quantity: 25 g, 50 g • Size: ~ 1.5 µm, ~ 18 µm, ~ 7 µm		Zinc Oxide Nanoparticles SKU: M2398 CAS Number: 1314-13-2 Options: • Quantity: 25 g, 50 g, 100 g
	Graphite Powder SKU: M2393 CAS Number: 7782-42-5 Options*: • Quantity: 5 g, 10 g, 25 g, 50 g, 100 g, 250 g		Antimony Oxide (Sb₂O₃) Powder SKU: M2399 CAS Number: 1309-64-4 Options: • Particle size: 20 – 30 nm, < 210 nm • Quantity: 50 g, 100 g, 250 g

* Some options have been omitted, please see our website to browse our full range

	Aluminum Oxide (Al₂O₃) Nanopowder SKU: M2400 CAS Number: 1344-28-1 Options: - Quantity: 50 g, 100 g, 250 g, 500 g - Specification: Particle size 20 nm; Purity ≥ 99.99%, Particles size ~ 18 nm; Purity > 99.5%		Lithium Nickel Manganese Cobalt Oxide Powder (NCM523) SKU: M2409 CAS Number: 346417-97-8 / 182442-95-1 Options*: Quantity: 50 g, 100 g, 250 g
	Lithium Manganese Oxide (LMO) Powder SKU: M2402 CAS Number: 12057-17-9 Options: - Particle size (d50): 14.8 μm, 23 ~ 27 μm - Quantity: 100 g, 250 g, 500 g		Graphene Dispersion SKU: M2448 CAS Number: 1034343-98-0 Options: - Quantity: 50 mL, 100 mL, 250 mL - Solvent base: NMP, Water
	Lithium Nickel Manganese Oxide (LNMO) Powder SKU: M2403 CAS Number: 12031-75-3 Options: Quantity: 50 g, 100 g, 250 g		Graphene Oxide Dispersion SKU: M2450 Options: Quantity: 100 mL, 250 mL, 1000 mL
	Lithium Nickel Cobalt Aluminum Oxide (NCA) Powder SKU: M2404 CAS Number: 193214-24-3 Options: Quantity: 10 g, 25 g, 50 g, 100 g		Silver Nanowire Dispersion SKU: M2456 CAS Number: 7440-22-4 Options: Quantity: 25 mL, 50 mL
	Lithium Nickel Manganese Cobalt Oxide (NMC811) Powder SKU: M2405 CAS Number: 179802-95-0 Options: Quantity: 50 g, 100 g		
	Holey Graphene SKU: M2406 CAS Number: 1034343-98-0 Options: Quantity: 100 mg, 250 mg		
	Bismuth Oxide (Bi₂O₃) Nanopowder SKU: M2407 CAS Number: 1304-76-3 Options: Quantity: 50 g, 100 g, 250 g		
	Lithium Iron Phosphate (LiFePO₄) Powder SKU: M2408 CAS Number: 15365-14-7 Options*: Quantity: 100 g, 250 g		

* Some options have been omitted, please see our website to browse our full range

Battery Materials

As the demand for clean and renewable energy sources continues to grow, so does the need for efficient electrical storage systems. Renewable energy must be stored and distributed to match energy demand effectively. Our range of lithium-ion battery materials includes the six primary types—NMC, NCA, NCM, LFP, LMO, and LCO—along with the newer LNMO. Additionally, our anode active materials, such as lithium titanate, carbon black, and graphite, can be coated onto copper foils for enhanced performance.



Popular Products



LNMO Powder

Cathode active,
high voltage and high stability

- Excellent thermal stability
- High voltage
- Good rate capability



LiFePO₄ Powder

Cathode active, long lifetime
and highly safe

- Prohibiting oxygen release
- Long Battery Life
- Stable phosphate group



LTO Powder

Anode active,
long lifetime and highly safe

- Long battery life
- High safety
- Cost Effective Value

To help you choose the right material, full technical specifications are available on our website.




100% Cobalt Oxide
powder

Solpro Business Park Winder...

May damage fertility or the embryo.
Wear protective gloves/goggles.
If concerned, get medical advice.



Cathode Active Materials

	Pyrene-4,5,9,10-tetraone SKU: B3721 CAS Number: 14727-71-0 Options: Quantity: 1 g, 250 mg, 5 g, 500 mg		Lithium Nickel Cobalt Aluminum Oxide (NCA) Powder SKU: M2404 CAS Number: 193214-24-3 Options: Quantity: 10 g, 25 g, 50 g, 100 g
	Titanium Dioxide (TiO2) Nanoparticles SKU: M2394 CAS Number: 1317-70-0 / 1317-80-2 / 13463-67-7 Options*: - Form: Powder, Solution in Water - Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml		Lithium Nickel Manganese Cobalt Oxide (NMC811) Powder SKU: M2405 CAS Number: 179802-95-0 Options: Quantity: 50 g, 100 g
	Titanium Dioxide Anatase Nanoparticles SKU: M2394 CAS Number: 1317-70-0 Options*: - Form: Powder, Solution in Water - Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml		Lithium Iron Phosphate (LiFePO4) Powder SKU: M2408 CAS Number: 15365-14-7 Options*: Quantity: 100 g, 250 g
	Rutile Titanium Dioxide (TiO2) Nanoparticles SKU: M2394 CAS Number: 1317-80-2 Options*: Quantity/volume: 50 g, 100 g, 250 g		Lithium Nickel Manganese Cobalt Oxide Powder (NCM523) SKU: M2409 CAS Number: 346417-97-8 / 182442-95-1 Options*: Quantity: 50 g, 100 g, 250 g
	Titanium Dioxide (TiO2) P25 Nanoparticles SKU: M2394 CAS Number: 13463-67-7 Options*: Quantity/volume: 50 g, 100 g, 250 g		Nickel (II) Oxide Powder SKU: M2441 CAS Number: 1313-99-1 Options: - Particle size: 8 - 18 nm - Quantity: 10 g, 25 g
	Lithium Cobalt Oxide (LiCoO2) Powder SKU: M2401 CAS Number: 12190-79-3 Options: - Particle size (d50): 13 µm - Quantity: 50 g, 100 g		Barium Titanate Powder SKU: M2443 CAS Number: 12047-27-7 Options: - Average particle size: 90 nm - Quantity: 25 g, 50 g, 100 g, 250 g
	Lithium Manganese Oxide (LMO) Powder SKU: M2402 CAS Number: 12057-17-9 Options: - Particle size (d50): 14.8 µm, 23 - 27 µm - Quantity: 100 g, 250 g, 500 g		Manganese (IV) Oxide Powder SKU: M2490 CAS Number: 1313-13-9 Options: Quantity: 100 g, 250 g
	Lithium Nickel Manganese Oxide (LNMO) Powder SKU: M2403 CAS Number: 12031-75-3 Options: Quantity: 50 g, 100 g, 250 g		Lithium Manganese Iron Phosphate (LMFP) Powder SKU: M2491 CAS Number: 1354700-46-1 Options: Quantity: 50 g, 100 g, 250 g, 500 g

* Some options have been omitted, please see our website to browse our full range

	Lithium Nickel Manganese Cobalt Oxide Powder (NMC622) SKU: M2546 CAS Number: 346417-97-8 / 182442-95-1 Options: • Quantity: 50 g, 100 g, 250 g
	Lithium Nickel Manganese Cobalt Oxide Powder (NCM111) SKU: M2547 CAS Number: 346417-97-8 / 182442-95-1 Options: • Quantity: 50 g, 100 g, 250 g
	Prussian Blue Analogue Powder SKU: M2551 CAS Number: 14038-43-8 Options: • Quantity: 10 g, 25 g
	Prussian White Powder SKU: M2552 CAS Number: 14038-43-8 Options: • Quantity: 10 g, 25 g, 50 g

Anode Active Materials

	Reduced Graphene Oxide Powders SKU: M0921 CAS Number: 1034343-98-0 Options*: - Amount: 1 g, 5 g - Flake size: 1 - 5 µm		Black Phosphorus Powder SKU: M2106 CAS Number: 7723-14-0 Options: - Packaging: Ampoule; Purity ≥99.999%, Vacuum sealed; Purity ≥99.995% - Quantity: 250 mg, 500 mg, 1 g
	Multi-Walled Carbon Nanotubes (MWCNT) SKU: M2008 CAS Number: 308068-56-6 Options*: - Purity: 96%, 99% - Quantity: 5 g, 10 g, 25 g, 50 g		Black Phosphorus Crystal SKU: M2106 CAS Number: 7723-14-0 Options*: Purity: ≥ 99.995%, ≥ 99.999%
	COOH Functionalized Multi-Walled Carbon Nanotubes SKU: M2010 CAS Number: 308068-56-6 Options*: - Purity: 96%, 99.9% - Quantity: 500 mg, 1 g, 5 g, 10 g		Molybdenum Diselenide (MoSe₂) Powder SKU: M2108 CAS Number: 12058-18-3 Options: Quantity: 500 mg, 1 g
	OH Functionalized Multi-Walled Carbon Nanotubes SKU: M2011 CAS Number: 308068-56-6 Options*: - Purity: 96%, 99.9% - Quantity: 500 mg, 1 g, 5 g, 10 g		Tin(IV) Sulfide (SnS₂) Powder and Crystal SKU: M2114 CAS Number: 1315-01-1 Options*: Weight/size: 500 mg, 1 g, Large (≥100 mm²), Medium (≥25 mm²), Small (≥10 mm²)
	Single-Walled Carbon Nanotubes (SWCNT) SKU: M2013 CAS Number: 308068-56-6 Options*: - Purity: > 65%, > 95%, > 96%, ≥ 95% - Quantity: 250 mg, 500 mg, 1 g		Molybdenum Sulfide Selenide (MoSSe) Powder SKU: M2143 CAS Number: 132004-88-7 Options: Quantity: 500 mg, 1 g
	COOH Functionalized Single-Walled Carbon Nanotubes SKU: M2014 CAS Number: 308068-56-6 Options*: - Purity: > 92% - Quantity: 250 mg, 500 mg, 1 g		Nickel Phosphorus Trisulfide (NiPS₃) Powder and Crystal SKU: M2211 Options*: Weight/size: 500 mg, 1 g, Medium (≥25 mm²), Small (≥10 mm²)
	OH Functionalized Single-Walled Carbon Nanotubes SKU: M2015 CAS Number: 308068-56-6 Options*: - Purity: 95% - Quantity: 100 mg, 250 mg, 500 mg, 1 g		Niobium Diselenide (NbSe₂) Powder and Crystals SKU: M2302 CAS Number: 12034-77-4 Options*: Weight/size: 500 mg, 1 g, Medium (≥25 mm²), Small (≥10 mm²)
	COOH Functionalized Double-Walled Carbon Nanotubes SKU: M2017 CAS Number: 308068-56-6 Options*: - Purity: > 60%, > 65% - Quantity: 250 mg, 500 mg, 1 g		Germanium Disulfide (GeS₂) Powder and Crystals SKU: M2305 CAS Number: 12025-34-2 Options*: Weight/size: 500 mg, 1 g, Medium (≥25 mm²), Small (≥10 mm²)

* Some options have been omitted, please see our website to browse our full range

	Graphene Sheet SKU: M2371 CAS Number: 1034343-98-0 Options: - Quantity: 10 Pieces, 5 Pieces - Size: 10×10 cm, 5×5 cm		Rutile Titanium Dioxide (TiO₂) Nanoparticles SKU: M2394 CAS Number: 1317-80-2 Options*: Quantity/volume: 50 g, 100 g, 250 g
	Boron-doped Graphene Powder SKU: M2372 CAS Number: 1034343-98-0 / 7440-42-8 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g		Titanium Dioxide (TiO₂) P25 Nanoparticles SKU: M2394 CAS Number: 13463-67-7 Options*: Quantity/volume: 50 g, 100 g, 250 g
	Graphitized Multi-Walled Carbon Nanotubes (G-MWCNT) SKU: M2374 CAS Number: 308068-56-6 Options*: Quantity: 1 g, 5 g, 10 g		Carbon Black Powder SKU: M2396 CAS Number: 1333-86-4 Options: - Average particle size: 30 nm, 4 μm - Quantity: 10 g, 25 g, 50 g, 100 g
	Graphene Nanoplatelet Powder SKU: M2392 CAS Number: 1034343-98-0 Options: - Quantity: 25 g, 50 g - Size: $\sim 1.5 \mu$m, $\sim 18 \mu$m, $\sim 7 \mu$m		Carbon Black Nanopowder SKU: M2396 CAS Number: 1333-86-4 Options: - Average particle size: 30 nm - Quantity: 25 g, 50 g, 100 g
	Graphite Powder SKU: M2393 CAS Number: 7782-42-5 Options*: Quantity: 5 g, 10 g, 25 g, 50 g, 100 g, 250 g		Antimony Oxide (Sb₂O₃) Powder SKU: M2399 CAS Number: 1309-64-4 Options: - Particle size: 20 $\sim$ 30 nm, < 210 nm - Quantity: 50 g, 100 g, 250 g
	Graphite Nanopowder SKU: M2393 CAS Number: 7782-42-5 Options*: Quantity: 5 g, 10 g, 25 g		Bismuth Oxide (Bi₂O₃) Nanopowder SKU: M2407 CAS Number: 1304-76-3 Options: Quantity: 50 g, 100 g, 250 g
	Titanium Dioxide (TiO₂) Nanoparticles SKU: M2394 CAS Number: 1317-70-0 / 1317-80-2 / 13463-67-7 Options*: - Form: Powder, Solution in Water - Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml		Ni-MWCNTs (Nickel Coated Multi-walled Carbon Nanotubes) SKU: M2412 CAS Number: 308068-56-6 / 7440-02-0 Options: - Quantity: 1 g, 2 g, 5 g - Size: 4 - 16 nm, 18 - 28 nm, 28 - 48 nm
	Titanium Dioxide Anatase Nanoparticles SKU: M2394 CAS Number: 1317-70-0 Options*: - Form: Powder, Solution in Water - Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml		Lithium Titanium Oxide Powder SKU: M2428 CAS Number: 12031-95-7 Options: Quantity: 25 g, 50 g, 100 g, 250 g

* Some options have been omitted, please see our website to browse our full range

	Aligned Multi-Walled Carbon Nanotubes SKU: M2438 CAS Number: 308068-56-6 Options: - Length: 9 - 14 μm, 25 - 95 μm - Quantity: 500 mg, 1 g, 2 g, 5 g		Antimony Tin Oxide Nanopowder SKU: M2487 CAS Number: 128221-48-7 Options: Quantity: 100 g, 250 g
	Iron Oxide (Fe_3O_4) Nanopowder SKU: M2439 CAS Number: 1317-61-9 Options: - Average particle size: 30 - 50 nm - Quantity: 100 g, 250 g, 500 g		Copper Foil SKU: S2011 Options*: - Size (mm): 148 x 210 - Pack of 25, 210 x 297 - Pack of 25 - Thickness: 6 μm, 9 μm, 12 μm, 15 μm
	Iron Oxide (Fe_2O_3) Nanopowder SKU: M2440 CAS Number: 1309-37-1 Options: - Crystalline phase: α-Fe_2O_3, γ-Fe_2O_3 - Quantity: 10 g, 25 g, 50 g		
	Silicon Powder SKU: M2442 CAS Number: 7440-21-3 Options: - Average particle size: 90 nm, 10 μm - Quantity: 5 g, 10 g, 25 g, 50 g		
	Silicon Nanopowder SKU: M2442 CAS Number: 7440-21-3 Options: - Average particle size: 90 nm - Quantity: 5 g, 10 g		
	Silver Nanopowder SKU: M2444 CAS Number: 7440-22-4 Options: - Average particle size: 20 nm - Quantity: 5 g, 10 g, 25 g		
	Silver Nanowire Dispersion SKU: M2456 CAS Number: 7440-22-4 Options: Quantity: 25 mL, 50 mL		
	Tin (IV) Oxide Powder SKU: M2484 CAS Number: 18282-10-5 Options: Quantity: 50 g, 100 g, 250 g		


* Some options have been omitted, please see our website to browse our full range

Battery Binders


	Polyvinylidene Fluoride (PVDF) SKU: M2548 CAS Number: 24937-79-9 Options: • Quantity: 50 g, 100 g, 250 g
	Polytetrafluoroethylene (PTFE) SKU: M2549 CAS Number: 9002-84-0 Options: • Quantity: 50 g, 100 g, 250 g
	Carboxymethylcellulose Sodium SKU: M2550 CAS Number: 9004-32-4 Options: • Quantity: 100 g, 250 g

Optoelectronic Materials

Optoelectronic materials have the unique ability to interact with and manipulate photons and electrical energy, resulting from their structures, which are often conjugated. These materials are at the functional center of many modern devices, including simple light-emitting diodes (LEDs), displays, and photovoltaics.




Popular Products



4'-Pentyl-4-biphenylcarbonitrile

Nematic liquid crystal at room temperature


- 5CB
- Cyano and pentyl substituents



HAT6

Discotic liquid crystal with triphenylene core

- Discotic Liquid Crystals
- Transparent hole transport material



Poly(pentafluorostyrene)

Low refractive index, high heat and chemical resistance

- Highly fluorinated
- Enables post-functionalization

To help you choose the right material, full technical specifications are available on our website.

Ossila

M2433A1-250mg

HAT6




2,3,6,7,10,11-Hexakis(hexyloxy)triphenylene

CAS# 70351-86-9

Purity (>98%)
MW 829.24 g/mol











Solpro Business P

Photonic and Optical Materials

	Poly(1,1,1,3,3,3-hexafluoroisopropyl methacrylate-co-glycidyl methacrylate) 90/10 SKU: M2368 Options: • Quantity: 500 mg, 1 g, 5 g
	Poly(pentafluorostyrene-co-glycidyl methacrylate) 50/50 SKU: M2369 CAS Number: 224962-38-3 Options: • Quantity: 500 mg, 1 g, 5 g
	Poly(pentafluorostyrene) SKU: M2370 CAS Number: 26838-55-1 Options*: • Quantity: 1 g, 5 g

* Some options have been omitted, please see our website to browse our full range

Liquid Crystals

	2,3,6,7,10,11-Hexamethoxytriphenylene SKU: B4221 CAS Number: 808-57-1 Options: • Quantity: 1 g, 5 g		4-Pentylphenyl 4-pentylbenzoate SKU: M2366 CAS Number: 74305-48-9 Options: • Quantity: 10 g, 25 g, 50 g, 100 g
	4-(trans-4-Pentylcyclohexyl)benzonitrile SKU: M2358 CAS Number: 61204-01-1 Options: • Size: 10 g, 25 g, 50 g		Phenyl 1-hydroxy-2-naphthoate SKU: M2367 CAS Number: 132-54-7 Options: • Quantity: 25 g, 50 g, 100 g
	4'-Heptyl-4-biphenylcarbonitrile SKU: M2359 CAS Number: 41122-71-8 Options: • Quantity: 5 g, 10 g, 25 g		2,3,6,7,10,11-Hexakis(hexyloxy)triphenylene - HAT6 SKU: M2433 CAS Number: 70351-86-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	4'-Octyloxy-4-biphenylcarbonitrile SKU: M2360 CAS Number: 52364-73-5 Options: • Quantity: 5 g, 10 g, 25 g		
	4'-Pentyl-4-biphenylcarbonitrile SKU: M2361 CAS Number: 40817-08-1 Options: • Quantity: 10 g, 25 g, 50 g		
	4-Pentylphenyl 4-Methoxybenzoate SKU: M2363 CAS Number: 38444-13-2 Options: • Quantity: 25 g, 50 g, 100 g		
	4-Pentylphenyl 4-Methylbenzoate SKU: M2364 CAS Number: 50649-59-7 Options: • Quantity: 10 g, 25 g, 50 g		
	4-Pentylphenyl 4-propylbenzoate SKU: M2365 CAS Number: 50649-60-0 Options: • Quantity: 10 g, 25 g, 50 g, 100 g		


Fullerenes & NFAs

Fullerenes, soluble fullerene derivatives, and fullerene replacements are available for use as electron acceptors in OPV devices and n-type interfaces in perovskite solar cells. Classic derivatives of C₆₀ and C₇₀, functionalized with groups such as phenyls and esters, improve solubility and processability.

Non-fullerene acceptors (NFAs) such as Y6 are a promising alternative, offering greater flexibility to tune optical properties and energy levels. ITIC-based solar cells deliver higher thermal and photochemical stability, longer device lifetimes, and improved power conversion efficiencies.




Popular Products



PCBM

Popular, high-purity, soluble fullerene

- Tailored High Purity
- Efficient charge transfer



Y6

High-purity and efficiency, with wide optical absorption

- High power conversion efficiencies
- Electron deficient



C₆₀ Fullerene

Nobel Prize-winning electron acceptor molecule

- N-type semiconductor
- High purity

To help you choose the right material, full technical specifications are available on our website.

Ossila

Solpro Business Park Windsor SL

M2200A1-50mg






>98%
mw 1451.93

Y6(BTP-4F)
















2,2'-((2Z,2'Z)-((12,13-bis(2-ethylhexyl)-3,9-diundecyl-12,13-dihydro-[1,2,5]thiadiazolo[3,4-e]thieno[2'',3'':4',5'])thieno[2',3':4,5]pyrrolo[3,2-g]thieno[2',3':4,5]thieno[3,2-b]indole-2,10-diyl)bis(methanylylidene))bis(5,6-difluoro-3-oxo-2,3-dihydro-1H-indene-2,1-diylidene))dimalononitrile

CAS# 2304444-49-1
















Fullerenes for Photovoltaics and Biotechnology


	ICBA (Indene-C60 Bisadduct) SKU: M0181 CAS Number: 1207461-57-1 Options: • Quantity: 100 mg, 250 mg, 500 mg
	PCBM SKU: M111 CAS Number: 160848-22-6 Options: • Purity: >99%, >99.5%, >99.9% • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g
	PC70BM (C70-PCBM), PC71BM SKU: M114 CAS Number: 609771-63-3 Options: • PC70bm: >99% • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	C70 Fullerene SKU: M115 CAS Number: 115383-22-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	C60 Fullerene SKU: M117 CAS Number: 99685-96-8 Options: • Purity: > 99.5%, > 99.9%, > 99.9% Solvent Free • Quantity: 1 g, 5 g, 10 g, 25 g
	bis-PCBM (PC60BM Bisadduct) SKU: M191 CAS Number: 1048679-01-1 Options: • Quantity: 200 mg, 500 mg
	Fullerenol C60 SKU: M2470 Options: • Purity: Fulleranol C60 99.5%, Fulleranol C60 99.9% Solvent Free • Quantity: 500 mg, 1 g, 5 g
	Fullerene C60 Malonic Acid SKU: M2471 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

Non-Fullerene Acceptors

	ITIC SKU: M1101 CAS Number: 1664293-06-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		o-IDTBR SKU: M2238 CAS Number: 2077945-91-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	ITIC-M, IT-M SKU: M2073 CAS Number: 2047352-80-5 / 2047352-83-8 / 2047352-86-1 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g		IHIC SKU: M2239 CAS Number: 2127354-15-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	ITIC-Th, IT-Th SKU: M2074 CAS Number: 1889344-13-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		IEICO-4F, IOTIC-4F SKU: M2240 CAS Number: 2089044-02-8 Options: • Quantity: 50 mg, 100 mg, 250 mg
	ITIC-4F, ITIC-2F, IT-4F SKU: M2075 CAS Number: 2097998-59-7 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g		IDIC-4F SKU: M2241 CAS Number: 2229934-71-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	ITIC-DM, IT-2M SKU: M2148 CAS Number: 2047352-92-9 Options: • Quantity: 100 mg		FBR SKU: M2242 CAS Number: 1644381-95-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	ITIC-4Cl, ITIC-DCI, IT-4Cl SKU: M2149 CAS Number: 2253663-81-7 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg		BTP-4F-12, Y6-BO (Y12) SKU: M2245 CAS Number: 2389125-23-7 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g
	Y6, BTP-4F SKU: M2200 CAS Number: 2304444-49-1 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g		COTIC-4F SKU: M2246 CAS Number: 2328118-91-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Y7, BTP-4Cl SKU: M2212 CAS Number: 2414918-25-3 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg		BTP-eC9 SKU: M2249 CAS Number: 2598965-39-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g


* Some options have been omitted, please see our website to browse our full range

	TPT10 SKU: M2250 CAS Number: 2414381-17-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		L8-BO, L8-BO-2F SKU: M2304 CAS Number: 2668341-40-8 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg
	IDT-2Br SKU: M2252 CAS Number: 2042521-91-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		L8-BO-F SKU: M2306 Options: • Quantity: 50 mg, 100 mg, 250 mg
	TPT10-C8C12 SKU: M2253 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		N3 - Y6N3 SKU: M2324 Options: • Quantity: 50 mg, 100 mg, 250 mg
	IDIC SKU: M2261 CAS Number: 1883441-92-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BTP-4Cl-12 SKU: M2327 CAS Number: 2447642-41-1 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg
	IDIC-4Cl SKU: M2262 CAS Number: 2361961-01-3 Options: • Quantity: 100 mg, 250 mg, 500 mg		6TIC SKU: M2329 CAS Number: 2244414-53-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DTY6 SKU: M2274 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg		COTIC-4Cl SKU: M2340 CAS Number: 2328092-51-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	IO-4Cl SKU: M2298 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		IEICO-4Cl SKU: M2341 CAS Number: 2240998-88-1 Options: • Quantity: 100 mg
	eC9-2Cl SKU: M2301 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		ZY-4Cl SKU: M2350 CAS Number: 2703920-28-7 Options: • Quantity: 50 mg, 100 mg, 250 mg


	EH-IDTBR SKU: M2413 CAS Number: 2102510-60-9 Options: • Quantity: 100 mg, 250 mg
---	--

Semiconducting Polymers

Semiconducting polymers for bulk heterojunction, OPV, OLED, OFET and perovskite interfaces and solar cell research. OFETs (organic field-effect transistors) and OLEDs (organic light-emitting diodes) are at the forefront of organic electronics. These stable polymer materials have semiconducting properties and are both solution-processable and morphologically advantageous. OPVs are distinguished by their flexibility, lightweight nature, and potential for low-cost production compared to traditional silicon-based solar cells.




Popular Products



PNDI(2OD)2T

High-mobility n-type polymer semiconductor


- Broad light absorption
- High electron affinity
- High purity via Soxhlet extraction



DPP-DTT

High-mobility p-type polymer with versatile applications

- High molecular weight
- High purity via Soxhlet extraction
- Batch-specific GPC data



P4T2F-HD

Strong interchain interaction for enhanced polymer crystallinity

- Designed for OSC active material packing
- High purity via Soxhlet extraction
- Batch-specific GPC data

To help you choose the right material, full technical specifications are available on our website.

Ossila

M0311A10-100mg

















| Mw 130,679

DPP-DTT
















Poly[2,5-(2-octyldodecyl)-3,6-diketopyrrolopyrrole-alt-5,5'-(2,5'-thieno [3,2-b]thiophene)]

CAS# 1260685-66-7















OPV Polymers

	PCDTBT SKU: M0131 CAS Number: 958261-50-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		PBDB-T (PCE12) SKU: M1001 CAS Number: 1415929-80-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PTB7 SKU: M0211 CAS Number: 1266549-31-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		P3HT SKU: M1011 CAS Number: 104934-50-1 Options*: • Quantity: 1 g, 5 g, 10 g, 25 g
	PTB7-Th, PCE10, PBDTTT-EFT SKU: M0261 CAS Number: 1469791-66-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PNDI(2OD)2T SKU: M1201 CAS Number: 1100243-40-0 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPP-DTT, PDPP2T-TT-OD SKU: M0311 CAS Number: 1260685-66-2 / 1444870-74-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		PBDTT-DPP SKU: M1401 CAS Number: 1380582-98-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PDBPyBT (DPPDPyBT, high-mobility n-type polymer) SKU: M0321 CAS Number: 1455028-36-0 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg		PBDD4T SKU: M2049 CAS Number: 1439937-09-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PDPP4T SKU: M0331 CAS Number: 1267540-03-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDD4T-2F SKU: M2050 CAS Number: 1890205-85-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	F8T2 SKU: M0501 CAS Number: 210347-56-1 Options*: • Quantity: 250 mg, 500 mg, 1 g		PDPP4T-2F SKU: M2051 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PDPP3T SKU: M0991 CAS Number: 1198291-01-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PNF222 SKU: M2052 CAS Number: 1773524-46-1 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g

* Some options have been omitted, please see our website to browse our full range
















	PBDTPD SKU: M2053 CAS Number: 1223479-75-1 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PNDI(2OD)T, PCE8 SKU: M2087 CAS Number: 1236111-76-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PBDTTTPD SKU: M2054 CAS Number: 1426534-44-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PSBTBT SKU: M2095 CAS Number: 1089687-02-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PJ71 SKU: M2056 CAS Number: 2035466-89-6 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g		PCPDTBT SKU: M2096 CAS Number: 920515-34-0 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TQ1, PTQ1 SKU: M2076 CAS Number: 565228-37-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDB-T-2Cl (PM7) SKU: M2134 CAS Number: 2239295-71-5 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PNDI(2HD)T SKU: M2081 CAS Number: 1675973-74-6 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDB-T-2F (PM6) SKU: M2150 CAS Number: 1802013-83-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PNDI(2HD)2T SKU: M2082 CAS Number: 1459168-70-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PFO-co-MEH-PPV SKU: M2197 CAS Number: 1383605-56-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg
	PffBT4T-C9C13 (PCE12) SKU: M2083 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PTEBS Water Soluble Polythiophenes SKU: M2197 CAS Number: 744198-90-1 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PNDI-T10, PCE9 SKU: M2086 CAS Number: 1977539-03-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PDCBT SKU: M2233 CAS Number: 1609536-17-5 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

* Some options have been omitted, please see our website to browse our full range








	D18 SKU: M2248 CAS Number: 2433725-54-1 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDB-TCI-T1 SKU: M2282 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PBBTCD SKU: M2256 CAS Number: 1334032-14-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PTQ10 SKU: M2283 CAS Number: 2270233-86-6 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PBBTsID SKU: M2258 CAS Number: 1334032-16-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		J52 SKU: M2288 CAS Number: 1887136-01-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PBDB-TF-T1 SKU: M2259 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		J52-Cl SKU: M2297 CAS Number: 1887136-01-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PTO2 SKU: M2260 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBQx-TF (PCE19) SKU: M2383 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PTQ11 SKU: M2277 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PCE11 (PffBT4T-2OD) SKU: M303 CAS Number: 1644164-62-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g
	PCITO2 SKU: M2278 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	D18-Cl SKU: M2281 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		

* Some options have been omitted, please see our website to browse our full range

OFET & OLED Polymer Materials

















	PCDTBT SKU: M0131 CAS Number: 958261-50-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		MEH-PPV SKU: M0541 CAS Number: 138184-36-8 Options: • Molecular weight and pDI: Mw = 204836 PDI = 6.96 • Quantity: 250 mg, 500 mg, 1 g
	PFO (F8) SKU: M0161 CAS Number: 19456-48-5 Options*: • Quantity: 500 mg, 1 g		P3HT SKU: M1011 CAS Number: 104934-50-1 Options*: • Quantity: 1 g, 5 g, 10 g, 25 g
	F8BT (PF8T) SKU: M0231 CAS Number: 210347-52-7 Options*: • Quantity: 250 mg, 500 mg, 1 g		PNDI(2OD)2T SKU: M1201 CAS Number: 1100243-40-0 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPP-DTT, PDPP2T-TT-OD SKU: M0311 CAS Number: 1260685-66-2 / 1444870-74-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		PTAA for OFET Applications SKU: M171
	PDBPyBT (DPPDPyBT, high-mobility n-type polymer) SKU: M0321 CAS Number: 1455028-36-0 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg		PNF222 SKU: M2052 CAS Number: 1773524-46-1 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g
	PDPP4T SKU: M0331 CAS Number: 1267540-03-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PNDI(2HD)T SKU: M2081 CAS Number: 1675973-74-6 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	F8T2 SKU: M0501 CAS Number: 210347-56-1 Options*: • Quantity: 250 mg, 500 mg, 1 g		PNDI(2HD)2T SKU: M2086 CAS Number: 1459168-70-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PolyTPD SKU: M0521 CAS Number: 472960-35-3 Options: • Grade: M0521A5 (Mw 63630; PDI 5.34) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PNDI-T10, PCE9 SKU: M2086 CAS Number: 1977539-03-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

* Some options have been omitted, please see our website to browse our full range


	PNDI(2OD)T, PCE8 SKU: M2087 CAS Number: 1236111-76-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PSBTBT SKU: M2095 CAS Number: 1089687-02-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PCPDTBT SKU: M2096 CAS Number: 920515-34-0 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PFO-co-MEH-PPV SKU: M2196 CAS Number: 1383605-56-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg
	P4T2F-HD SKU: M2348 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Poly(9-vinylcarbazole) (PVK) SKU: M2384 CAS Number: 25067-59-8 Options: • Batch specification: M2384A6 (Mw = 235700; PDI = 2.22) • Quantity: 5 g, 10 g
	PF8Cz SKU: M2421 Options: • Molecular weight: M2421A2 • Quantity: 100 mg, 250 mg, 500 mg, 1 g

* Some options have been omitted, please see our website to browse our full range

Interface Polymers

	PCDTBT SKU: M0131 CAS Number: 958261-50-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		PDPP3T SKU: M0991 CAS Number: 1198291-01-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPP-DTT, PDPP2T-TT-OD SKU: M0311 CAS Number: 1260685-66-2 / 1444870-74-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		P3HT SKU: M1011 CAS Number: 104934-50-1 Options*: • Quantity: 1 g, 5 g, 10 g, 25 g
	PDBPyBT (DPPDPyBT, high-mobility n-type polymer) SKU: M0321 CAS Number: 1455028-36-0 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg		PFN-FP SKU: M2078 CAS Number: 673474-79-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg
	PDPP4T SKU: M0331 CAS Number: 1267540-03-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PFN-I SKU: M2079 CAS Number: 1687752-60-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg
	F8T2 SKU: M0501 CAS Number: 210347-56-1 Options*: • Quantity: 250 mg, 500 mg, 1 g		PFN-FP-I SKU: M2080 CAS Number: 1026670-72-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg
	PTAA for Perovskite Applications SKU: M0511 CAS Number: 1333317-99-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PTEBS Water Soluble Polythiophenes SKU: M2197 CAS Number: 744198-90-1 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PolyTPD SKU: M0521 CAS Number: 472960-35-3 Options: • Grade: M0521A5 (Mw 63630; PDI 5.34) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PFN, PFN-DOF SKU: M222 CAS Number: 673474-74-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TFB SKU: M0981 CAS Number: 220797-16-0 Options*: • Quantity: 100 mg, 250 mg, 500 mg		PFN-bromide, PFN-Br SKU: M2230 CAS Number: 889672-99-5 Options*: • Quantity: 100 mg, 250 mg, 500 mg

* Some options have been omitted, please see our website to browse our full range



PF8Cz


SKU: M2421

Options:


- Molecular weight: M2421A2
- Quantity: 100 mg, 250 mg, 500 mg, 1 g

Luminosyn™ Polymers

The Luminosyn™ range by Ossila offers high-purity, batch-specific polymer semiconductors designed to give researchers consistent performance and reproducible results across experiments. Every material is subjected to a rigorous Soxhlet extraction process, removing impurities and ensuring exceptional quality. This level of consistency and purity makes Luminosyn™ polymers a dependable choice for cutting-edge research in electronics and optoelectronics.




Popular Products



D18

Narrow bandgap copolymer
with high hole mobility


- High purity via Soxhlet extraction
- Can be processed in non-halogenated solvents



F8T2

Popular, high-purity material
for organic electronics

- Absorption in the blue region
- Polymer donor
- High purity via Soxhlet extraction



J52

Medium bandgap
semiconducting copolymer
donor

- Batch-specific GPC data
- Maximum absorption around 450 - 600 nm
- High purity via Soxhlet extraction

To help you choose the right material, full technical specifications are available on our website.

Ossila

M2283A1-100mg

PTQ10

Poly[[[6,7-difluoro[[(2-hexyldecyloxy)-5,8-quinoline-2,5-thiophenediyl]]

CAS# 2270233-86-6

Solpro Bu

Solpro Business

















Ossila

M2200A1-100mg

















Y6(TTPPTI-4F)

2,2'-(2,2,2,2-tetrafluoro-1,1'-bis(2-ethyl-3,9-diundecyl-1,3,4,5-tetrahydro[1,2,5]thiadiazol[4,5-c]thienyl)-4,4'-diyl)-3,3'-bipyridine-2,10-diol bis(methanesulfonyl)dimaleonitrile










Luminosyn™ Polymers

	PCDTBT SKU: M0131 CAS Number: 958261-50-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		PNDI(2OD)2T SKU: M1201 CAS Number: 1100243-40-0 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PTB7 SKU: M0211 CAS Number: 1266549-31-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDTT-DPP SKU: M1401 CAS Number: 1380582-98-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PTB7-Th, PCE10, PBDTTT-EFT SKU: M0261 CAS Number: 1469791-66-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDD4T SKU: M2049 CAS Number: 1439937-09-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPP-DTT, PDPP2T-TT-OD / 1444870-74-9 SKU: M0311 CAS Number: 1260685-66-2 / 1444870-74-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		PBDD4T-2F SKU: M2050 CAS Number: 1890205-82-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PDPP4T SKU: M0331 CAS Number: 1267540-03-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PNF222 SKU: M2052 CAS Number: 1773524-46-1 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g
	F8T2 SKU: M0501 CAS Number: 210347-56-1 Options*: • Quantity: 250 mg, 500 mg, 1 g		PBDTTPD SKU: M2053 CAS Number: 1223479-75-1 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PDPP3T SKU: M0991 CAS Number: 1198291-01-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDTTTPD SKU: M2054 CAS Number: 1426534-44-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PBDB-T (PCE12) SKU: M1001 CAS Number: 1415929-80-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PJ71 SKU: M2056 CAS Number: 2035466-89-6 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g

* Some options have been omitted, please see our website to browse our full range

	TQ1, PTQ1 SKU: M2076 CAS Number: 965228-37-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBBT5ID SKU: M2258 CAS Number: 1334032-16-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PNDI(2HD)T SKU: M2081 CAS Number: 1675973-74-6 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDB-TF-T1 SKU: M2259 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PNDI-T10, PCE9 SKU: M2086 CAS Number: 1977539-03-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PTO2 SKU: M2260 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PBDB-T-2Cl (PM7) SKU: M2134 CAS Number: 2239295-71-5 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PTQ11 SKU: M2277 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PBDB-T-2F (PM6) SKU: M2150 CAS Number: 1802013-83-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PCITO2 SKU: M2278 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PDCBT SKU: M2233 CAS Number: 1609536-17-5 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		D18-Cl SKU: M2281 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	D18 SKU: M2248 CAS Number: 2433725-54-1 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PBDB-TCI-T1 SKU: M2282 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PBBTCD SKU: M2256 CAS Number: 1334032-14-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PTQ10 SKU: M2283 CAS Number: 2270233-86-6 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

* Some options have been omitted, please see our website to browse our full range

	J52 SKU: M2288 CAS Number: 1887136-01-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		P(g2T-T) SKU: M2477 CAS Number: 1974310-62-7 Options*: • Quantity: 100 mg
	J52-Cl SKU: M2297 CAS Number: 1887136-01-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	P4T2F-HD SKU: M2348 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	P5TCN-F25 SKU: M2349 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	PBQx-TF (PCE19) SKU: M2383 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	Poly(9-vinylcarbazole) (PVK) SKU: M2384 CAS Number: 25067-59-8 Options: • Batch specification: M2384A6 (Mw = 235700; PDI = 2.22) • Quantity: 5 g, 10 g		
	PF8Cz SKU: M2421 Options: • Molecular weight: M2421A2 • Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	P(g2T-TT) SKU: M2475 CAS Number: 2058240-28-9 Options*: • Quantity: 100 mg, 250 mg		

* Some options have been omitted, please see our website to browse our full range

Perovskite Materials

Perovskite materials are defined by their unique ABX₃ crystal structure, where "A" and "B" are cations and "X" is an anion. These materials often incorporate multiple A, B, and X components within a single crystal to optimize properties such as bandgap tuning and stability. The composition and stoichiometry are critical, significantly influencing device performance. Perovskites can be tailored for various applications, with current uses including solar cells, LEDs, sensors, and photodetectors. Their potential applications are vast and continue to be actively explored.




Popular Products



Perovskite Inks

For processing in air or nitrogen with triple cation options available


- Specially formulated
- For fabrication of photovoltaic devices



Methylammonium Iodide (MAI)

High-purity, low-price precursor for the synthesis of perovskites

- For synthesis of MAPbI₃ based solar cells
- Recrystallised variants for added purity



Formamidinium Iodide (FAI)

High-purity precursor for optimum solar conversion efficiencies

- For synthesis of FAPbI₃ based solar cells
- Recrystallised variants for added purity

To help you choose the right material, full technical specifications are available on our website.



Sheffield S4 799B

GER

Harmful if swallowed. May cause eye irritation. Harmful to the unborn child. May cause delayed or repeated effects. Keep away from children. - No smoking. Avoid breathing surfaces/spray. Avoid breathing. - SWALLOWED: Call a doctor. IF ON SKIN: Wash well.

Sheffield S4 799B

GER

Harmful if swallowed. May cause eye irritation. Harmful to the unborn child. May cause delayed or repeated effects. Keep away from children. - No smoking. Avoid breathing surfaces/spray. Avoid breathing. - SWALLOWED: Call a doctor. IF ON SKIN: Wash well.

Sheffield S4 799B

GER

Harmful if swallowed. May cause eye irritation. Harmful to the unborn child. May cause delayed or repeated effects. Keep away from children. - No smoking. Avoid breathing surfaces/spray. Avoid breathing. - SWALLOWED: Call a doctor. IF ON SKIN: Wash well.

Sheffield S4 799B

GER

Harmful if swallowed. May cause eye irritation. Harmful to the unborn child. May cause delayed or repeated effects. Keep away from children. - No smoking. Avoid breathing surfaces/spray. Avoid breathing. - SWALLOWED: Call a doctor. IF ON SKIN: Wash well.

799B, UK

Sheffield S4 799B

GER

Harmful if swallowed. May cause eye irritation. Harmful to the unborn child. May cause delayed or repeated effects. Keep away from children. - No smoking. Avoid breathing surfaces/spray. Avoid breathing. - SWALLOWED: Call a doctor. IF ON SKIN: Wash well.

Sheffield S4 799B

GER

Harmful if swallowed. May cause eye irritation. Harmful to the unborn child. May cause delayed or repeated effects. Keep away from children. - No smoking. Avoid breathing surfaces/spray. Avoid breathing. - SWALLOWED: Call a doctor. IF ON SKIN: Wash well.

Sheffield S4 799B

GER

Harmful if swallowed. May cause eye irritation. Harmful to the unborn child. May cause delayed or repeated effects. Keep away from children. - No smoking. Avoid breathing surfaces/spray. Avoid breathing. - SWALLOWED: Call a doctor. IF ON SKIN: Wash well.

Perovskite Inks



Perovskite Precursor Ink for Air Processing
SKU: I101
Options:
• Quantity: 10 x 0.5 ml (5 ml total), 3 x 10 ml (30 ml total)








Perovskite Precursor Ink for Nitrogen Processing
SKU: I201
Options:
• Quantity: 10 x 0.5 ml (5 ml total), 3 x 10ml (30 ml total)



















Triple Cation Perovskite Precursor Ink (for Nitrogen Processing)
SKU: I301
Options:
• Size: 10 x 0.5 ml (5 ml total), 3 x 10 ml (30 ml total)

Perovskite Precursor Materials

















	n-Butylammonium Iodide SKU: M2057 CAS Number: 36945-08-1 Options: • Quantity: 10 g		Benzylammonium Iodide (BzAI) SKU: M2065 CAS Number: 45579-91-7 Options: • Quantity: 10 g
	n-Butylammonium Bromide SKU: M2058 CAS Number: 15567-09-6 Options: • Quantity: 25 g		Ethylammonium Iodide, EAI SKU: M2066 CAS Number: 506-58-1 Options: • Quantity: 10 g, 25 g
	iso-Butylammonium Iodide SKU: M2059 CAS Number: 205508-75-4 Options: • Quantity: 10 g		Ethylammonium Bromide, EABr SKU: M2067 CAS Number: 593-55-5 Options: • Quantity: 25 g
	iso-Butylammonium Bromide SKU: M2060 CAS Number: 74098-36-5 Options: • Quantity: 10 g		Dimethylammonium Iodide (DMAI) SKU: M2068 CAS Number: 51066-74-1 Options: • Quantity: 10 g
	t-Butylammonium Iodide SKU: M2061 CAS Number: 39557-45-4 Options: • Quantity: 5 g, 10 g		Dimethylammonium Bromide (DMABr) SKU: M2069 CAS Number: 6912-12-5 Options: • Quantity: 25 g
	t-Butylammonium Bromide SKU: M2062 CAS Number: 60469-70-7 Options: • Quantity: 25 g		Imidazolium Iodide SKU: M2070 CAS Number: 68007-08-9 Options: • Quantity: 10 g
	n-Propylammonium Iodide SKU: M2063 CAS Number: 14488-45-0 Options: • Quantity: 5 g, 10 g		Imidazolium Bromide SKU: M2071 CAS Number: 101023-55-6 Options: • Quantity: 10 g, 25 g
	Acetamidinium Iodide SKU: M2064 CAS Number: 1452099-14-7 Options: • Quantity: 10 g		Guanidinium Iodide SKU: M2072 CAS Number: 19227-70-4 Options: • Quantity: 25 g

	Guanidinium Bromide SKU: M2077 CAS Number: 19244-98-5 Options: • Quantity: 25 g
	Methylammonium Iodide (MAI) SKU: M271 CAS Number: 14965-49-2 Options: • Quantity: 5 g, 10 g, 25 g, 100 g, 500 g • Quantity/grade: 99.9% purity, >99.9% purity - recrystallised
	Formamidinium Iodide (FAI) SKU: M551 CAS Number: 879643-71-7 Options: • Grade: 99.995% purity, >99.995% purity - recrystallised • Quantity: 5 g, 10 g, 25 g
	Formamidinium Bromide (FABr) SKU: M561 CAS Number: 146958-06-7 Options: • Grade: 99% purity, >99% purity - recrystallised • Quantity: 5 g, 10 g, 25 g
	Methylammonium Bromide (MABr) SKU: M571 CAS Number: 6876-37-5 Options: • Grade: 99.995% purity, >99.995% purity - recrystallised • Quantity: 5 g, 10 g, 25 g
















Perovskite Interface Materials














	Bathocuproine (BCP) SKU: B231 CAS Number: 4733-39-5 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 1 g, 5 g		TFB SKU: M0981 CAS Number: 220797-16-0 Options*: • Quantity: 100 mg, 250 mg, 500 mg
	PCDTBT SKU: M0131 CAS Number: 958261-50-2 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 2 g		P3HT SKU: M1011 CAS Number: 104934-50-1 Options*: • Quantity: 1 g, 5 g, 10 g, 25 g
	PFO (F8) SKU: M0161 CAS Number: 19456-48-5 Options*: • Quantity: 500 mg, 1 g		FK102-Co(II)PF6 Salt SKU: M2007 CAS Number: 1392221-69-0 Options: • Quantity: 10 g
	F8BT (PFBT) SKU: M0231 CAS Number: 210347-52-7 Options*: • Quantity: 250 mg, 500 mg, 1 g		FK102-Co(II)TFSI Salt SKU: M2019 Options: • Quantity: 10 g
	TPD - N,N'-Bis(3-methylphenyl)-N,N'-diphenylbenzidine SKU: M0401 CAS Number: 65181-78-4 Options: • Quantity: 1 g, 5 g		FK102-Co(III)PF6 Salt SKU: M2020 CAS Number: 1346416-71-4 Options: • Quantity: 10 g
	PTAA for Perovskite Applications SKU: M0511 CAS Number: 1333317-99-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		FK102-Co(III)TFSI Salt SKU: M2021 CAS Number: 2057441-94-6 Options: • Quantity: 10 g
	PolyTPD SKU: M0521 CAS Number: 472960-35-3 Options: • Grade: M0521A5 (Mw 63630; PDI 5.34) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TQ1, PTQ1 SKU: M0276 CAS Number: 565228-37-7 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4-tert-Butylpyridine SKU: M0591 CAS Number: 3978-81-2 Options: • Quantity: 25 g, 50 g, 100 g		PFN-FP SKU: M2078 CAS Number: 673474-79-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg




* Some options have been omitted, please see our website to browse our full range

	PFN-I SKU: M2079 CAS Number: 1687752-60-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg		B4PyPPm SKU: M2176 CAS Number: 1097652-83-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PFN-FP-I SKU: M2080 CAS Number: 1026670-72-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg		Spiro-TTB SKU: M2192 CAS Number: 515834-67-0 Options: • Quantity: 500 mg, 1 g, 5 g
	MADN SKU: M2092 CAS Number: 804560-00-7 Options: • Quantity: 1 g, 5 g		PFN, PFN-DOF SKU: M222 CAS Number: 673474-74-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PSBTBT SKU: M2205 CAS Number: 1089687-02-4 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PFN-bromide, PFN-Br SKU: M2230 CAS Number: 889672-99-5 Options*: • Quantity: 100 mg, 250 mg, 500 mg
	PCPDTBT SKU: M2096 CAS Number: 920515-34-0 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Guanidinium Tetrafluoroborate SKU: M2422 CAS Number: 36595-00-3 Options: • Quantity: 5 g, 10 g, 25 g, 50 g
	TTPA SKU: M2129 CAS Number: 177799-16-5 Options: • Quantity: 250 mg, 500 mg, 1 g		2PACz SKU: M2425 CAS Number: 20999-38-6 Options: • Quantity: 500 mg, 1 g
	B3PymPm SKU: M2174 CAS Number: 925425-96-3 Options: • Purity: Sublimed >99.0% purity, Unsublimed >98.0% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g		Br-2PACz SKU: M2427 CAS Number: 2762888-11-7 Options: • Quantity: 500 mg, 1 g
	B4PymPm SKU: M2175 CAS Number: 1030380-51-8 Options: • Quantity: 250 mg, 500 mg, 1 g		MeO-2PACz SKU: M2427 CAS Number: 2377770-18-6 Options: • Quantity: 500 mg, 1 g

















* Some options have been omitted, please see our website to browse our full range










	2,3,6,7,10,11-Hexakis(hexyloxy)triphenylene – HAT6 SKU: M2433 CAS Number: 70351-86-9 Options: • Quantity: 250 mg, 500 mg, 1 g		Me-4PACz SKU: M2454 CAS Number: 2747959-96-0 Options: • Quantity: 500 mg, 1 g, 5 g
	Nickel (II) Oxide Powder SKU: M2441 CAS Number: 1313-99-1 Options: • Particle size: 8 - 18 nm • Quantity: 10 g, 25 g		Br-2EPT SKU: M2455 CAS Number: 2826271-17-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	DMACPA SKU: M2445 CAS Number: 2971088-37-4 Options: • Quantity: 250 mg, 500 mg, 1 g		TDPA-Cl SKU: M2457 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4PADCB SKU: M2446 CAS Number: 2882156-63-8 Options: • Quantity: 250 mg, 500 mg, 1 g		2BrPXZPA SKU: M2458 Options: • Quantity: 250 mg, 500 mg, 1 g
	MPA-CPA SKU: M2447 CAS Number: 2212003-31-9 Options: • Quantity: 250 mg, 500 mg, 1 g		2BrPTZPA SKU: M2459 CAS Number: 2913184-53-7 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	Cl-2PACz SKU: M2451 CAS Number: 3036926-72-1 Options: • Quantity: 500 mg, 1 g, 5 g		3PATAT-C3 SKU: M2460 CAS Number: 3048690-10-1 Options: • Quantity: 250 mg, 500 mg, 1 g
	4PACz SKU: M2452 CAS Number: 20999-36-4 Options: • Quantity: 500 mg, 1 g, 5 g		9CAA SKU: M2461 CAS Number: 524-80-1 Options: • Quantity: 500 mg, 1 g, 5 g
	MeO-4PACz SKU: M2453 CAS Number: 2922526-56-3 Options: • Quantity: 500 mg, 1 g, 5 g		DCB-BPA SKU: M2462 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g

	9CPA SKU: M2463 CAS Number: 6622-54-4 Options: • Quantity: 500 mg, 1 g, 5 g		Tantalum (V) Oxide Powder SKU: M2489 CAS Number: 1314-61-0 Options: • Quantity: 10 g, 25 g, 50 g
	3-BPIC-F SKU: M2464 CAS Number: 3038707-50-2 Options: • Quantity: 250 mg, 500 mg, 1 g		Spiro-OMeTAD (Spiro-MeOTAD) SKU: M291 CAS Number: 207739-72-8 Options: • Grade: Sublimed (>99.5% purity), Unsublimed (>99.0% purity) • Quantity: 500 mg, 1 g, 5 g, 10 g
	3-BPIC SKU: M2465 CAS Number: 3038707-49-9 Options: • Quantity: 250 mg, 500 mg, 1 g		F4TCNQ SKU: M351 CAS Number: 29261-33-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g
	Me-PhpPACz SKU: M2468 CAS Number: 3046309-16-1 Options: • Quantity: 250 mg, 500 mg, 1 g		NPB (NPD) SKU: M361 CAS Number: 123847-85-8 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 1 g, 5 g
	Phosphorylcholine chloride SKU: M2468 CAS Number: 107-73-3 Options: • Quantity: 1 g, 5 g, 10 g, 25 g		CBP, 4,4'-Bis(N-carbazolyl)-1,1'-biphenyl SKU: M391 CAS Number: 58328-31-7 Options: • Grade: Sublimed Grade (>99.5% purity) • Quantity: 1 g, 5 g
	Ph-4PACz SKU: M2481 CAS Number: 2814500-04-2 Options: • Quantity: 250 mg, 500 mg, 1 g		LiTFSI SKU: M531 CAS Number: 90076-65-6 Options: • Quantity: 50 g
	Magnesium Oxide Powder SKU: M2483 CAS Number: 1309-48-4 Options: • Quantity: 100 g, 250 g, 500 g		TAZ, 3-(Biphenyl-4-yl)-5-(4-tert-butylphen- enyl)-4-phenyl-4H-1,2,4-triazole SKU: M601 CAS Number: 150405-69-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	Tin (IV) Oxide Powder SKU: M2484 CAS Number: 18282-10-5 Options: • Quantity: 50 g, 100 g, 250 g		m-MTDATA SKU: M621 CAS Number: 124729-98-2 Options: • Quantity: 250 mg, 500 mg, 1 g

	TPBi SKU: M651 CAS Number: 192198-85-9 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 250 mg, 500 mg, 1 g
	HATCN SKU: M791 CAS Number: 105598-27-4 Options: • Grade: Sublimed (>99.0%), Unsublimed (>98.0%) • Quantity: 500 mg, 1 g, 5 g, 10 g
	BPhen (Bathophenanthroline) SKU: M961 CAS Number: 1662-01-7 Options: • Quantity: 1 g, 5 g


Self-Assembled Monolayers (SAMs)

	2PACz SKU: M2425 CAS Number: 20999-38-6 Options: • Quantity: 500 mg, 1 g		MeO-4PACz SKU: M2453 CAS Number: 2922526-56-3 Options: • Quantity: 500 mg, 1 g, 5 g
	Br-2PACz SKU: M2426 CAS Number: 2762888-11-7 Options: • Quantity: 500 mg, 1 g		Me-4PACz SKU: M2454 CAS Number: 2747959-96-0 Options: • Quantity: 500 mg, 1 g, 5 g
	MeO-2PACz SKU: M2427 CAS Number: 2377770-18-6 Options: • Quantity: 500 mg, 1 g		Br-2EPT SKU: M2455 CAS Number: 2826271-17-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	DMACPA SKU: M2445 CAS Number: 2971088-37-4 Options: • Quantity: 250 mg, 500 mg, 1 g		TDPA-Cl SKU: M2457 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	4PADCB SKU: M2446 CAS Number: 2882156-63-8 Options: • Quantity: 250 mg, 500 mg, 1 g		2BrPXZPA SKU: M2458 Options: • Quantity: 250 mg, 500 mg, 1 g
	MPA-CPA SKU: M2447 CAS Number: 2212003-31-9 Options: • Quantity: 250 mg, 500 mg, 1 g		2BrPTZPA SKU: M2459 CAS Number: 2913184-53-7 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g
	Cl-2PACz SKU: M2451 CAS Number: 3036926-72-1 Options: • Quantity: 500 mg, 1 g, 5 g		3PATAT-C3 SKU: M2460 CAS Number: 3048690-10-1 Options: • Quantity: 250 mg, 500 mg, 1 g
	4PACz SKU: M2452 CAS Number: 20999-36-4 Options: • Quantity: 500 mg, 1 g, 5 g		9CAA SKU: M2461 CAS Number: 524-80-1 Options: • Quantity: 500 mg, 1 g, 5 g


	DCB-BPA SKU: M2462 Options: Quantity: 250 mg, 500 mg, 1 g, 5 g
	Ph-4PACz SKU: M2481 CAS Number: 2814500-04-2 Options: Quantity: 250 mg, 500 mg, 1 g
	9CPA SKU: M2463 CAS Number: 6622-54-4 Options: Quantity: 500 mg, 1 g, 5 g
	3-BPIC-F SKU: M2464 CAS Number: 3038707-50-2 Options: Quantity: 250 mg, 500 mg, 1 g
	3-BPIC SKU: M2465 CAS Number: 3038707-49-9 Options: Quantity: 250 mg, 500 mg, 1 g
	Me-PhpPACz SKU: M2467 CAS Number: 3046309-16-1 Options: Quantity: 250 mg, 500 mg, 1 g
	Phosphorylcholine chloride SKU: M2468 CAS Number: 107-73-3 Options: Quantity: 1 g, 5 g, 10 g, 25 g
	PTZ-CPA SKU: M2474 Options: Quantity: 250 mg, 500 mg, 1 g
	Py3 SKU: M2478 CAS Number: 1520027-62-6 Options: Quantity: 250 mg, 500 mg, 1 g

DSSC Dyes & Electrolytes

Dye-sensitized solar cells (DSSCs) are thin-film solar cells designed around organic dyes in electrochemical cells. Considered to be the third generation of solar cell, DSSCs are comparatively cheap and straightforward to produce and do not require the use of toxic chemicals. Our range of DSSC materials include the foundation dye N3, black dye N749, and cobalt electrolytes.




Popular Products



N749 Dye

High-purity ruthenium terpyridine photosensitizer


- Absorbs light over a very broad spectrum



Z907 Dye

High-purity ruthenium bipyridine photosensitizer

- Amphiphilic
- Excellent thermal stability
- High efficiency



FK102-Co(III)TFSI Salt

High-purity cobalt pyrazolopyridine complex

- DSSC Electrolyte
- P-type dopant
- Large stabilizing TFSI⁻ anion

To help you choose the right material, full technical specifications are available on our website.

Ossila

M291-500mg

Spiro-OMeTAD












2,2',7,7'-Tetrakis[N,N-di(4-methoxyphenyl)amino]-
9,9'-spirobifluorene

CAS# 207739-72-8

Purity: >99.0%
MW 12


Solpro Business Park Winc

DSSC Dyes & Electrolytes


	N749 Black Dye SKU: M2001 CAS Number: 359415-47-7 Options: • Quantity: 100 mg, 250 mg		FK102-Co(III)PF6 Salt SKU: M2020 CAS Number: 1346416-71-4 Options: • Quantity: 10 g
	Z907 Dye SKU: M2002 CAS Number: 502693-09-6 Options: • Quantity: 1 g		FK102-Co(III)TFSI Salt SKU: M2021 CAS Number: 2057441-94-6 Options: • Quantity: 10 g
	N719 Dye SKU: M2003 CAS Number: 207347-46-4 Options: • Quantity: 500 mg, 1 g		Silver Nanopowder SKU: M2444 CAS Number: 7440-22-4 Options: • Average particle size: 20 nm • Quantity: 5 g, 10 g, 25 g
	C106 Dye SKU: M2004 CAS Number: 1152310-69-4 Options: • Quantity: 250 mg		Spiro-OMeTAD (Spiro-MeOTAD) SKU: M291 CAS Number: 207739-72-8 Options: • Grade: Sublimed (>99.5% purity), Unsublimed (>99.0% purity) • Quantity: 500 mg, 1 g, 5 g, 10 g
	K19 Dye SKU: M2005 CAS Number: 847665-45-6 Options: • Quantity: 250 mg		
	N3 Dye SKU: M2006 CAS Number: 141460-19-7 Options: • Quantity: 1 g		
	FK102-Co(II)PF6 Salt SKU: M2007 CAS Number: 1392221-69-0 Options: • Quantity: 10 g		
	FK102-Co(II)TFSI Salt SKU: M2019 Options: • Quantity: 10 g		

OLED Materials

Organic light-emitting diodes (OLEDs) consist of layers of organic materials with charge transport or emissive properties sandwiched between conductive electrodes. OLEDs emit light at a specific color and wavelength, depending on the combination of organic material components. By selecting appropriate host and dopant OLED materials, the emission can be tailored to achieve desired results. We offer high-purity OLED materials, including fluorescent, phosphorescent, and TADF host materials, as well as a range of compatible dopant materials and charge transport layer materials, such as electron- and hole-injecting materials.




Popular Products



TAPC

High-purity, high-hole mobility, transport layer and host material.


- Exciplex host
- HOMO = 5.5 eV
- LUMO = 2.0 eV



Alq3

High thermal stability with high quantum yield of fluorescence

- High electron-transport
- HOMO = 5.62 eV
- LUMO = 2.85 eV



4CzIPN

High PLQY and metal-free organophotocatalyst

- Donor-acceptor fluorophore
- Fully substituted benzene ring

To help you choose the right material, full technical specifications are available on our website.



Ossila

M291-500mg

Spiro-OMeTAD

2,2',7,7'-Tetrakis(N,N'-di(4-methoxyphenylamino)-9,9'-spirobifluorene)

















CAS# 207739-72-8

Solpro Business Park

















Purity: >99.0%















MW: 424.4















Charge Transport Layer Materials

















	Bathocuproine (BCP) SKU: B231 CAS Number: 4733-39-5 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 1 g, 5 g		2,4,6-tris(biphenyl-3-yl)-1,3,5-triazine (T2T) SKU: M2103 CAS Number: 1201800-83-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	TPD - N,N'-Bis(3-methylphenyl)-N,N'-diphenylbenzidine SKU: M0401 CAS Number: 65181-78-4 Options: • Quantity: 1 g, 5 g		CzSi SKU: M2104 CAS Number: 898546-82-2 Options: • Quantity: 100 mg, 250 mg, 500 mg
	PTAA for Perovskite Applications SKU: M0511 CAS Number: 1333317-99-9 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DPEPO SKU: M2105 CAS Number: 808142-23-6 Options: • Quantity: 250 mg, 500 mg, 1 g
	TAPC SKU: M0811 CAS Number: 58473-78-2 Options: • Quantity: 250 mg, 500 mg, 1 g		B3PyPB SKU: M2127 CAS Number: 1030380-38-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	3TPYMB SKU: M2088 CAS Number: 929203-02-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		TTPA SKU: M2129 CAS Number: 177799-16-5 Options: • Quantity: 250 mg, 500 mg, 1 g
	ZADN SKU: M2091 CAS Number: 561064-11-7 Options: • Quantity: 500 mg, 1 g		DPPS SKU: M2130 CAS Number: 1152162-74-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	MADN SKU: M2092 CAS Number: 804560-00-7 Options: • Quantity: 1 g, 5 g		UGH-2 SKU: M2131 CAS Number: 40491-34-7 / 18856-08-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TmPyPB SKU: M2093 CAS Number: 921205-03-0 Options: • Quantity: 500 mg, 1 g		PPT, PO15 SKU: M2132 CAS Number: 1019842-99-9 Options: • Quantity: 250 mg, 500 mg, 1 g










* Some options have been omitted, please see our website to browse our full range

	B2PymPm SKU: M2173 CAS Number: 1266181-51-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Tris-PCz SKU: M2191 CAS Number: 1141757-83-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	B3PymPm SKU: M2174 CAS Number: 925425-96-3 Options: • Purity: Sublimed >99.0% purity, Unsublimed >98.0% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g		Spiro-TTB SKU: M2192 CAS Number: 515834-67-0 Options: • Quantity: 500 mg, 1 g, 5 g
	B4PymPm SKU: M2175 CAS Number: 1030380-51-8 Options: • Quantity: 250 mg, 500 mg, 1 g		DNTPD SKU: M2193 CAS Number: 199121-98-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	B4PyPPm SKU: M2176 CAS Number: 1097652-83-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TSP01 SKU: M2199 CAS Number: 1286708-86-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BTB SKU: M2179 CAS Number: 266349-83-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PO-T2T SKU: M2266 CAS Number: 1646906-26-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g, 25 g
	NPNPB SKU: M2182 CAS Number: 936355-01-0 Options: • Quantity: 500 mg, 1 g		DPPyA SKU: M2356 CAS Number: 1257879-34-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	NPB-DPA SKU: M2185 CAS Number: 910058-11-6 Options: • Grade: Sublimed, Unsublimed • Quantity: 500 mg, 1 g, 5 g		α-NPD SKU: M2411 CAS Number: 495416-60-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	PPDN SKU: M2187 CAS Number: 215611-93-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		2PACz SKU: M2425 CAS Number: 20999-38-6 Options: • Quantity: 500 mg, 1 g

















	Br-2PACz SKU: M2426 CAS Number: 2762888-11-7 Options: • Quantity: 500 mg, 1 g		MPA-CPA SKU: M2447 CAS Number: 2212003-31-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	MeO-2PACz SKU: M2427 CAS Number: 2377770-18-6 Options: • Quantity: 500 mg, 1 g		Cl-2PACz SKU: M2451 CAS Number: 3036926-72-1 Options: • Quantity: 500 mg, 1 g, 5 g
	2,3,6,7,10,11-Hexakis(hexyloxy)triphenylene - HAT6 SKU: M2433 CAS Number: 70351-86-9 Options: • Quantity: 250 mg, 500 mg, 1 g		4PACz SKU: M2452 CAS Number: 20999-36-4 Options: • Quantity: 500 mg, 1 g, 5 g
	PCzAc SKU: M2434 CAS Number: 1705584-08-2 Options: • Quantity: 250 mg, 500 mg, 1 g		MeO-4PACz SKU: M2453 CAS Number: 2922526-56-3 Options: • Quantity: 500 mg, 1 g, 5 g
	BPBPA SKU: M2436 CAS Number: 164724-35-0 Options: • Quantity: 500 mg, 1 g, 5 g		Me-4PACz SKU: M2454 CAS Number: 2747959-96-0 Options: • Quantity: 500 mg, 1 g, 5 g
	Nickel (II) Oxide Powder SKU: M2441 CAS Number: 1313-99-1 Options: • Particle size: 8 - 18 nm • Quantity: 10 g, 25 g		Br-2EPT SKU: M2455 CAS Number: 2826271-17-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	DMACPA SKU: M2445 CAS Number: 2971088-37-4 Options: • Quantity: 250 mg, 500 mg, 1 g		Silver Nanowire Dispersion SKU: M2456 CAS Number: 7440-22-4 Options: • Quantity: 25 mL, 50 mL
	4PADCB SKU: M2446 CAS Number: 2882156-63-8 Options: • Quantity: 250 mg, 500 mg, 1 g		TDPA-Cl SKU: M2457 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g

















	2BrPXZPA SKU: M2458 Options: • Quantity: 250 mg, 500 mg, 1 g		Me-PhpPACz SKU: M2467 CAS Number: 3046309-16-1 Options: • Quantity: 250 mg, 500 mg, 1 g
	2BrPTZPA SKU: M2459 CAS Number: 2913184-53-7 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		Benzyl Viologen Dichloride SKU: M2469 CAS Number: 1102-19-8 Options: • Quantity: 5 g, 10 g, 25 g
	3PATAT-C3 SKU: M2460 CAS Number: 3048690-10-1 Options: • Quantity: 250 mg, 500 mg, 1 g		PTZ-CPA SKU: M2474 Options: • Quantity: 250 mg, 500 mg, 1 g
	9CAA SKU: M2461 CAS Number: 524-80-1 Options: • Quantity: 500 mg, 1 g, 5 g		Py3 SKU: M2478 CAS Number: 1520027-62-6 Options: • Quantity: 250 mg, 500 mg, 1 g
	DCB-BPA SKU: M2462 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		Ph-4PACz SKU: M2481 CAS Number: 2814500-04-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	9CPA SKU: M2463 CAS Number: 6622-54-4 Options: • Quantity: 500 mg, 1 g, 5 g		Tin (IV) Oxide Powder SKU: M2484 CAS Number: 18282-10-5 Options: • Quantity: 50 g, 100 g, 250 g
	3-BPIC-F SKU: M2464 CAS Number: 3038707-50-2 Options: • Quantity: 250 mg, 500 mg, 1 g		Tungsten Oxide Powder SKU: M2485 CAS Number: 1314-35-8 Options: • Quantity: 100 g, 250 g, 500 g
	3-BPIC SKU: M2465 CAS Number: 3038707-49-9 Options: • Quantity: 250 mg, 500 mg, 1 g		Spiro-OMeTAD (Spiro-MeOTAD) SKU: M291 CAS Number: 207739-72-8 Options: • Grade: Sublimed (>99.5% purity), Unsublimed (>99.0% purity) • Quantity: 500 mg, 1 g, 5 g, 10 g
















	CuPc, Copper(II) phthalocyanine SKU: M341 CAS Number: 147-14-8 Options: Quantity: 1 g, 5 g		OXD-7 SKU: M451 CAS Number: 138372-67-5 Options: Quantity: 250 mg, 500 mg, 1 g
	F4TCNQ SKU: M351 CAS Number: 29261-33-4 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g		TCTA SKU: M471 CAS Number: 139092-78-7 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) - Quantity: 500 mg, 1 g
	NPB (NPD) SKU: M361 CAS Number: 123847-85-8 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) - Quantity: 1 g, 5 g		Pentacene SKU: M581 CAS Number: 135-48-8 Options: Quantity: 1 g, 5 g
	mCP SKU: M371 CAS Number: 550378-78-4 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) - Quantity: 1 g, 5 g		TAZ, 3-(Biphenyl-4-yl)-5-(4-tert-butylphenyl)-4-phenyl-4H-1,2,4-triazole SKU: M601 CAS Number: 150405-69-9 Options: Quantity: 250 mg, 500 mg, 1 g
	Alq3 SKU: M381 CAS Number: 2085-33-8 Options: Quantity: 5 g, 10 g		2-TNATA SKU: M611 CAS Number: 185690-41-9 Options: - Grade: Sublimed (>99.0% purity), Unsublimed (>98.0% purity) - Quantity: 500 mg, 1 g, 5 g
	Zn(BTZ)2 SKU: M411 CAS Number: 58280-31-2 Options: Quantity: 1 g		m-MTDATA SKU: M621 CAS Number: 124729-98-2 Options: Quantity: 250 mg, 500 mg, 1 g
	BAlq SKU: M421 CAS Number: 146162-54-1 Options: Quantity: 1 g		Cobalt(II) Phthalocyanine, CoPc SKU: M631 CAS Number: 3317-67-7 Options: Quantity: 500 mg, 1 g
	DPVBi SKU: M431 CAS Number: 142289-08-5 Options: Quantity: 1 g		TPBi SKU: M651 CAS Number: 192198-85-9 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) - Quantity: 250 mg, 500 mg, 1 g

















	TCNQ SKU: M681 CAS Number: 1518-16-7 Options: Quantity: 1 g, 5 g		BPhen (Bathophenanthroline) SKU: M961 CAS Number: 1662-01-7 Options: Quantity: 1 g, 5 g
	Bepp2 SKU: M691 CAS Number: 220694-90-6 Options: Quantity: 1 g		
	Titanyl phthalocyanine, TiOPc SKU: M701 CAS Number: 26201-32-1 Options: - Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) - Quantity: 500 mg, 1 g, 5 g		
	Ir(ppz)3 SKU: M721 CAS Number: 562824-31-1 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	8-Quinolinolato lithium (Liq) SKU: M731 CAS Number: 25387-93-3 Options: - Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) - Quantity: 1 g, 5 g		
	Bebq2 SKU: M741 CAS Number: 148896-39-3 Options: Quantity: 1 g		
	tCP, 1,3,5-Tris(carbazol-9-yl)benzene SKU: M751 CAS Number: 148044-07-9 Options: - Grade: Sublimed (>99.0% purity), Unsublimed (>98.0% purity) - Quantity: 1 g, 5 g		
	HATCN SKU: M791 CAS Number: 105598-27-4 Options: - Grade: Sublimed (>99.0%), Unsublimed (>98.0%) - Quantity: 500 mg, 1 g, 5 g, 10 g		


Host Materials

	TPD - N,N'-Bis(3-methylphenyl)-N,N'-diphenylbenzidine SKU: M0401 CAS Number: 65181-78-4 Options: • Quantity: 1 g, 5 g		CzSi SKU: M2104 CAS Number: 898546-82-2 Options: • Quantity: 100 mg, 250 mg, 500 mg
	26DCzPPy SKU: M0771 CAS Number: 1013405-24-7 Options: • Grade: Sublimed Grade (>99% Purity), Unsublimed Grade (>99% Purity) • Quantity: 250 mg, 500 mg, 1 g		DPEPO SKU: M2105 CAS Number: 808142-23-6 Options: • Quantity: 250 mg, 500 mg, 1 g
	TAPC SKU: M0811 CAS Number: 58473-78-2 Options: • Quantity: 250 mg, 500 mg, 1 g		BCzPh SKU: M2119 CAS Number: 57102-62-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	3TPYMB SKU: M0771 CAS Number: 929203-02-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		mCPSOB SKU: M2105 CAS Number: 1374770-41-8 Options: • Quantity: 250 mg, 500 mg, 1 g
	MADN SKU: M2092 CAS Number: 804560-00-7 Options: • Quantity: 1 g, 5 g		PXZ-DPS SKU: M2122 CAS Number: 1477511-57-1 Options: • Quantity: 100 mg, 250 mg, 500 mg
	TmPyPB SKU: M2093 CAS Number: 921205-03-0 Options: • Quantity: 500 mg, 1 g		UGH-2 SKU: M2131 CAS Number: 40491-34-7 / 18856-08-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PYD-2Cz (PYD2) SKU: M2102 CAS Number: 168127-49-9 Options: • Quantity: 250 mg, 500 mg, 1 g		PPT, PO15 SKU: M2132 CAS Number: 1019842-99-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	2,4,6-tris(biphenyl-3-yl)-1,3,5-triazine (T2T) SKU: M2103 CAS Number: 1201800-83-0 Options: • Quantity: 250 mg, 500 mg, 1 g		BTB SKU: M2179 CAS Number: 266349-83-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g


	mCBP SKU: M2186 CAS Number: 342638-54-4 Options: • Grade: Sublimed (>99%), Unsublimed (>98%) • Quantity: 500 mg, 1 g, 5 g		PXZ-TRZ SKU: M2268 CAS Number: 1411910-25-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	CDBP SKU: M2188 CAS Number: 120260-01-7 Options: • Quantity: 250 mg, 500 mg, 1 g		UGH3 SKU: M2270 CAS Number: 18920-16-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Tris-PCz SKU: M2191 CAS Number: 1141757-83-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCPCN SKU: M2272 CAS Number: 1392506-99-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BCBP SKU: M2194 CAS Number: 858131-70-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		3P-T2T SKU: M2273 CAS Number: 352196-01-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BCPO SKU: M2198 CAS Number: 1233407-28-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		3N-T2T SKU: M2275 CAS Number: 939430-26-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TSP01 SKU: M2199 CAS Number: 1286708-86-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		CbBPCb SKU: M2289 CAS Number: 1469997-91-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMAC-BPP SKU: M2247 CAS Number: 1836192-40-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		FATPA SKU: M2290 CAS Number: 1131007-94-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PO-T2T SKU: M2266 CAS Number: 1646906-24-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g, 25 g		SimCP2 SKU: M2291 CAS Number: 944465-42-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

	BSB SKU: M2293 CAS Number: 18826-13-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BCzTPA SKU: M2334 CAS Number: 1032174-52-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	SimCP SKU: M2294 CAS Number: 850221-63-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCPPO1 SKU: M2335 CAS Number: 1296229-26-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PPF SKU: M2295 CAS Number: 911397-27-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		POSTF SKU: M2337 CAS Number: 1647050-25-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mSiTrz SKU: M2309 CAS Number: 2252416-90-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCBP-CN SKU: M2335 CAS Number: 1327163-09-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPTPCz SKU: M2313 CAS Number: 1313391-57-9 Options: • Quantity: 250 mg, 500 mg, 1 g		mCPBC SKU: M2353 CAS Number: 1799525-58-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TCPZ SKU: M2314 CAS Number: 890148-68-2 Options: • Quantity: 250 mg, 500 mg, 1 g		46DCzPPm (CzPhPy) SKU: M2354 CAS Number: 1262678-77-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TBA DN SKU: M2314 CAS Number: 274905-73-6 Options: • Quantity: 250 mg, 500 mg, 1 g		PPO27 SKU: M2376 CAS Number: 1299463-56-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	SPPO13 SKU: M2333 CAS Number: 1234510-13-4 Options: • Quantity: 250 mg, 500 mg, 1 g		SiCzCz SKU: M2376 CAS Number: 2665768-26-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g


	SiTrzCz2 SKU: M2377 CAS Number: 2422045-58-5 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g		CzPA SKU: M2506 CAS Number: 910647-23-3 Options: Quantity: 250 mg, 500 mg, 1 g
	Poly(9-vinylcarbazole) (PVK) SKU: M2384 CAS Number: 25067-59-8 Options: - Batch specification: M2384A6 (Mw = 235700; PDI = 2.22) - Quantity: 5 g, 10 g		DIC-TRZ SKU: M2507 CAS Number: 1024598-01-3 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PCzAc SKU: M2434 CAS Number: 1705584-08-2 Options: Quantity: 250 mg, 500 mg, 1 g		DBF-TRZ SKU: M2520 CAS Number: 2040300-54-5 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PhCzBCz SKU: M2479 CAS Number: 2692683-03-5 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCP SKU: M371 CAS Number: 550378-78-4 Options: - Grade: Sublimed Grade ($>99.5\%$ purity), Unsublimed Grade ($>98.0\%$ purity) - Quantity: 1 g, 5 g
	SF3-TRZ SKU: M2480 CAS Number: 1639424-37-5 Options: - Purity: Sublimed $> 99\%$, Unsublimed $>98\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g		Alq3 SKU: M381 CAS Number: 2085-33-8 Options: Quantity: 5 g, 10 g
	mMDBA-DI SKU: M2496 CAS Number: 2770257-85-5 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g		CBP, 4,4'-Bis(N-carbazolyl)-1,1'-biphenyl SKU: M391 CAS Number: 58328-31-7 Options: - Grade: Sublimed Grade ($>99.5\%$ purity) - Quantity: 1 g, 5 g
	3CzBN SKU: M2502 CAS Number: 1799923-03-7 Options: - Grade: Sublimed, Unsublimed - Quantity: 250 mg, 500 mg, 1 g		Zn(BTZ)2 SKU: M411 CAS Number: 58280-31-2 Options: Quantity: 1 g
	DMIC-TRZ SKU: M2504 CAS Number: 1257248-13-7 Options: Quantity: 250 mg, 500 mg, 1 g		BAiq SKU: M421 CAS Number: 146162-54-1 Options: Quantity: 1 g




DPVBi
SKU: M431 | CAS Number: 142289-08-5
Options:
• Quantity: 1 g




35DCzPPy
SKU: M781 | CAS Number: 1013405-25-8
Options:
• Quantity: 1 g




OXD-7
SKU: M451 | CAS Number: 138372-67-5
Options:
• Quantity: 250 mg, 500 mg, 1 g




ADN - 9,10-Bis(2-naphthyl)anthracene
SKU: M461 | CAS Number: 122648-99-1
Options:
• Grade: Sublimed (>99.0% purity),
Unsublimed (>98.0% purity)
• Quantity: 1 g, 5 g




TCTA
SKU: M471 | CAS Number: 139092-78-7
Options:
• Grade: Sublimed Grade (>99.5% purity),
Unsublimed Grade (>98.0% purity)
• Quantity: 500 mg, 1 g




TAZ, 3-(Biphenyl-4-yl)-5-(4-tert-butylphenyl)-4-phenyl-4H-1,2,4-triazole
SKU: M601 | CAS Number: 150405-69-9
Options:
• Quantity: 250 mg, 500 mg, 1 g



m-MTDATA
SKU: M621 | CAS Number: 124729-98-2
Options:
• Quantity: 250 mg, 500 mg, 1 g


































TPBi
SKU: M651 | CAS Number: 192198-85-9
Options:
• Grade: Sublimed Grade (>99.5% purity),
Unsublimed Grade (>98.0% purity)
• Quantity: 250 mg, 500 mg, 1 g



















tCP, 1,3,5-Tris(carbazol-9-yl)benzene
SKU: M751 | CAS Number: 148044-07-9
Options:
• Grade: Sublimed (>99.0% purity),
Unsublimed (>98.0% purity)
• Quantity: 1 g, 5 g


Dopant Materials

	FirPic SKU: M0711 CAS Number: 376367-93-0 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DMAC-DPS SKU: M2121 CAS Number: 1477512-32-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMAC-BP SKU: M2090 CAS Number: 1685287-55-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		PXZ-DPS SKU: M2122 CAS Number: 1477511-57-1 Options: • Quantity: 100 mg, 250 mg, 500 mg
	MADN SKU: M2092 CAS Number: 804560-00-7 Options: • Quantity: 1 g, 5 g		Cab-Ph-TRZ SKU: M2123 CAS Number: 440354-93-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	C545T, Coumarin 545T SKU: M2094 CAS Number: 155306-71-1 Options: • Quantity: 500 mg, 1 g		TBRb SKU: M2124 CAS Number: 682806-51-5 Options: • Quantity: 250 mg, 500 mg, 1 g
	4CzIPN SKU: M2100 CAS Number: 1416881-52-1 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Fir6 SKU: M2126 CAS Number: 664374-03-2 Options: • Quantity: 100 mg, 250 mg
	DBP SKU: M2101 CAS Number: 175606-05-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TXO-TPA SKU: M2128 CAS Number: 1623010-63-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TBPe SKU: M2117 CAS Number: 80663-92-9 Options: • Quantity: 250 mg, 500 mg, 1 g		TTPA SKU: M2129 CAS Number: 177799-16-5 Options: • Quantity: 250 mg, 500 mg, 1 g
	ACSRA SKU: M2118 CAS Number: 1206626-95-0 Options: • Quantity: 1 g, 250 mg, 500 mg		DCJTb SKU: M2177 CAS Number: 200052-70-6 Options: • Quantity: 100 mg, 250 mg, 500 mg


	2CzPN SKU: M2181 CAS Number: 1416881-50-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		5CzBN SKU: M2219 CAS Number: 1469700-24-0 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(MDQ)2(acac) SKU: M2184 CAS Number: 536755-34-7 Options: • Quantity: 100 mg, 250 mg		5TCzBN SKU: M2220 CAS Number: 1469700-26-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPAVBi SKU: M2189 CAS Number: 119586-44-6 Options: • Quantity: 500 mg, 1 g		4CzPN SKU: M2221 CAS Number: 1416881-51-0 Options: • Grade: Sublimed > 99% purity, Unsublimed > 98% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(mppy)3 SKU: M2195 CAS Number: 149005-33-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BD-1 SKU: M2222 CAS Number: 1252914-27-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzTPN SKU: M2214 CAS Number: 1416881-53-2 Options: • Quantity: 250 mg, 500 mg, 1 g		Ir(dpm)PQ2 SKU: M2231 CAS Number: 713079-03-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzTPN-Ph SKU: M2216 CAS Number: 1416881-55-4 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 250 mg, 500 mg, 1 g		DMAC-BPP SKU: M2247 CAS Number: 1836192-40-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzBN SKU: M2214 CAS Number: 1996609-93-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzTPN-tBu SKU: M2251 CAS Number: 2153433-46-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	4TCzBN SKU: M2218 CAS Number: 2055722-93-3 Options: • Quantity: 100 mg, 250 mg, 500 mg		Phen-TRZ SKU: M2264 CAS Number: 1357066-21-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

	TXO-PhCz SKU: M2265 CAS Number: 1623010-64-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TCzTRZ SKU: M2332 CAS Number: 1808158-40-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PXZ-TRZ SKU: M2268 CAS Number: 1411910-25-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TmCzTRZ SKU: M2336 CAS Number: 1808158-41-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DDCzTRZ SKU: M2286 CAS Number: 1685282-47-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DCzTRZ SKU: M2338 CAS Number: 1106730-48-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Eu(dbm)3(Phen) SKU: M2296 CAS Number: 17904-83-5 Options: • Quantity: 250 mg, 500 mg, 1 g		DMAC-TRZ SKU: M2343 CAS Number: 1628752-98-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	FK306 SKU: M2308 CAS Number: 1421058-47-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DMPPP SKU: M2344 CAS Number: 1036404-84-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(dpm)(piq)2 SKU: M2310 CAS Number: 1202867-58-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzPN-Me SKU: M2378 CAS Number: 1469700-28-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	Hex-Ir(phq)3 SKU: M2312 CAS Number: 1268460-37-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzPN-tBu SKU: M2379 CAS Number: 1469705-93-8 Options: • Grade: Sublimed > 99% purity, Unsublimed > 98% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzFCN SKU: M2331 CAS Number: 1819362-10-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzIPN-Me SKU: M2380 CAS Number: 1469703-61-4 Options: • Quantity: 250 mg, 500 mg, 1 g


* Some options have been omitted, please see our website to browse our full range




4CzIPN-tBu
SKU: M2381 | CAS Number: 1630263-99-8
Options:
• Quantity: 250 mg, 500 mg, 1 g




TPA-DCPP
SKU: M2418 | CAS Number: 1803287-94-6
Options:
• Quantity: 250 mg, 500 mg, 1 g




4CzIPN-Ph
SKU: M2382 | CAS Number: 1469705-37-0
Options:
• Quantity: 250 mg, 500 mg, 1 g




TmPPPyTz
SKU: M2419 | CAS Number: 939430-31-6
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g



Cz-PS
SKU: M2386 | CAS Number: 1396165-20-0
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g



t-DABNA
SKU: M2420 | CAS Number: 1802003-07-1
Options:
• Grade: Sublimed > 99.5%, Unsublimed >99%
• Quantity: 100 mg, 250 mg, 500 mg, 1 g




AQ(PhDPA)₂
SKU: M2387 | CAS Number: 1640978-33-1
Options:
• Purity grade: Sublimed 99%, Unsublimed 98%
• Quantity: 100 mg, 250 mg, 500 mg, 1 g



v-DABNA
SKU: M2423 | CAS Number: 2251782-13-3
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g




Px-VPN
SKU: M2388 | CAS Number: 1784766-39-7
Options:
• Quantity: 100 mg, 250 mg, 500 mg




DABNA-1
SKU: M2429 | CAS Number: 1689552-89-3
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g




Ac-CNP
SKU: M2389 | CAS Number: 1883400-34-7
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g




CzBN
SKU: M2430 | CAS Number: 2170487-27-9
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g



DACT-II
SKU: M2390 | CAS Number: 1613074-59-1
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g



















3tPAB
SKU: M2431 | CAS Number: 2056877-25-7
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g











CzDBA
SKU: M2417 | CAS Number: 2171334-43-1
Options:
• Quantity: 100 mg, 250 mg


































t-Bu-v-DABNA
SKU: M2432 | CAS Number: 2914070-28-1
Options:
• Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g
















	SpiroAC-TRZ SKU: M2435 CAS Number: 1980037-96-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		2PTZBN SKU: M2512 CAS Number: 2250280-17-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DtBuCzB SKU: M2493 CAS Number: 2170487-30-4 Options: • Purity: Sublimed > 99%, Unsublimed > 97% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DiKta SKU: M2513 CAS Number: 32081-22-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mMDBA-DI SKU: M2496 CAS Number: 2770257-85-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		3DPA-DiKta SKU: M2514 Options: • Quantity: 250 mg, 500 mg, 1 g
	PPCzTRz SKU: M2499 CAS Number: 2329651-89-8 Options: • Grade: Sublimed ≥99.0%, Unsublimed ≥98.0% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TB-tPCz SKU: M2515 Options: • Quantity: 250 mg, 500 mg, 1 g
	5Cz-TRZ SKU: M2500 CAS Number: 2043597-42-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		tDPA-DtCzB SKU: M2540 CAS Number: 2760492-67-7 Options: • Grade: Sublimed ≥99% • Quantity: 250 mg, 500 mg, 1 g
	3CzBN SKU: M2502 CAS Number: 1799923-03-7 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g		Alq3 SKU: M381 CAS Number: 2085-33-8 Options: • Quantity: 5 g, 10 g
	TDBA-SAF SKU: M2509 CAS Number: 2643306-28-7 Options: • Quantity: 250 mg, 500 mg, 1 g		Rubrene SKU: M441 CAS Number: 517-51-1 Options: • Quantity: 500 mg, 1 g
	TB-P3Cz SKU: M2511 CAS Number: 2411720-50-6 Options: • Quantity: 250 mg, 500 mg, 1 g		Ir(ppy)3 SKU: M481 CAS Number: 94928-86-6 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed (>98.0% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g

















	Ir(btpy)3 SKU: M491 CAS Number: 405289-74-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(PIQ)3 SKU: M641 CAS Number: 435293-93-9 Options: • Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(ppy)2(acac) SKU: M661 CAS Number: 337526-85-9 Options: • Grade: Sublimed (>99.5% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(PIQ)2(acac) SKU: M671 CAS Number: 435294-03-4 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>99.0% purity) • Quantity: 250 mg, 500 mg, 1 g
	Ir(ppz)3 SKU: M721 CAS Number: 562824-31-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(ppy)2(bmp) SKU: M761 CAS Number: 1215692-34-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(dmpq)2(acac) SKU: M801 CAS Number: 1056874-46-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMQA SKU: M971 CAS Number: 19205-19-7 Options: • Quantity: 1 g, 5 g
















TADF Materials















	26DCzPPy SKU: M0771 CAS Number: 1013405-24-7 Options: • Grade: Sublimed Grade (>99% Purity), Unsublimed Grade (>99% Purity) • Quantity: 250 mg, 500 mg, 1 g		2,4,6-tris(biphenyl-3-yl)-1,3,5-triazine (T2T) SKU: M2103 CAS Number: 1201800-83-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	3TPYMB SKU: M2088 CAS Number: 929203-02-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		CzSi SKU: M2104 CAS Number: 898546-82-2 Options: • Quantity: 100 mg, 250 mg, 500 mg
	DMAC-BP SKU: M2090 CAS Number: 1685287-55-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		DPEPO SKU: M2105 CAS Number: 808142-23-6 Options: • Quantity: 250 mg, 500 mg, 1 g
	ZADN SKU: M2091 CAS Number: 561064-11-7 Options: • Quantity: 500 mg, 1 g		TBPe SKU: M2117 CAS Number: 80663-92-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	TmPyPB SKU: M2093 CAS Number: 921205-03-0 Options: • Quantity: 500 mg, 1 g		ACRSA SKU: M2118 CAS Number: 1206626-95-0 Options: • Quantity: 1 g, 250 mg, 500 mg
	4CzIPN SKU: M2100 CAS Number: 1416881-52-1 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BCzPh SKU: M2119 CAS Number: 57102-62-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	DBP SKU: M2101 CAS Number: 175606-05-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCPSOB SKU: M2120 CAS Number: 1374770-41-8 Options: • Quantity: 250 mg, 500 mg, 1 g
	PYD-2Cz (PYD2) SKU: M2102 CAS Number: 168127-49-9 Options: • Quantity: 250 mg, 500 mg, 1 g		DMAC-DPS SKU: M2121 CAS Number: 1477512-32-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

	PXZ-DPS SKU: M2122 CAS Number: 1477511-57-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		PPT, PO15 SKU: M2132 CAS Number: 1019842-99-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	Cab-Ph-TRZ SKU: M2123 CAS Number: 440354-93-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		B2PymPm SKU: M2173 CAS Number: 1266181-51-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TBRb SKU: M2125 CAS Number: 682806-51-5 Options: • Quantity: 250 mg, 500 mg, 1 g		B4PymPm SKU: M2175 CAS Number: 1030380-51-8 Options: • Quantity: 250 mg, 500 mg, 1 g
	B3PyPB SKU: M2127 CAS Number: 1030380-38-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		B4PyPPm SKU: M2173 CAS Number: 1097652-83-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TXO-TPA SKU: M2128 CAS Number: 1623010-63-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DCJTb SKU: M2177 CAS Number: 200052-70-6 Options: • Quantity: 100 mg, 250 mg, 500 mg
	TTPA SKU: M2129 CAS Number: 177799-16-5 Options: • Quantity: 250 mg, 500 mg, 1 g		BP4mPy SKU: M2178 CAS Number: 1009033-94-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPPS SKU: M2130 CAS Number: 1152162-74-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BTB SKU: M2179 CAS Number: 266349-83-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	UGH-2 SKU: M2131 CAS Number: 40491-34-7 / 18856-08-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		2CzPN SKU: M2181 CAS Number: 1416881-50-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g















	NBPhen SKU: M2183 CAS Number: 1174006-43-9 Options: • Quantity: 250 mg, 500 mg, 1 g		TSP01 SKU: M2199 CAS Number: 1286708-86-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	NPB-DPA SKU: M2185 CAS Number: 910058-11-6 Options: • Grade: Sublimed, Unsublimed • Quantity: 500 mg, 1 g, 5 g		TSBF SKU: M2213 CAS Number: 518997-91-6 Options: • Quantity: 500 mg, 1 g
	mCBP SKU: M2186 CAS Number: 342638-54-4 Options: • Grade: Sublimed (>99%), Unsublimed (>98%) • Quantity: 500 mg, 1 g, 5 g		4CzTPN SKU: M2214 CAS Number: 1416881-53-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	PPDN SKU: M2187 CAS Number: 215611-93-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzTPN-Me SKU: M2215 CAS Number: 1416881-54-3 Options: • Purity: Sublimed grade A1 (> 99.0%), Sublimed grade A2 (> 99.0%), Unsublimed grade (> 98%) • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	CDBP SKU: M2188 CAS Number: 120260-01-7 Options: • Quantity: 250 mg, 500 mg, 1 g		4CzTPN-Ph SKU: M2216 CAS Number: 1416881-55-4 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 250 mg, 500 mg, 1 g
	Tris-PCz SKU: M2191 CAS Number: 1141757-28-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzBN SKU: M2217 CAS Number: 1996609-93-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BCBP SKU: M2194 CAS Number: 858131-70-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4TCzBN SKU: M2218 CAS Number: 2055722-93-3 Options: • Quantity: 100 mg, 250 mg, 500 mg
	BCPO SKU: M2198 CAS Number: 1233407-28-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		5CzBN SKU: M2219 CAS Number: 1469700-24-0 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g


	STCzBN SKU: M2220 CAS Number: 1469700-26-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		UGH3 SKU: M2270 CAS Number: 18920-16-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzPN SKU: M2221 CAS Number: 1416881-51-0 Options: • Grade: Sublimed > 99% purity, Unsublimed > 98% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g		CzAcSF SKU: M2271 CAS Number: 1792173-34-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMAC-BPP SKU: M2247 CAS Number: 1836192-40-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCPCN SKU: M2272 CAS Number: 1392506-99-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzTPN-tBu SKU: M2251 CAS Number: 2153433-46-4 Options: • Quantity: 250 mg, 500 mg, 1 g		3P-T2T SKU: M2273 CAS Number: 352196-01-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TXO-PhCz SKU: M2265 CAS Number: 1623010-64-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		3N-T2T SKU: M2275 CAS Number: 939430-26-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PO-T2T SKU: M2266 CAS Number: 1646906-26-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g, 25 g		DDCzTRZ SKU: M2286 CAS Number: 1685282-47-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DCzDBT SKU: M2267 CAS Number: 913738-04-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		CbBPCb SKU: M2289 CAS Number: 1469997-91-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PXZ-TRZ SKU: M2268 CAS Number: 1411910-25-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		FATPA SKU: M2290 CAS Number: 1131007-94-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

	SimCP2 SKU: M2291 CAS Number: 944465-42-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TCPZ SKU: M2314 CAS Number: 890148-68-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	SPPO1 SKU: M2292 CAS Number: 1125547-88-7 Options: • Quantity: 250 mg, 500 mg		4CzFCN SKU: M2331 CAS Number: 1819362-10-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BSB SKU: M2293 CAS Number: 18826-13-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TCzTRZ SKU: M2332 CAS Number: 1808158-40-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	SimCP SKU: M2294 CAS Number: 850221-63-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SPPO13 SKU: M2333 CAS Number: 1234510-13-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	mSiTrz SKU: M2309 CAS Number: 2252416-90-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BCzTPA SKU: M2334 CAS Number: 1032174-52-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	Ir(dpm)(piq)2 SKU: M2310 CAS Number: 1202867-58-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCPPO1 SKU: M2335 CAS Number: 1296229-26-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TSBPA SKU: M2309 CAS Number: 205327-13-5 Options: • Quantity: 250 mg, 500 mg, 1 g		POSTF SKU: M2337 CAS Number: 1647050-25-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPTCz SKU: M2313 CAS Number: 1313391-57-9 Options: • Quantity: 250 mg, 500 mg, 1 g		DCzTRZ SKU: M2338 CAS Number: 1106730-48-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g


	BPy-OXD SKU: M2339 CAS Number: 866117-19-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzPN-Ph SKU: M2375 CAS Number: 1469707-47-8 Options: • Quantity: 250 mg, 500 mg, 1 g
	DTAF SKU: M2342 CAS Number: 159526-57-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SiCzCz SKU: M2376 CAS Number: 2665768-26-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMAC-TRZ SKU: M2343 CAS Number: 1628752-98-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SiTrzCz2 SKU: M2377 CAS Number: 2422045-58-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMPPP SKU: M2344 CAS Number: 1036404-84-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzPN-Me SKU: M2378 CAS Number: 1469700-28-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	mCBP-CN SKU: M2345 CAS Number: 1327163-09-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzPN-tBu SKU: M2379 CAS Number: 1469705-93-8 Options: • Grade: Sublimed > 99% purity, Unsublimed > 98% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	46DCzPPm (CzPhPy) SKU: M2354 CAS Number: 1262678-77-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzIPN-Me SKU: M2380 CAS Number: 1469703-61-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	DPPyA SKU: M2356 CAS Number: 1257879-34-7 Options: • Quantity: 250 mg, 500 mg, 1 g		4CzIPN-tBu SKU: M2381 CAS Number: 1630263-99-8 Options: • Quantity: 250 mg, 500 mg, 1 g
	PPO27 SKU: M2357 CAS Number: 1299463-56-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzIPN-Ph SKU: M2382 CAS Number: 1469705-37-0 Options: • Quantity: 250 mg, 500 mg, 1 g

* Some options have been omitted, please see our website to browse our full range


	Cz-DPS SKU: M2385 CAS Number: 733038-89-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		v-DABNA SKU: M2423 CAS Number: 2251782-13-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Cz-PS SKU: M2386 CAS Number: 1396165-20-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DABNA-1 SKU: M2429 CAS Number: 1689552-89-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	AQ(PhDPA)₂ SKU: M2387 CAS Number: 1640978-33-1 Options: • Purity grade: Sublimed 99%, Unsublimed 98% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		CzBN SKU: M2430 CAS Number: 2170487-27-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Px-VPN SKU: M2389 CAS Number: 1784766-39-7 Options: • Quantity: 100 mg, 250 mg, 500 mg		3tPAB SKU: M2431 CAS Number: 2056877-25-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ac-CNP SKU: M2389 CAS Number: 1883400-34-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		t-Bu-v-DABNA SKU: M2432 CAS Number: 2914070-28-1 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g
	DACT-II SKU: M2390 CAS Number: 1613074-59-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PCzAc SKU: M2434 CAS Number: 1705584-08-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	CzDBA SKU: M2417 CAS Number: 2171334-43-1 Options: • Quantity: 100 mg, 250 mg		SpiroAC-TRZ SKU: M2493 CAS Number: 1980037-96-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	t-DABNA SKU: M2420 CAS Number: 1802003-07-1 Options: • Grade: Sublimed > 99.5%, Unsublimed > 99% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DtBuCzB SKU: M2493 CAS Number: 2170487-30-4 Options: • Purity: Sublimed > 99%, Unsublimed > 97% • Quantity: 100 mg, 250 mg, 500 mg, 1 g




3CzClIPN
SKU: M2498 | CAS Number: 1469704-61-7
Options:
• Quantity: 250 mg, 500 mg, 1 g, 5 g




4DPAIPN
SKU: M2508 | CAS Number: 1846598-27-3
Options:
• Grade: Sublimed, Unsublimed
• Quantity: 250 mg, 500 mg, 1 g




PPCzTRz
SKU: M2499 | CAS Number: 2329651-89-8
Options:
• Grade: Sublimed $\geq 99.0\%$, Unsublimed $\geq 98.0\%$
• Quantity: 100 mg, 250 mg, 500 mg, 1 g




TDBA-SAF
SKU: M2509 | CAS Number: 2643306-28-7
Options:
• Quantity: 250 mg, 500 mg, 1 g




5Cz-TRZ
SKU: M2500 | CAS Number: 2043597-42-6
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g




TB-P3Cz
SKU: M2511 | CAS Number: 2411720-50-6
Options:
• Quantity: 250 mg, 500 mg, 1 g




3DPAFIPN
SKU: M2501 | CAS Number: 2260543-73-3
Options:
• Grade: Sublimed, Unsublimed
• Quantity: 250 mg, 500 mg, 1 g




2PTZBN
SKU: M2512 | CAS Number: 2250280-17-0
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g




3CzBN
SKU: M2502 | CAS Number: 1799923-03-7
Options:
• Grade: Sublimed, Unsublimed
• Quantity: 250 mg, 500 mg, 1 g




DiKta
SKU: M2513 | CAS Number: 32081-22-4
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g




3DPA2FBN
SKU: M2503 | CAS Number: 1403850-00-9
Options:
• Grade: Sublimed, Unsublimed
• Quantity: 250 mg, 500 mg, 1 g




3DPA-DiKta
SKU: M2514
Options:
• Quantity: 250 mg, 500 mg, 1 g




DMIC-TRZ
SKU: M2504 | CAS Number: 1257248-13-7
Options:
• Quantity: 250 mg, 500 mg, 1 g



















TB-tPCz
SKU: M2515
Options:
• Quantity: 250 mg, 500 mg, 1 g


















DIC-TRZ
SKU: M2507 | CAS Number: 1024598-01-3
Options:
• Quantity: 100 mg, 250 mg, 500 mg, 1 g



































Mes3DiKta
SKU: M2516
Options:
• Quantity: 250 mg, 500 mg, 1 g















	3TPA-DiKta SKU: M2517 Options: Quantity: 250 mg, 500 mg, 1 g		p-bPPhenB SKU: M2541 CAS Number: 2135536-38-6 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	T3PyTRZ SKU: M2519 CAS Number: 2703850-25-1 Options: - Grade: Sublimed $\geq 99\%$, Unsublimed $\geq 98\%$ - Quantity: 250 mg, 500 mg, 1 g		SiDBFCz SKU: M2543 CAS Number: 1234996-95-2 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DBF-TRZ SKU: M2520 CAS Number: 2040300-54-5 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g		DCz-BTP SKU: M2545 CAS Number: 1235883-57-4 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DOBNA-OAr SKU: M2536 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 250 mg, 500 mg, 1 g		mCP SKU: M371 CAS Number: 550378-78-4 Options: - Grade: Sublimed Grade ($>99.5\%$ purity), Unsublimed Grade ($>98.0\%$ purity) - Quantity: 1 g, 5 g
	4DBFHPB SKU: M2537 CAS Number: 2280942-02-9 Options: - Grade: Sublimed $\geq 99\%$, Unsublimed $\geq 98\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g		CBP, 4,4'-Bis(N-carbazolyl)-1,1'-biphenyl SKU: M391 CAS Number: 58328-31-7 Options: - Grade: Sublimed Grade ($>99.5\%$ purity) - Quantity: 1 g, 5 g
	BN-CP1 SKU: M2538 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 250 mg, 500 mg, 1 g		TCTA SKU: M471 CAS Number: 139092-78-7 Options: - Grade: Sublimed Grade ($>99.5\%$ purity), Unsublimed Grade ($>98.0\%$ purity) - Quantity: 500 mg, 1 g
	BPBiPA SKU: M2539 CAS Number: 1256723-60-0 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 500 mg, 1 g		TPBi SKU: M651 CAS Number: 192198-85-9 Options: - Grade: Sublimed Grade ($>99.5\%$ purity), Unsublimed Grade ($>98.0\%$ purity) - Quantity: 250 mg, 500 mg, 1 g
	tDPA-DtCzB SKU: M2540 CAS Number: 2760492-67-7 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 250 mg, 500 mg, 1 g		35DCzPPy SKU: M781 CAS Number: 1013405-25-8 Options: Quantity: 1 g
















High Purity Sublimed Materials

















	Bathocuproine (BCP) SKU: B231 CAS Number: 4733-39-5 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 1 g, 5 g		TmPyPB SKU: M2093 CAS Number: 921205-03-0 Options: • Quantity: 500 mg, 1 g
	TPD - N,N'-Bis(3-methylphenyl)-N,N'-diphenylbenzidine SKU: M0401 CAS Number: 65181-78-4 Options: • Quantity: 1 g, 5 g		C545T, Coumarin 545T SKU: M2094 CAS Number: 155306-71-1 Options: • Quantity: 500 mg, 1 g
	FIrPic SKU: M0711 CAS Number: 376367-93-0 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzIPN SKU: M2100 CAS Number: 1416881-52-1 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TAPC SKU: M0811 CAS Number: 58473-78-2 Options: • Quantity: 250 mg, 500 mg, 1 g		DBP SKU: M2101 CAS Number: 175606-05-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	3TPYMB SKU: M2088 CAS Number: 929203-02-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		PYD-2Cz (PYD2) SKU: M2102 CAS Number: 168127-49-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	DMAC-BP SKU: M2090 CAS Number: 1685287-55-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		2,4,6-tris(biphenyl-3-yl)-1,3,5-triazine (T2T) SKU: M2103 CAS Number: 1201800-83-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	ZADN SKU: M2091 CAS Number: 561064-11-7 Options: • Quantity: 500 mg, 1 g		CzSi SKU: M2104 CAS Number: 898546-82-2 Options: • Quantity: 100 mg, 250 mg, 500 mg
	MADN SKU: M2092 CAS Number: 804560-00-7 Options: • Quantity: 1 g, 5 g		DPEO SKU: M2105 CAS Number: 808142-23-6 Options: • Quantity: 250 mg, 500 mg, 1 g















	TBPe SKU: M2117 CAS Number: 80663-92-9 Options: • Quantity: 250 mg, 500 mg, 1 g		B3PyPB SKU: M2127 CAS Number: 1030380-38-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	ACRSA SKU: M2118 CAS Number: 1206626-95-0 Options: • Quantity: 1 g, 250 mg, 500 mg		TXO-TPA SKU: M2128 CAS Number: 1623010-63-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BCzPh SKU: M2119 CAS Number: 57102-62-2 Options: • Quantity: 250 mg, 500 mg, 1 g		TTPA SKU: M2129 CAS Number: 177799-16-5 Options: • Quantity: 250 mg, 500 mg, 1 g
	mCPSOB SKU: M2121 CAS Number: 1374770-41-8 Options: • Quantity: 250 mg, 500 mg, 1 g		DPPS SKU: M2130 CAS Number: 1152162-74-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMAC-DPS SKU: M2121 CAS Number: 1477512-32-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		UGH-2 SKU: M2131 CAS Number: 40491-34-7 / 18856-08-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PXZ-DPS SKU: M2122 CAS Number: 1477511-57-1 Options: • Quantity: 100 mg, 250 mg, 500 mg		PPT, PO15 SKU: M2173 CAS Number: 1019842-99-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	Cab-Ph-TRZ SKU: M2121 CAS Number: 440354-93-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		B2PymPm SKU: M2174 CAS Number: 1266181-51-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TBRb SKU: M2125 CAS Number: 682806-51-5 Options: • Quantity: 250 mg, 500 mg, 1 g		B3PymPm SKU: M2174 CAS Number: 925425-96-3 Options: • Purity: Sublimed >99.0% purity, Unsublimed >98.0% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g

	B4PymPm SKU: M2175 CAS Number: 1030380-51-8 Options: • Quantity: 250 mg, 500 mg, 1 g		mCBP SKU: M2186 CAS Number: 342638-54-4 Options: • Grade: Sublimed (>99%), Unsublimed (>98%) • Quantity: 500 mg, 1 g, 5 g
	B4PyPPm SKU: M2176 CAS Number: 1097652-83-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PPDN SKU: M2187 CAS Number: 215611-93-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BP4mPy SKU: M2178 CAS Number: 1009033-94-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		CDBP SKU: M2188 CAS Number: 120260-01-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	BTB SKU: M2179 CAS Number: 266349-83-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DPAVBi SKU: M2189 CAS Number: 119586-44-6 Options: • Quantity: 500 mg, 1 g
	2CzPN SKU: M2181 CAS Number: 1416881-50-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Tris-PCz SKU: M2191 CAS Number: 1141757-83-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	NPNPB SKU: M2182 CAS Number: 936355-01-0 Options: • Quantity: 500 mg, 1 g		Spiro-TTB SKU: M2192 CAS Number: 515834-67-0 Options: • Quantity: 500 mg, 1 g, 5 g
	NBPhen SKU: M2183 CAS Number: 1174006-43-9 Options: • Quantity: 250 mg, 500 mg, 1 g		DNTPD SKU: M2193 CAS Number: 199121-98-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	NPB-DPA SKU: M2185 CAS Number: 910058-11-6 Options: • Grade: Sublimed, Unsublimed • Quantity: 500 mg, 1 g, 5 g		BCBP SKU: M2194 CAS Number: 858131-70-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g














	Ir(mppy)3 SKU: M2195 CAS Number: 149005-33-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4TCzBN SKU: M2218 CAS Number: 2055722-93-3 Options: • Quantity: 100 mg, 250 mg, 500 mg
	BCPO SKU: M2198 CAS Number: 1233407-28-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		5CzBN SKU: M2219 CAS Number: 1469700-24-0 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TSP01 SKU: M2199 CAS Number: 1286708-86-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		5TCzBN SKU: M2220 CAS Number: 1469700-26-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TSBF SKU: M2213 CAS Number: 518997-91-6 Options: • Quantity: 500 mg, 1 g		BD-1 SKU: M2222 CAS Number: 1252914-27-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzTPN SKU: M2214 CAS Number: 1416881-53-2 Options: • Quantity: 250 mg, 500 mg, 1 g		Ir(dpm)PQ2 SKU: M2231 CAS Number: 713079-03-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzTPN-Me SKU: M2215 CAS Number: 1416881-54-3 Options: • Purity: Sublimed grade A1 (> 99.0%), Sublimed grade A2 (> 99.0%), Unsublimed grade (> 98%) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DMAC-BPP SKU: M2247 CAS Number: 1836192-40-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzTPN-Ph SKU: M2216 CAS Number: 1416881-55-4 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 250 mg, 500 mg, 1 g		4CzTPN-tBu SKU: M2251 CAS Number: 2153433-46-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	4CzBN SKU: M2217 CAS Number: 1996609-93-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Phen-TRZ SKU: M2264 CAS Number: 1357066-21-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

















	TXO-PhCz SKU: M2265 CAS Number: 1623010-64-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		3N-T2T SKU: M2275 CAS Number: 939430-26-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PO-T2T SKU: M2266 CAS Number: 1646906-26-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g, 25 g		DDCzTRZ SKU: M2286 CAS Number: 1685282-47-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DCzDBT SKU: M2267 CAS Number: 913738-04-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		CbBPCb SKU: M2289 CAS Number: 1469997-91-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PXZ-TRZ SKU: M2268 CAS Number: 1411910-25-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		FATPA SKU: M2290 CAS Number: 1131007-94-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	UGH3 SKU: M2270 CAS Number: 18920-16-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SimCP2 SKU: M2291 CAS Number: 944465-42-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	CzAcSF SKU: M2271 CAS Number: 1792173-34-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SPPO1 SKU: M2292 CAS Number: 1125547-88-7 Options: • Quantity: 250 mg, 500 mg
	mCPCN SKU: M2272 CAS Number: 1392506-99-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BSB SKU: M2293 CAS Number: 18826-13-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	3P-T2T SKU: M2273 CAS Number: 352196-01-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SimCP SKU: M2294 CAS Number: 850221-63-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

















	PPF SKU: M2295 CAS Number: 911397-27-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TCPZ SKU: M2314 CAS Number: 890148-68-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	Eu(dbm)3(Phen) SKU: M2296 CAS Number: 17904-83-5 Options: • Quantity: 250 mg, 500 mg, 1 g		TBADN SKU: M2330 CAS Number: 274905-73-6 Options: • Quantity: 250 mg, 500 mg, 1 g
	FK306 SKU: M2308 CAS Number: 1421058-47-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzFCN SKU: M2331 CAS Number: 1819362-10-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mSiTrz SKU: M2309 CAS Number: 2252416-90-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TCzTRZ SKU: M2332 CAS Number: 1808158-40-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(dpm)(piq)2 SKU: M2310 CAS Number: 1202867-58-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BCZTPA SKU: M2334 CAS Number: 1032174-52-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	TSBPA SKU: M2311 CAS Number: 205327-13-5 Options: • Quantity: 250 mg, 500 mg, 1 g		mCPPO1 SKU: M2335 CAS Number: 1296229-26-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Hex-Ir(phq)3 SKU: M2312 CAS Number: 1268460-37-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TmCzTRZ SKU: M2336 CAS Number: 1808158-41-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPTPCz SKU: M2313 CAS Number: 1313391-57-9 Options: • Quantity: 250 mg, 500 mg, 1 g		POSTF SKU: M2337 CAS Number: 1647050-25-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g















	DCzTRZ SKU: M2338 CAS Number: 1106730-48-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DPPyA SKU: M2356 CAS Number: 1257879-34-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	BPy-OXD SKU: M2339 CAS Number: 866117-19-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PPO27 SKU: M2357 CAS Number: 1299463-56-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DTAF SKU: M2342 CAS Number: 159526-57-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SiCzCz SKU: M2376 CAS Number: 2665768-26-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMAC-TRZ SKU: M2343 CAS Number: 1628752-98-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SiTrzCz2 SKU: M2377 CAS Number: 2422045-58-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMPPP SKU: M2344 CAS Number: 1036404-84-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Cz-DPS SKU: M2385 CAS Number: 733038-89-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mCBP-CN SKU: M2345 CAS Number: 1327163-09-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Cz-PS SKU: M2386 CAS Number: 1396165-20-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mCPBC SKU: M2354 CAS Number: 1799525-58-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		AQ(PhDPA)₂ SKU: M2387 CAS Number: 1640978-33-1 Options: • Purity grade: Sublimed 99%, Unsublimed 98% • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	46DzPPm (CzPhPy) SKU: M2354 CAS Number: 1262678-77-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Px-VPN SKU: M2388 CAS Number: 1784766-39-7 Options: • Quantity: 100 mg, 250 mg, 500 mg

















* Some options have been omitted, please see our website to browse our full range

	Ac-CNP SKU: M2389 CAS Number: 1883400-34-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		v-DABNA SKU: M2423 CAS Number: 2251782-13-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DACT-II SKU: M2390 CAS Number: 1613074-59-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DABNA-1 SKU: M2429 CAS Number: 1689552-89-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	α-NPD SKU: M2411 CAS Number: 495416-60-9 Options: • Quantity: 250 mg, 500 mg, 1 g		CzBN SKU: M2430 CAS Number: 2170487-27-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	ANT-BIZ SKU: M2415 CAS Number: 1254961-38-0 Options: • Quantity: 250 mg, 500 mg, 1 g		3tPAB SKU: M2431 CAS Number: 2056877-25-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TaTm SKU: M2416 CAS Number: 952431-34-4 Options: • Quantity: 250 mg, 500 mg, 1 g		t-Bu-v-DABNA SKU: M2432 CAS Number: 2914070-28-1 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g
	CzDBA SKU: M2417 CAS Number: 2171334-43-1 Options: • Quantity: 100 mg, 250 mg		PCzAc SKU: M2434 CAS Number: 1705584-08-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	TmPPPyTz SKU: M2419 CAS Number: 939430-31-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SpiroAC-TRZ SKU: M2435 CAS Number: 1980037-96-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	t-DABNA SKU: M2420 CAS Number: 1802003-07-1 Options: • Grade: Sublimed > 99.5%, Unsublimed >99% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BPBPA SKU: M2436 CAS Number: 164724-35-0 Options: • Quantity: 500 mg, 1 g, 5 g

	PhCzBCz SKU: M2479 CAS Number: 2692683-03-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DIC-TRZ SKU: M2507 CAS Number: 1024598-01-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DtBuCzB SKU: M2493 CAS Number: 2170487-30-4 Options: • Purity: Sublimed > 99%, Unsublimed > 97% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4DPAIPN SKU: M2508 CAS Number: 1846598-27-3 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g
	mMDBA-DI SKU: M2496 CAS Number: 2770257-85-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TDBA-SAF SKU: M2509 CAS Number: 2643306-28-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	3DPAFIPN SKU: M2501 CAS Number: 2260543-73-3 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g		TB-P3Cz SKU: M2511 CAS Number: 2411720-50-6 Options: • Quantity: 250 mg, 500 mg, 1 g
	3CzBN SKU: M2502 CAS Number: 1799923-03-7 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g		2PTZBN SKU: M2512 CAS Number: 2250280-17-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	3DPA2FBN SKU: M2503 CAS Number: 1403850-00-9 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g		3DPA-DiKta SKU: M2514 Options: • Quantity: 250 mg, 500 mg, 1 g
	DMIC-TRZ SKU: M2504 CAS Number: 1257248-13-7 Options: • Quantity: 250 mg, 500 mg, 1 g		TB-tPCz SKU: M2515 Options: • Quantity: 250 mg, 500 mg, 1 g
	CzPA SKU: M2506 CAS Number: 910647-23-3 Options: • Quantity: 250 mg, 500 mg, 1 g		Mes3DiKta SKU: M2516 Options: • Quantity: 250 mg, 500 mg, 1 g


	3TPA-DiKtA SKU: M2517 Options: Quantity: 250 mg, 500 mg, 1 g		tDPA-DtCzB SKU: M2540 CAS Number: 2760492-67-7 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 250 mg, 500 mg, 1 g
	T3PyTRZ SKU: M2519 CAS Number: 2703850-25-1 Options: - Grade: Sublimed $\geq 99\%$, Unsublimed $\geq 98\%$ - Quantity: 250 mg, 500 mg, 1 g		p-bPPhenB SKU: M2541 CAS Number: 2135536-38-6 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DBF-TRZ SKU: M2520 CAS Number: 2040300-54-5 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g		SiDBFCz SKU: M2543 CAS Number: 1234996-95-2 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	1,5-DTNA SKU: M2522 CAS Number: 2762325-17-5 Options: - Purity: Sublimed $\geq 99.0\%$, Unsublimed $\geq 98.0\%$ - Quantity: 250 mg, 500 mg, 1 g		DCz-BTP SKU: M2545 CAS Number: 1235883-57-4 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DOBNA-OAr SKU: M2536 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 250 mg, 500 mg, 1 g		Spiro-OMeTAD (Spiro-MeOTAD) SKU: M291 CAS Number: 207739-72-8 Options: - Grade: Sublimed ($>99.5\%$ purity), Unsublimed ($>99.0\%$ purity) - Quantity: 500 mg, 1 g, 5 g, 10 g
	4DBFHPB SKU: M2537 CAS Number: 2280942-02-9 Options: - Grade: Sublimed $\geq 99\%$, Unsublimed $\geq 98\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g		CuPc, Copper(II) phthalocyanine SKU: M341 CAS Number: 147-14-8 Options: Quantity: 1 g, 5 g
	BN-CP1 SKU: M2538 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 250 mg, 500 mg, 1 g		F4TCNQ SKU: M351 CAS Number: 29261-33-4 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g
	BPBiPA SKU: M2539 CAS Number: 1256723-60-0 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 500 mg, 1 g		NPB (NPD) SKU: M361 CAS Number: 123847-85-8 Options: - Grade: Sublimed Grade ($>99.5\%$ purity), Unsublimed Grade ($>98.0\%$ purity) - Quantity: 1 g, 5 g

	mCP SKU: M371 CAS Number: 550378-78-4 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) - Quantity: 1 g, 5 g		TCTA SKU: M471 CAS Number: 139092-78-7 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) - Quantity: 500 mg, 1 g
	Alq3 SKU: M381 CAS Number: 2085-33-8 Options: Quantity: 5 g, 10 g		Ir(ppy)3 SKU: M481 CAS Number: 94928-86-6 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed (>98.0% purity) - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	CBP, 4,4'-Bis(N-carbazolyl)-1,1'-biphenyl SKU: M391 CAS Number: 58328-31-7 Options: - Grade: Sublimed Grade (>99.5% purity) - Quantity: 1 g, 5 g		Ir(btpp)3 SKU: M491 CAS Number: 405289-74-9 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Zn(BTZ)2 SKU: M421 CAS Number: 58280-31-2 Options: Quantity: 1 g		Pentacene SKU: M431 CAS Number: 135-48-8 Options: Quantity: 1 g, 5 g
	BAiq SKU: M421 CAS Number: 146162-54-1 Options: Quantity: 1 g		TAZ, 3-(Biphenyl-4-yl)-5-(4-tert-butylphenyl)-4-phenyl-4H-1,2,4-triazole SKU: M601 CAS Number: 150405-69-9 Options: Quantity: 250 mg, 500 mg, 1 g
	DPVBi SKU: M431 CAS Number: 142289-08-5 Options: Quantity: 1 g		2-TNATA SKU: M611 CAS Number: 185690-41-9 Options: - Grade: Sublimed (>99.0% purity), Unsublimed (>98.0% purity) - Quantity: 500 mg, 1 g, 5 g
	OXD-7 SKU: M461 CAS Number: 138372-67-5 Options: Quantity: 250 mg, 500 mg, 1 g		m-MTDATA SKU: M621 CAS Number: 124729-98-2 Options: Quantity: 250 mg, 500 mg, 1 g
	ADN - 9,10-Bis(2-naphthyl)anthracene SKU: M461 CAS Number: 122648-99-1 Options: - Grade: Sublimed (>99.0% purity), Unsublimed (>98.0% purity) - Quantity: 1 g, 5 g		Cobalt(II) Phthalocyanine, CoPc SKU: M631 CAS Number: 3317-67-7 Options: Quantity: 500 mg, 1 g


	Ir(piq)3 SKU: M641 CAS Number: 435293-93-9 Options: • Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		8-Quinololato lithium (Liq) SKU: M731 CAS Number: 25387-93-3 Options: • Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) • Quantity: 1 g, 5 g
	TPBi SKU: M651 CAS Number: 192198-85-9 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 250 mg, 500 mg, 1 g		Bebp2 SKU: M741 CAS Number: 148896-39-3 Options: • Quantity: 1 g
	Ir(ppy)2(acac) SKU: M661 CAS Number: 337526-85-9 Options: • Grade: Sublimed (>99.5% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		tCP, 1,3,5-Tris(carbazol-9-yl)benzene SKU: M751 CAS Number: 148044-07-9 Options: • Grade: Sublimed (>99.0% purity), Unsublimed (>98.0% purity) • Quantity: 1 g, 5 g
	Ir(piq)2(acac) SKU: M671 CAS Number: 435294-03-4 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>99.0% purity) • Quantity: 250 mg, 500 mg, 1 g		Ir(ppy)2(bmp) SKU: M761 CAS Number: 1215692-34-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TCNQ SKU: M681 CAS Number: 1518-16-7 Options: • Quantity: 1 g, 5 g		HATCN SKU: M791 CAS Number: 105598-27-4 Options: • Grade: Sublimed (>99.0%), Unsublimed (>98.0%) • Quantity: 500 mg, 1 g, 5 g, 10 g
	Bepp2 SKU: M691 CAS Number: 220694-90-6 Options: • Quantity: 1 g		Ir(dmpq)2(acac) SKU: M801 CAS Number: 1056874-46-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Titanyl phthalocyanine, TiOPc SKU: M701 CAS Number: 26201-32-1 Options: • Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) • Quantity: 500 mg, 1 g, 5 g		BPhen (Bathophenanthroline) SKU: M961 CAS Number: 1662-01-7 Options: • Quantity: 1 g, 5 g
	Ir(ppz)3 SKU: M721 CAS Number: 562824-31-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DMQA SKU: M971 CAS Number: 19205-19-7 Options: • Quantity: 1 g, 5 g

Organic Conductors

Organic conductors are an environmentally friendly alternative to traditional inorganic conductors. They are composed of carbon-based molecules that facilitate the flow of electricity through their molecular structure due to the delocalization of electrons.




Popular Products



TIPS-Pentacene

High-purity, air stable, p-type organic semiconductor


- Excellent hole mobility
- Solution processable
- OFET applications



NDI-N

High-purity, n-type semiconductor with high electron affinity

- Cathode interlayer molecule
- π - π stacking self-organization
- polymer building block



PDIN

High-purity with favourable out-of-plane π - π stacking

- Improved surface morphology
- Cathode interlayer molecules
- Precipitates DNA

To help you choose the right material, full technical specifications are available on our website.

Ossila

M2283A1-100mg

PTQ10

Poly[[6,7-difluoro[(2-hexyldodecyloxy)-5,8-quinoline
diyl]-2,5-thiophenediyl]]











CAS# 2270233-86-6

Solpro Business Park

Solpro Business Park












undecyl
lene(2',3':
4',5')
arylidene)
indene-2,1-

Cathode Interlayer Materials (CIMs)

	PFN-FP SKU: M2078 CAS Number: 673474-79-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg		NDI-N SKU: M2237 CAS Number: 3436-54-2 Options: • Quantity: 1 g, 5 g, 10 g
	PFN-I SKU: M2079 CAS Number: 1687752-60-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg		Phen-NaDPO SKU: M2510 CAS Number: 1480371-38-7 Options: • Purity: Unsublimed > 97%, Unsublimed > 99% • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PFN-FP-I SKU: M2080 CAS Number: 1026670-72-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg		
	PFN, PFN-DOF SKU: M222 CAS Number: 673474-74-3 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		
	PFN-bromide, PFN-Br SKU: M2230 CAS Number: 889672-99-5 Options*: • Quantity: 100 mg, 250 mg, 500 mg		
	PDIN SKU: M2232 CAS Number: 117901-97-0 Options: • Quantity: 1 g, 5 g, 10 g		
	PDINN SKU: M2234 CAS Number: 1020180-01-1 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		
	PDINO SKU: M2236 CAS Number: 1558023-86-1 Options: • Quantity: 250 mg, 500 mg, 1 g		

* Some options have been omitted, please see our website to browse our full range

Green Energy Materials

	Y6, BTP-4F SKU: M2200 CAS Number: 2304444-49-1 Options*: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g		L8-BO, L8-BO-2F SKU: M2304 CAS Number: 2668341-40-8 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg
	BTP-4F-12, Y6-BO (Y12) SKU: M2245 CAS Number: 2389125-23-7 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g		L8-BO-F SKU: M2306 Options: • Quantity: 50 mg, 100 mg, 250 mg
	BTP-eC9 SKU: M2249 CAS Number: 2598965-39-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BTP-4Cl-12 SKU: M2327 CAS Number: 2447642-41-1 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg
	TPT10-C8C12 SKU: M2253 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Lithium Iron Phosphate (LiFePO4) Powder SKU: M2402 CAS Number: 15365-14-7 Options*: • Quantity: 100 g, 250 g
	PTO2 SKU: M2260 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Lithium Nickel Manganese Cobalt Oxide Powder (NCM523) SKU: M2409 CAS Number: 346417-97-8 / 182442-95-1 Options*: • Quantity: 50 g, 100 g, 250 g
	DTY6 SKU: M2274 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg		Hydroxyapatite Powder SKU: M2437 CAS Number: 12167-74-7 Options: • Average particle size: 3 µm, 10 - 20 µm • Quantity: 10 g, 25 g, 50 g
	PCITO2 SKU: M2278 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Iron Oxide (Fe₃O₄) Nanopowder SKU: M2439 CAS Number: 1317-61-9 Options: • Average particle size: 30 - 50 nm • Quantity: 100 g, 250 g, 500 g
	eC9-2Cl SKU: M2301 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Iron Oxide (Fe₂O₃) Nanopowder SKU: M2440 CAS Number: 1309-37-1 Options: • Crystalline phase: α-Fe ₂ O ₃ , γ-Fe ₂ O ₃ • Quantity: 10 g, 25 g, 50 g

* Some options have been omitted, please see our website to browse our full range



Silicon Powder
SKU: M2442 | CAS Number: 7440-21-3
Options:
• Average particle size: 90 nm, 10 μ m
• Quantity: 5 g, 10 g, 25 g, 50 g


































Silicon Nanopowder
SKU: M2442 | CAS Number: 7440-21-3
Options:
• Average particle size: 90 nm
• Quantity: 5 g, 10 g


































Lithium Manganese Iron Phosphate (LMFP) Powder
SKU: M2491 | CAS Number: 1354700-46-1
Options:
• Quantity: 50 g, 100 g, 250 g, 500 g















Semiconducting Molecules
















	Bathocuproine (BCP) SKU: B231 CAS Number: 4733-39-5 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 1 g, 5 g		TAPC SKU: M0811 CAS Number: 58473-78-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	4,4,5,5-Tetramethyl-2-[4-(triphenylvinyl)-phenyl]-1,3,2-dioxaborolane SKU: B4061 CAS Number: 1260865-91-5 Options: • Quantity: 500 mg, 1 g		DRCN5T SKU: M2084 CAS Number: 1674394-69-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DTT, Dithieno[3,2-b:2',3'-d]thiophene SKU: B931 CAS Number: 3593-75-7 Options: • Quantity: 1 g, 5 g		DTS(FBTTh₂)₂ SKU: M2085 CAS Number: 1402460-84-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TIPS-Pentacene SKU: M0151 CAS Number: 373596-08-8 Options: • Quantity: 100 mg, 500 mg, 1 g, 5 g		3TPYMB SKU: M2088 CAS Number: 929203-02-1 Options: • Quantity: 100 mg, 250 mg, 500 mg
	TPD, N,N'-Bis(3-methylphenyl)-N,N'-diphenylbenzidine SKU: M0401 CAS Number: 65181-78-4 Options: • Quantity: 1 g, 5 g		DMAC-BP SKU: M2090 CAS Number: 1685287-55-1 Options: • Quantity: 100 mg, 250 mg, 500 mg
	4-tert-Butylpyridine SKU: M0591 CAS Number: 3978-81-2 Options: • Quantity: 25 g, 50 g, 100 g		ZADN SKU: M2091 CAS Number: 561064-11-7 Options: • Quantity: 500 mg, 1 g
	FirPic SKU: M0711 CAS Number: 376367-93-0 Options: • Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		MADN SKU: M2092 CAS Number: 804560-00-7 Options: • Quantity: 1 g, 5 g
	26Dc2PPy SKU: M0771 CAS Number: 1013405-24-7 Options: • Grade: Sublimed Grade (>99% Purity), Unsublimed Grade (>99% Purity) • Quantity: 250 mg, 500 mg, 1 g		TmPyPB SKU: M2093 CAS Number: 921205-03-0 Options: • Quantity: 500 mg, 1 g
















	C545T, Coumarin 545T SKU: M2094 CAS Number: 155306-71-1 Options: • Quantity: 500 mg, 1 g		ACRSA SKU: M2118 CAS Number: 1206626-95-0 Options: • Quantity: 1 g, 250 mg, 500 mg
	4CzIPN SKU: M2100 CAS Number: 1416881-52-1 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BCzPh SKU: M2119 CAS Number: 57102-62-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	DBP SKU: M2101 CAS Number: 175606-05-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCPSOB SKU: M2120 CAS Number: 1374770-41-8 Options: • Quantity: 250 mg, 500 mg, 1 g
	PYD-2Cz (PYD2) SKU: M2102 CAS Number: 168127-49-9 Options: • Quantity: 250 mg, 500 mg, 1 g		DMAC-DPS SKU: M2121 CAS Number: 1477512-32-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	2,4,6-tris(biphenyl-3-yl)-1,3,5-triazine (T2T) SKU: M2103 CAS Number: 1201800-83-0 Options: • Quantity: 250 mg, 500 mg, 1 g		PXZ-DPS SKU: M2122 CAS Number: 1477511-57-1 Options: • Quantity: 100 mg, 250 mg, 500 mg
	CzSi SKU: M2104 CAS Number: 898546-82-2 Options: • Quantity: 100 mg, 250 mg, 500 mg		Cab-Ph-TRZ SKU: M2123 CAS Number: 440354-93-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPEO SKU: M2105 CAS Number: 808142-23-6 Options: • Quantity: 250 mg, 500 mg, 1 g		TBRb SKU: M2125 CAS Number: 682806-51-5 Options: • Quantity: 250 mg, 500 mg, 1 g
	TBPe SKU: M2117 CAS Number: 80663-92-9 Options: • Quantity: 250 mg, 500 mg, 1 g		Flr6 SKU: M2126 CAS Number: 664374-03-2 Options: • Quantity: 100 mg, 250 mg

	NBPhen SKU: M2183 CAS Number: 1174006-43-9 Options: • Quantity: 250 mg, 500 mg, 1 g		Spiro-TTB SKU: M2192 CAS Number: 515834-67-0 Options: • Quantity: 500 mg, 1 g, 5 g
	Ir(MDQ)2(acac) SKU: M2184 CAS Number: 536755-34-7 Options: • Quantity: 100 mg, 250 mg		DNTPD SKU: M2193 CAS Number: 199121-98-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	NPB-DPA SKU: M2185 CAS Number: 910058-11-6 Options: • Grade: Sublimed, Unsublimed • Quantity: 500 mg, 1 g, 5 g		BCBP SKU: M2194 CAS Number: 858131-70-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mCBP SKU: M2186 CAS Number: 342638-54-4 Options: • Grade: Sublimed (>99%), Unsublimed (>98%) • Quantity: 500 mg, 1 g, 5 g		Ir(mppy)3 SKU: M2195 CAS Number: 149005-33-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PPDN SKU: M2187 CAS Number: 215611-93-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BCPO SKU: M2196 CAS Number: 1233407-28-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	CDBP SKU: M2188 CAS Number: 120260-01-7 Options: • Quantity: 250 mg, 500 mg, 1 g		TSPO1 SKU: M2199 CAS Number: 1286708-86-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPAVBi SKU: M2189 CAS Number: 119586-44-6 Options: • Quantity: 500 mg, 1 g		TSBF SKU: M2213 CAS Number: 518997-91-6 Options: • Quantity: 500 mg, 1 g
	Tris-PCz SKU: M2191 CAS Number: 1141757-83-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzTPN SKU: M2214 CAS Number: 1416881-53-2 Options: • Quantity: 250 mg, 500 mg, 1 g

















	B3PyPB SKU: M2127 CAS Number: 1030380-38-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		B4PymPm SKU: M2175 CAS Number: 1030380-51-8 Options: • Quantity: 250 mg, 500 mg, 1 g
	TXO-TPA SKU: M2128 CAS Number: 1623010-63-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		B4PyPPm SKU: M2176 CAS Number: 1097652-83-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TTPA SKU: M2129 CAS Number: 177799-16-5 Options: • Quantity: 250 mg, 500 mg, 1 g		DCJTb SKU: M2177 CAS Number: 200052-70-6 Options: • Quantity: 100 mg, 250 mg, 500 mg
	DPPS SKU: M2130 CAS Number: 1152162-74-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BP4mPy SKU: M2178 CAS Number: 1009033-94-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	UGH-2 SKU: M2131 CAS Number: 40491-34-7 / 18856-08-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BTB SKU: M2179 CAS Number: 266349-83-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PPT, PO15 SKU: M2132 CAS Number: 1019842-99-9 Options: • Quantity: 250 mg, 500 mg, 1 g		Picene SKU: M2180 CAS Number: 213-46-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	B2PymPm SKU: M2131 CAS Number: 1266181-51-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		2CzPN SKU: M2182 CAS Number: 1416881-50-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	B3PymPm SKU: M2174 CAS Number: 925425-96-3 Options: • Purity: Sublimed >99.0% purity, Unsublimed >98.0% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g		NPNPB SKU: M2182 CAS Number: 936355-01-0 Options: • Quantity: 500 mg, 1 g















	DCzDBT SKU: M2267 CAS Number: 913738-04-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DDcZTRZ SKU: M2286 CAS Number: 1685282-47-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PXZ-TRZ SKU: M2268 CAS Number: 1411910-25-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		CbBPCb SKU: M2289 CAS Number: 1469997-91-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	UGH3 SKU: M2270 CAS Number: 18920-16-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		FATPA SKU: M2290 CAS Number: 1131007-94-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	CzAcSF SKU: M2271 CAS Number: 1792173-34-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SimCP2 SKU: M2291 CAS Number: 944465-42-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mCPCN SKU: M2272 CAS Number: 1392506-99-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SPPO1 SKU: M2292 CAS Number: 1125547-88-7 Options: • Quantity: 250 mg, 500 mg
	3P-T2T SKU: M2273 CAS Number: 352196-01-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BSB SKU: M2293 CAS Number: 18826-13-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	3N-T2T SKU: M2275 CAS Number: 939430-26-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SimCP SKU: M2294 CAS Number: 850221-63-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	ZR1 SKU: M2285 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PPF SKU: M2295 CAS Number: 911397-27-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g















	4CzTPN-Me SKU: M2215 CAS Number: 1416881-54-3 Options: • Purity: Sublimed grade A1 (> 99.0%), Sublimed grade A2 (> 99.0%), Unsublimed grade (> 98%) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Ir(dpm)PQ2 SKU: M2231 CAS Number: 713079-03-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzTPN-Ph SKU: M2216 CAS Number: 1416881-55-4 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 250 mg, 500 mg, 1 g		DMAC-BPP SKU: M2247 CAS Number: 1836192-40-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzBN SKU: M2217 CAS Number: 1996609-93-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzTPN-tBu SKU: M2251 CAS Number: 2153433-46-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	4TCzBN SKU: M2218 CAS Number: 2055722-93-3 Options: • Quantity: 100 mg, 250 mg, 500 mg		BTR SKU: M2255 CAS Number: 2041283-06-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	5CzBN SKU: M2219 CAS Number: 1469700-24-0 Options: • Grade: Sublimed (>99% purity), Unsublimed (>98% purity) • Quantity: 100 mg, 250 mg, 500 mg, 1 g		BTR-Cl SKU: M2257 CAS Number: 2410661-17-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	5TCzBN SKU: M2220 CAS Number: 1469700-26-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		Phen-TRZ SKU: M2264 CAS Number: 1357066-21-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzPN SKU: M2221 CAS Number: 1416881-51-0 Options: • Grade: Sublimed > 99% purity, Unsublimed > 98% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TXO-PhCz SKU: M2265 CAS Number: 1623010-64-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BD-1 SKU: M2222 CAS Number: 1252914-27-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PO-T2T SKU: M2266 CAS Number: 1646906-26-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g, 25 g
















	DCzTRZ SKU: M2338 CAS Number: 1106730-48-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		46DCzPPm (CzPhPy) SKU: M2354 CAS Number: 1262678-77-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BPy-OXD SKU: M2339 CAS Number: 866117-19-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DPPyA SKU: M2356 CAS Number: 1257879-34-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	DTAF SKU: M2342 CAS Number: 159526-57-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		PPO27 SKU: M2357 CAS Number: 1299463-56-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DMAC-TRZ SKU: M2343 CAS Number: 1628752-98-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzPN-Ph SKU: M2376 CAS Number: 1469707-47-8 Options: • Quantity: 250 mg, 500 mg, 1 g
	DMPPP SKU: M2344 CAS Number: 1036404-84-8 Options*: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SiCzCz SKU: M2376 CAS Number: 2665768-26-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mCBP-CN SKU: M2345 CAS Number: 1327163-09-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SiTrzCz2 SKU: M2377 CAS Number: 2422045-58-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	2-(2-Thienylmethylene)indane-1,3-dione SKU: M2352 CAS Number: 24760-36-9 Options: • Quantity: 1 g, 5 g, 10 g, 25 g		4CzPN-Me SKU: M2378 CAS Number: 1469700-28-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	mCPBC SKU: M2353 CAS Number: 1799525-58-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzPN-tBu SKU: M2379 CAS Number: 1469705-93-8 Options: • Grade: Sublimed > 99% purity, Unsublimed > 98% purity • Quantity: 100 mg, 250 mg, 500 mg, 1 g

* Some options have been omitted, please see our website to browse our full range
















	Eu(dbm)3(Phen) SKU: M2296 CAS Number: 17904-83-5 Options: • Quantity: 250 mg, 500 mg, 1 g		TBADN SKU: M2330 CAS Number: 274905-73-6 Options: • Quantity: 250 mg, 500 mg, 1 g
	FK306 SKU: M2308 CAS Number: 1421058-47-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4CzFCN SKU: M2331 CAS Number: 1819362-10-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	mSiTrz SKU: M2309 CAS Number: 2252416-90-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TCzTRZ SKU: M2332 CAS Number: 1808158-40-8 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(dpm)(piq)2 SKU: M2310 CAS Number: 1202867-58-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		SPPO13 SKU: M2333 CAS Number: 1234510-13-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	TSBPA SKU: M2311 CAS Number: 205327-13-5 Options: • Quantity: 250 mg, 500 mg, 1 g		BCzTPA SKU: M2334 CAS Number: 1032174-52-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	Hex-Ir(phq)3 SKU: M2312 CAS Number: 1268460-37-2 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		mCPPO1 SKU: M2335 CAS Number: 1296229-26-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	DPTPCz SKU: M2313 CAS Number: 1313391-57-9 Options: • Quantity: 250 mg, 500 mg, 1 g		TmCzTRZ SKU: M2336 CAS Number: 1808158-41-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TCPZ SKU: M2314 CAS Number: 890148-68-2 Options: • Quantity: 250 mg, 500 mg, 1 g		POSTF SKU: M2337 CAS Number: 1647050-25-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g

















	TmPPPyTz SKU: M2419 CAS Number: 939430-31-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		3tPAB SKU: M2431 CAS Number: 2056877-25-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	t-DABNA SKU: M2420 CAS Number: 1802003-07-1 Options: • Grade: Sublimed > 99.5%, Unsublimed >99% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		t-Bu-v-DABNA SKU: M2432 CAS Number: 2914070-28-1 Options: • Quantity: 50 mg, 100 mg, 250 mg, 500 mg, 1 g
	v-DABNA SKU: M2423 CAS Number: 2251782-13-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		2,3,6,7,10,11-Hexakis(hexyloxy)triphenylene - HAT6 SKU: M2433 CAS Number: 70351-86-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	2PACz SKU: M2425 CAS Number: 20999-38-6 Options: • Quantity: 500 mg, 1 g		PCzAc SKU: M2434 CAS Number: 1705584-08-2 Options: • Quantity: 250 mg, 500 mg, 1 g
	Br-2PACz SKU: M2426 CAS Number: 2762888-11-7 Options: • Quantity: 500 mg, 1 g		SpiroAC-TRZ SKU: M2435 CAS Number: 1980037-96-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	MeO-2PACz SKU: M2427 CAS Number: 2377770-18-6 Options: • Quantity: 500 mg, 1 g		BPBPA SKU: M2436 CAS Number: 164724-35-0 Options: • Quantity: 500 mg, 1 g, 5 g
	DABNA-1 SKU: M2429 CAS Number: 1689552-89-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DMAcPA SKU: M2445 CAS Number: 2971088-37-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	CzBN SKU: M2430 CAS Number: 2170487-27-9 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4PADCB SKU: M2446 CAS Number: 2882156-63-8 Options: • Quantity: 250 mg, 500 mg, 1 g
















	4CzIPN-Me SKU: M2380 CAS Number: 1469703-61-4 Options: • Quantity: 250 mg, 500 mg, 1 g		DACT-II SKU: M2390 CAS Number: 1613074-59-1 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	4CzIPN-tBu SKU: M2381 CAS Number: 1630263-99-8 Options: • Quantity: 250 mg, 500 mg, 1 g		1,4-Di-tert-butyl-2,5-dimethoxybenzene SKU: M2391 CAS Number: 7323-63-9 Options: • Quantity: 25 g, 100 g, 250 g, 500 g, 1 kg
	4CzIPN-Ph SKU: M2382 CAS Number: 1469705-37-0 Options: • Quantity: 250 mg, 500 mg, 1 g		TMTES-Pentacene SKU: M241 Options: • Quantity: 200 mg
	Cz-DPS SKU: M2385 CAS Number: 733038-89-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		α-NPD SKU: M2411 CAS Number: 495416-60-9 Options: • Quantity: 250 mg, 500 mg, 1 g
	Cz-PS SKU: M2386 CAS Number: 1396165-20-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		ANT-BIZ SKU: M2415 CAS Number: 1254961-38-0 Options: • Quantity: 250 mg, 500 mg, 1 g
	AQ(PhDPA)₂ SKU: M2387 CAS Number: 1640978-33-1 Options: • Purity grade: Sublimed 99%, Unsublimed 98% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TaTm SKU: M2416 CAS Number: 952431-34-4 Options: • Quantity: 250 mg, 500 mg, 1 g
	Px-VPN SKU: M2388 CAS Number: 1784766-39-7 Options: • Quantity: 100 mg, 250 mg, 500 mg		CzDBA SKU: M2417 CAS Number: 2171334-43-1 Options: • Quantity: 100 mg, 250 mg
	Ac-CNP SKU: M2389 CAS Number: 1883400-34-7 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TPA-DCPP SKU: M2418 CAS Number: 1803287-94-6 Options: • Quantity: 250 mg, 500 mg, 1 g









	Benzyl Viologen Dichloride SKU: M2469 CAS Number: 1102-19-8 Options: • Quantity: 5 g, 10 g, 25 g		3CzBN SKU: M2502 CAS Number: 1799923-03-7 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g
	Indium Oxide Powder SKU: M2488 CAS Number: 1312-43-2 Options: • Quantity: 5 g, 10 g, 25 g		3DPA2FBN SKU: M2503 CAS Number: 1403850-00-9 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g
	DtBuCzB SKU: M2493 CAS Number: 2170487-30-4 Options: • Purity: Sublimed > 99%, Unsublimed > 97% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DMIC-TRZ SKU: M2504 CAS Number: 1257248-13-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	mMDBA-DI SKU: M2496 CAS Number: 2770257-85-5 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		CzPA SKU: M2506 CAS Number: 910647-23-3 Options: • Quantity: 250 mg, 500 mg, 1 g
	3CzCIIPN SKU: M2507 CAS Number: 1469704-61-7 Options: • Quantity: 250 mg, 500 mg, 1 g, 5 g		DIC-TRZ SKU: M2507 CAS Number: 1024598-01-3 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	PPCzTRz SKU: M2499 CAS Number: 2329651-89-8 Options: • Grade: Sublimed ≥99.0%, Unsublimed ≥98.0% • Quantity: 100 mg, 250 mg, 500 mg, 1 g		4DPAIPN SKU: M2508 CAS Number: 1846598-27-3 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g
	5Cz-TRZ SKU: M2500 CAS Number: 2043597-42-6 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		TDBA-SAF SKU: M2509 CAS Number: 2643306-28-7 Options: • Quantity: 250 mg, 500 mg, 1 g
	3DPAFIPN SKU: M2501 CAS Number: 2260543-73-3 Options: • Grade: Sublimed, Unsublimed • Quantity: 250 mg, 500 mg, 1 g		Phen-NaDPO SKU: M2510 CAS Number: 1480371-38-7 Options: • Purity: Unsublimed > 97%, Unsublimed > 99% • Quantity: 100 mg, 250 mg, 500 mg, 1 g

www.ossila.com info@ossila.com +44 (0)114 2999 180

	SiDBFCz SKU: M2543 CAS Number: 1234996-95-2 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g		CBP, 4,4'-Bis(N-carbazolyl)-1,1'-biphenyl SKU: M391 CAS Number: 58328-31-7 Options: - Grade: Sublimed Grade ($>99.5\%$ purity) - Quantity: 1 g, 5 g
	DCz-BTP SKU: M2545 CAS Number: 1235883-57-4 Options: - Grade: Sublimed $\geq 99\%$ - Quantity: 100 mg, 250 mg, 500 mg, 1 g		Zn(BTZ)2 SKU: M411 CAS Number: 58280-31-2 Options: Quantity: 1 g
	Spiro-OMeTAD (Spiro-MeOTAD) SKU: M291 CAS Number: 207739-72-8 Options: - Grade: Sublimed ($>99.5\%$ purity), Unsublimed ($>99.0\%$ purity) - Quantity: 500 mg, 1 g, 5 g, 10 g		BAlq SKU: M421 CAS Number: 146162-54-1 Options: Quantity: 1 g
	CuPc, Copper(II) phthalocyanine SKU: M341 CAS Number: 147-14-8 Options: Quantity: 1 g, 5 g		DPVBi SKU: M431 CAS Number: 142289-08-5 Options: Quantity: 1 g
	F4TCNQ SKU: M351 CAS Number: 29261-33-4 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g, 5 g, 10 g		Rubrene SKU: M441 CAS Number: 517-51-1 Options: Quantity: 500 mg, 1 g
	NPB (NPD) SKU: M361 CAS Number: 123847-85-8 Options: - Grade: Sublimed Grade ($>99.5\%$ purity), Unsublimed Grade ($>98.0\%$ purity) - Quantity: 1 g, 5 g		OXD-7 SKU: M451 CAS Number: 138372-67-5 Options: Quantity: 250 mg, 500 mg, 1 g
	mCP SKU: M371 CAS Number: 550378-78-4 Options: - Grade: Sublimed Grade ($>99.5\%$ purity), Unsublimed Grade ($>98.0\%$ purity) - Quantity: 1 g, 5 g		ADN - 9,10-Bis(2-naphthyl)anthracene SKU: M461 CAS Number: 122648-99-1 Options: - Grade: Sublimed ($>99.0\%$ purity), Unsublimed ($>98.0\%$ purity) - Quantity: 1 g, 5 g
	Alq3 SKU: M381 CAS Number: 2085-33-8 Options: Quantity: 5 g, 10 g		TCTA SKU: M471 CAS Number: 139092-78-7 Options: - Grade: Sublimed Grade ($>99.5\%$ purity), Unsublimed Grade ($>98.0\%$ purity) - Quantity: 500 mg, 1 g

	TB-P3Cz SKU: M2511 CAS Number: 2411720-50-6 Options: • Quantity: 250 mg, 500 mg, 1 g		DBF-TRZ SKU: M2520 CAS Number: 2040300-54-5 Options: • Grade: Sublimed $\geq 99\%$ • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	2PTZBN SKU: M2512 CAS Number: 2250280-17-0 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		1,5-DTNA SKU: M2522 CAS Number: 2762325-17-5 Options: • Purity: Sublimed $\geq 99.0\%$, Unsublimed $\geq 98.0\%$ • Quantity: 250 mg, 500 mg, 1 g
	DiKTa SKU: M2513 CAS Number: 32081-22-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g		DOBNA-OAr SKU: M2536 Options: • Grade: Sublimed $\geq 99\%$ • Quantity: 250 mg, 500 mg, 1 g
	3DPA-DiKTa SKU: M2514 Options: • Quantity: 250 mg, 500 mg, 1 g		4DBFHPB SKU: M2537 CAS Number: 2280942-02-9 Options: • Grade: Sublimed $\geq 99\%$, Unsublimed $\geq 98\%$ • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	TB-tPCz SKU: M2515 Options: • Quantity: 250 mg, 500 mg, 1 g		BN-CP1 SKU: M2538 Options: • Grade: Sublimed $\geq 99\%$ • Quantity: 250 mg, 500 mg, 1 g
	Mes3DiKTa SKU: M2516 Options: • Quantity: 250 mg, 500 mg, 1 g		BPBiPA SKU: M2539 CAS Number: 1256723-68-0 Options: • Grade: Sublimed $\geq 99\%$ • Quantity: 500 mg, 1 g
	3TPA-DiKTa SKU: M2517 Options: • Quantity: 250 mg, 500 mg, 1 g		tDPA-DtCzB SKU: M2540 CAS Number: 2760492-67-7 Options: • Grade: Sublimed $\geq 99\%$ • Quantity: 250 mg, 500 mg, 1 g
	T3PyTRZ SKU: M2519 CAS Number: 2703850-25-1 Options: • Grade: Sublimed $\geq 99\%$, Unsublimed $\geq 98\%$ • Quantity: 250 mg, 500 mg, 1 g		p-bPPhenB SKU: M2541 CAS Number: 2135536-38-6 Options: • Grade: Sublimed $\geq 99\%$ • Quantity: 100 mg, 250 mg, 500 mg, 1 g

	Ir(ppy)3 SKU: M481 CAS Number: 94928-86-6 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed (>98.0% purity) - Quantity: 100 mg, 250 mg, 500 mg, 1 g		TPBi SKU: M651 CAS Number: 192198-85-9 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) - Quantity: 250 mg, 500 mg, 1 g
	Ir(btpp)3 SKU: M491 CAS Number: 405289-74-9 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g		Ir(ppy)2(acac) SKU: M661 CAS Number: 337526-85-9 Options: - Grade: Sublimed (>99.5% purity) - Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Pentacene SKU: M581 CAS Number: 135-48-8 Options: Quantity: 1 g, 5 g		Ir(piq)2(acac) SKU: M671 CAS Number: 435294-03-4 Options: - Grade: Sublimed Grade (>99.5% purity), Unsublimed Grade (>98.0% purity) - Quantity: 250 mg, 500 mg, 1 g
	TAZ, 3-(Biphenyl-4-yl)-5-(4-tert-butylphenyl)-4-phenyl-4H-1,2,4-triazole SKU: M601 CAS Number: 150405-69-9 Options: Quantity: 250 mg, 500 mg, 1 g		TCNQ SKU: M681 CAS Number: 1518-16-7 Options: Quantity: 1 g, 5 g
	2-TNATA SKU: M611 CAS Number: 185690-41-9 Options: - Grade: Sublimed (>99.0% purity), Unsublimed (>98.0% purity) - Quantity: 500 mg, 1 g, 5 g		Bepp2 SKU: M691 CAS Number: 220694-90-6 Options: Quantity: 1 g
	m-MTDATA SKU: M621 CAS Number: 124729-98-2 Options: Quantity: 250 mg, 500 mg, 1 g		Titanyl phthalocyanine, TiOPc SKU: M701 CAS Number: 26201-32-1 Options: - Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) - Quantity: 500 mg, 1 g, 5 g
	Cobalt(II) Phthalocyanine, CoPc SKU: M631 CAS Number: 3317-67-7 Options: Quantity: 500 mg, 1 g		Ir(ppz)3 SKU: M721 CAS Number: 562824-31-1 Options: Quantity: 100 mg, 250 mg, 500 mg, 1 g
	Ir(piq)3 SKU: M641 CAS Number: 435293-93-9 Options: - Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) - Quantity: 100 mg, 250 mg, 500 mg, 1 g		8-Quinololinato lithium (Liq) SKU: M731 CAS Number: 25387-93-3 Options: - Grade: Sublimed Grade (>99% purity), Unsublimed Grade (>98% purity) - Quantity: 1 g, 5 g

	Bebq2 SKU: M741 CAS Number: 148896-39-3 Options: • Quantity: 1 g
	tCP, 1,3,5-Tris(carbazol-9-yl)benzene SKU: M751 CAS Number: 148044-07-9 Options: • Grade: Sublimed (>99.0% purity), Unsublimed (>98.0% purity) • Quantity: 1 g, 5 g
	Ir(ppy)2(bmp) SKU: M761 CAS Number: 1215692-34-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	35DCzPPy SKU: M781 CAS Number: 1013405-25-8 Options: • Quantity: 1 g
	HATCN SKU: M791 CAS Number: 105598-27-4 Options: • Grade: Sublimed (>99.0%), Unsublimed (>98.0%) • Quantity: 500 mg, 1 g, 5 g, 10 g
	Ir(dmpq)2(acac) SKU: M801 CAS Number: 1056874-46-4 Options: • Quantity: 100 mg, 250 mg, 500 mg, 1 g
	BPhen (Bathophenanthroline) SKU: M961 CAS Number: 1662-01-7 Options: • Quantity: 1 g, 5 g
	DMQA SKU: M971 CAS Number: 19205-19-7 Options: • Quantity: 1 g, 5 g

PEDOT & PEDOT:PSS

PEDOT:PSS is a benchmark material in thin-film electronics, renowned for its versatility and wide range of applications. Its properties can be tuned by adjusting the ratio of conductive PEDOT to insulating PSS, as well as controlling the microstructure of the resulting film. This makes it possible to achieve precise control over conductivity and work function, enabling use in everything from transparent electrodes to surface modification layers. At Ossila, we offer a broad range of PEDOT:PSS formulations for different research needs, giving you the flexibility to select the material best suited to your project.



Popular Products



PH 1000 PEDOT:PSS

Popular high conductivity PEDOT:PSS formulation

- Optical Transparency
- Versatile Use
- High Quality



HTL Solar 3 PEDOT Complex

Toulene-based formulation for OPVs & perovskite solar cells

- Low conductivity
- Suitable for water sensitive environments



S V4 STAB PEDOT:PSS

Ideal for transparent electrodes and pressure sensors

- Blue stabilized dispersion/paste
- Highly transparent
- Tailored conductivity

To help you choose the right material, full technical specifications are available on our website.



Cossita

M1214-25ml

PEDOT:Complex
(HTL Solar 4)

poly(3,4-ethylenedioxythiophene);
complex in non-aqueous solution
(amine based)
Keep refrigerated - store above
0°C
CAS# N/A

Solpro Business Park

Highly flammable liquid and
vapour. May cause respiratory irritation
of the organs (lungs) through
inhalation. Keep away from heat,
smoking, fire or open flame.
Do not mix with water. If inhaled,
call a POISON CENTER or doctor
immediately and keep comfortable for
transportation. In case of fire use dry powder or
water to cool.

FREZETM



>0° C / 32° F =



<0° C / 32° F =



100ml

PEDOT:PSS (PH 1000)

poly(3,4-ethylenedioxythiophene);poly(styrene)
sulfonate in aqueous solution

Keep refrigerated - store above 0°C

155590-43-8

Solpro Business Park Windsor Street

PSS (S V4 STAB)

poly(3,4-ethylenedioxythiophene)sulfonate
(amine based)

Contains isocyanates. May
contain epoxy constituents.

PEDOT Polymer Blends

	AI 4083 PEDOT:PSS SKU: M121 CAS Number: 155090-83-8 Options: • Quantity: 100 ml, 1000 ml		P JET X N PEDOT:PSS SKU: M1217 CAS Number: 155090-83-8 Options: • Volume: 25 ml, 50 ml, 100 ml
	HIL 8 PEDOT Complex SKU: M1210 Options: • Title: Default Title		PH 1000 PEDOT:PSS SKU: M1217 CAS Number: 155090-83-8 Options: • Title: Default Title
	F 020 PEDOT:PSS SKU: M1211 CAS Number: 155090-83-8 Options: • Quantity: 100 ml, 500 ml, 1000 ml		HTL Solar PEDOT:PSS SKU: M124 CAS Number: 155090-83-8 Options: • Title: Default Title
	F ET PEDOT:PSS SKU: M1212 CAS Number: 155090-83-8 Options: • Quantity: 100 ml		HTL Solar 3 PEDOT Complex SKU: M125 Options: • Title: Default Title
	P T4 PEDOT:PSS SKU: M1213 CAS Number: 155090-83-8 Options: • Quantity: 100 ml, 500 ml, 1000 ml		F HC Solar PEDOT:PSS SKU: M126 CAS Number: 155090-83-8 Options: • Title: Default Title
	HTL Solar 4 PEDOT Complex SKU: M1214 Options: • Quantity: 25 ml, 50 ml		S V4 STAB PEDOT:PSS SKU: M128 CAS Number: 155090-83-8 Options: • Title: Default Title
	F 100T PEDOT:PSS SKU: M1215 CAS Number: 155090-83-8 Options: • Quantity: 100 mL, 1 L		P JET (OLED) PEDOT:PSS SKU: M129 CAS Number: 155090-83-8 Options: • Title: Default Title
	HY E PEDOT:PSS - AgNW/PEDOT:PSS SKU: M1216 CAS Number: 155090-83-8 / 7440-22-4 Options: • Quantity: 10 mL, 25 mL, 50 mL		

Inorganic Electronic Materials

Inorganic electronic materials include metals, metal oxides, carbon based nanomaterials and more. They possess unique electrical, magnetic, or optical properties, that make them suitable for a wide range of electronic, optoelectronic, and energy-related devices. These materials typically feature ionic, covalent, or metallic bonding and often exhibit high thermal stability, mechanical strength, and well-defined electronic behaviors.



Popular Products



Titanium Dioxide (TiO₂) Nanoparticles

High purity titanium dioxide nanoparticles for multiple applications

- Anatase, rutile, P25 phase
- Photocatalysts
- Wide bandgap



Silver Nanowire Dispersion

>99.9% purity, 5 mg/ml aqueous silver nanowire dispersion

- Easy handling
- Length: >20 μm
- Diameter: <30 nm



Aluminum Oxide (Al₂O₃) Powder

High purity (>99.5 %) gamma-phase aluminum oxide nanoparticles

- ~ 18 nm ave. particle size
- Mechanical & thermal stability

To help you choose the right material, full technical specifications are available on our website.



Solpro Business Park Windsor Street, Sheffield S10 2TE

Zinc Oxide (Y2O3) Nanopowder

Purity (by ICP-MS): 99.999%
MW: 77.547 g/mol

Average Particle Size = 20-30 nm

www.ossila.com

Solpro Business Park Windsor Street, Sheffield S10 2TE

TiO₂ Anatase Nanoparticles

Purity (by ICP-MS): 99.999%
MW: 79.866 g/mol

Average Particle Size = 20-30 nm

Ossila

M2398B1-25g

Zinc Oxide (ZnO)

Purity (by ICP-MS): 99.999%
MW: 81.376 g/mol

Average Particle Size = 20-30 nm

CAS# 1314-13-2

Metal Oxides

	Titanium Dioxide (TiO₂) Nanoparticles SKU: M2394 CAS Number: 1317-70-0 / 1317-80-2 / 13463-67-7 Options*: • Form: Powder, Solution in Water • Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml		Lithium Cobalt Oxide (LiCoO₂) Powder SKU: M2401 CAS Number: 12190-79-3 Options: • Particle size (d ₅₀): 13 µm • Quantity: 50 g, 100 g
	Titanium Dioxide Anatase Nanoparticles SKU: M2394 CAS Number: 1317-70-0 Options*: • Form: Powder, Solution in Water • Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml		Lithium Manganese Oxide (LMO) Powder SKU: M2402 CAS Number: 12057-17-9 Options: • Particle size (d ₅₀): 14.8 µm, 23 – 27 µm • Quantity: 100 g, 250 g, 500 g
	Rutile Titanium Dioxide (TiO₂) Nanoparticles SKU: M2394 CAS Number: 1317-80-2 Options*: • Quantity/volume: 50 g, 100 g, 250 g		Lithium Nickel Manganese Oxide (LNMO) Powder SKU: M2403 CAS Number: 12031-75-3 Options: • Quantity: 50 g, 100 g, 250 g
	Titanium Dioxide (TiO₂) P25 Nanoparticles SKU: M2394 CAS Number: 13463-67-7 Options*: • Quantity/volume: 50 g, 100 g, 250 g		Lithium Nickel Cobalt Aluminum Oxide (NCA) Powder SKU: M2404 CAS Number: 193214-24-3 Options: • Quantity: 10 g, 25 g, 50 g, 100 g
	Yttrium Oxide (Y₂O₃) Nanopowder SKU: M2397 CAS Number: 1314-36-9 Options: • Quantity: 50 g, 100 g, 250 g		Lithium Nickel Manganese Cobalt Oxide (NMC811) Powder SKU: M2405 CAS Number: 179802-95-0 Options: • Quantity: 50 g, 100 g
	Zinc Oxide Nanoparticles SKU: M2398 CAS Number: 1314-13-2 Options: • Quantity: 25 g, 50 g, 100 g		Bismuth Oxide (Bi₂O₃) Nanopowder SKU: M2407 CAS Number: 1304-76-3 Options: • Quantity: 50 g, 100 g, 250 g
	Antimony Oxide (Sb₂O₃) Powder SKU: M2399 CAS Number: 1309-64-4 Options: • Particle size: 20 – 30 nm, < 210 nm • Quantity: 50 g, 100 g, 250 g		Lithium Nickel Manganese Cobalt Oxide Powder (NCM523) SKU: M2409 CAS Number: 346417-97-8 / 182442-95-1 Options*: • Quantity: 50 g, 100 g, 250 g
	Aluminum Oxide (Al₂O₃) Nanopowder SKU: M2400 CAS Number: 1344-28-1 Options: • Quantity: 50 g, 100 g, 250 g, 500 g • Specification: Particle size 20 nm; Purity ≥ 99.99%, Particles size < 18 nm; Purity > 99.5%		Lithium Titanium Oxide Powder SKU: M2428 CAS Number: 12031-95-7 Options: • Quantity: 25 g, 50 g, 100 g, 250 g

* Some options have been omitted, please see our website to browse our full range

	Iron Oxide (Fe₃O₄) Nanopowder SKU: M2439 CAS Number: 1317-61-9 Options: - Average particle size: 30 – 50 nm - Quantity: 100 g, 250 g, 500 g		Antimony Tin Oxide Nanopowder SKU: M2487 CAS Number: 128221-48-7 Options: Quantity: 100 g, 250 g
	Iron Oxide (Fe₂O₃) Nanopowder SKU: M2440 CAS Number: 1309-37-1 Options: - Crystalline phase: α-Fe₂O₃, γ-Fe₂O₃ - Quantity: 10 g, 25 g, 50 g		Indium Oxide Powder SKU: M2488 CAS Number: 1312-43-2 Options: Quantity: 5 g, 10 g, 25 g
	Nickel (II) Oxide Powder SKU: M2441 CAS Number: 1313-99-1 Options: - Particle size: 8 - 18 nm - Quantity: 10 g, 25 g		Tantalum (V) Oxide Powder SKU: M2489 CAS Number: 1314-61-0 Options: Quantity: 10 g, 25 g, 50 g
	Barium Titanate Powder SKU: M2443 CAS Number: 12047-27-7 Options: - Average particle size: 90 nm - Quantity: 25 g, 50 g, 100 g, 250 g		Manganese (IV) Oxide Powder SKU: M2490 CAS Number: 1313-13-9 Options: Quantity: 100 g, 250 g
	Magnesium Oxide Powder SKU: M2483 CAS Number: 1309-48-4 Options: Quantity: 100 g, 250 g, 500 g		Lithium Nickel Manganese Cobalt Oxide Powder (NMC622) SKU: M2546 CAS Number: 346417-97-8 / 182442-95-1 Options: Quantity: 50 g, 100 g, 250 g
	Tin (IV) Oxide Powder SKU: M2484 CAS Number: 18282-10-5 Options: Quantity: 50 g, 100 g, 250 g		Lithium Nickel Manganese Cobalt Oxide Powder (NCM111) SKU: M2547 CAS Number: 346417-97-8 / 182442-95-1 Options: Quantity: 50 g, 100 g, 250 g
	Tungsten Oxide Powder SKU: M2485 CAS Number: 1314-35-8 Options: Quantity: 100 g, 250 g, 500 g		
	Aluminum Zinc Oxide Powder SKU: M2486 CAS Number: 952185-00-1 Options: Quantity: 50 g, 100 g		

Advanced Ceramic Materials

	Hexagonal Boron Nitride (h-BN) Crystals SKU: M2133 CAS Number: 10043-11-5 Options: • Package/substrate: Pack of 10 crystals, Pack of 5 crystals		Zinc Oxide Nanoparticles SKU: M2398 CAS Number: 1314-13-2 Options: • Quantity: 25 g, 50 g, 100 g
	Hexagonal Boron Nitride (h-BN) Powder SKU: M2133 CAS Number: 10043-11-5 Options: • Particle size: 1 – 2 µm, 100 – 200 nm, 30 – 40 µm, 40 – 50 µm, 790 nm, 800 nm • Quantity: 10 g, 25 g, 50 g, 100 g, 250 g, 500 g		Antimony Oxide (Sb2O3) Powder SKU: M2399 CAS Number: 1309-64-4 Options: • Particle size: 20 – 30 nm, < 210 nm • Quantity: 50 g, 100 g, 250 g
	Hexagonal Boron Nitride Nanopowders SKU: M2133 CAS Number: 10043-11-5 Options: • Particle size: 100 – 200 nm, 790 nm, 800 nm • Quantity: 10 g, 25 g, 50 g, 100 g, 250 g, 500 g		Aluminum Oxide (Al2O3) Nanopowder SKU: M2400 CAS Number: 1344-28-1 Options: • Quantity: 50 g, 100 g, 250 g, 500 g • Specification: Particle size 20 nm; Purity ≥ 99.99%, Particles size ~ 18 nm; Purity > 99.5%
	Titanium Dioxide (TiO2) Nanoparticles SKU: M2394 CAS Number: 1317-70-0 / 1317-80-2 / 13463-67-7 Options*: • Form: Powder, Solution in Water • Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml		Lithium Cobalt Oxide (LiCoO2) Powder SKU: M2401 CAS Number: 12190-79-3 Options: • Particle size (d50): 13 µm • Quantity: 50 g, 100 g
	Titanium Dioxide Anatase Nanoparticles SKU: M2394 CAS Number: 1317-70-0 Options*: • Form: Powder, Solution in Water • Quantity/volume: 50 g, 100 g, 250 g, 100 ml, 250 ml		Lithium Manganese Oxide (LMO) Powder SKU: M2402 CAS Number: 12057-17-9 Options: • Particle size (d50): 14.8 µm, 23 – 27 µm • Quantity: 100 g, 250 g, 500 g
	Rutile Titanium Dioxide (TiO2) Nanoparticles SKU: M2394 CAS Number: 1317-80-2 Options*: • Quantity/volume: 50 g, 100 g, 250 g		Lithium Nickel Manganese Oxide (LNMO) Powder SKU: M2403 CAS Number: 12031-75-3 Options: • Quantity: 50 g, 100 g, 250 g
	Titanium Dioxide (TiO2) P25 Nanoparticles SKU: M2394 CAS Number: 13463-67-7 Options*: • Quantity/volume: 50 g, 100 g, 250 g		Lithium Nickel Cobalt Aluminum Oxide (NCA) Powder SKU: M2404 CAS Number: 193214-24-3 Options: • Quantity: 10 g, 25 g, 50 g, 100 g
	Yttrium Oxide (Y2O3) Nanopowder SKU: M2397 CAS Number: 1314-36-9 Options: • Quantity: 50 g, 100 g, 250 g		Lithium Nickel Manganese Cobalt Oxide (NMC811) Powder SKU: M2405 CAS Number: 179802-95-0 Options: • Quantity: 50 g, 100 g

* Some options have been omitted, please see our website to browse our full range

	Bismuth Oxide (Bi₂O₃) Nanopowder SKU: M2407 CAS Number: 1304-76-3 Options: • Quantity: 50 g, 100 g, 250 g		Barium Titanate Powder SKU: M2443 CAS Number: 12047-27-7 Options: • Average particle size: 90 nm • Quantity: 25 g, 50 g, 100 g, 250 g
	Lithium Iron Phosphate (LiFePO₄) Powder SKU: M2408 CAS Number: 15365-14-7 Options*: • Quantity: 100 g, 250 g		Magnesium Oxide Powder SKU: M2483 CAS Number: 1309-48-4 Options: • Quantity: 100 g, 250 g, 500 g
	Lithium Nickel Manganese Cobalt Oxide Powder (NCM523) SKU: M2409 CAS Number: 346417-97-8 / 182442-95-1 Options*: • Quantity: 50 g, 100 g, 250 g		Tin (IV) Oxide Powder SKU: M2484 CAS Number: 18282-10-5 Options: • Quantity: 50 g, 100 g, 250 g
	Lithium Titanium Oxide Powder SKU: M2428 CAS Number: 12031-95-7 Options: • Quantity: 25 g, 50 g, 100 g, 250 g		Tungsten Oxide Powder SKU: M2485 CAS Number: 1314-35-8 Options: • Quantity: 100 g, 250 g, 500 g
	Hydroxyapatite Powder SKU: M2437 CAS Number: 12167-74-7 Options: • Average particle size: 3 µm, 10 - 20 µm • Quantity: 10 g, 25 g, 50 g		Aluminum Zinc Oxide Powder SKU: M2486 CAS Number: 952185-00-1 Options: • Quantity: 50 g, 100 g
	Iron Oxide (Fe₃O₄) Nanopowder SKU: M2439 CAS Number: 1317-61-9 Options: • Average particle size: 30 - 50 nm • Quantity: 100 g, 250 g, 500 g		Antimony Tin Oxide Nanopowder SKU: M2487 CAS Number: 128221-48-7 Options: • Quantity: 100 g, 250 g
	Iron Oxide (Fe₂O₃) Nanopowder SKU: M2440 CAS Number: 1309-37-1 Options: • Crystalline phase: α-Fe ₂ O ₃ , γ-Fe ₂ O ₃ • Quantity: 10 g, 25 g, 50 g		Indium Oxide Powder SKU: M2488 CAS Number: 1312-43-2 Options: • Quantity: 5 g, 10 g, 25 g
	Nickel (II) Oxide Powder SKU: M2441 CAS Number: 1313-99-1 Options: • Particle size: 8 - 18 nm • Quantity: 10 g, 25 g		Tantalum (V) Oxide Powder SKU: M2489 CAS Number: 1314-61-0 Options: • Quantity: 10 g, 25 g, 50 g

* Some options have been omitted, please see our website to browse our full range

	Manganese (IV) Oxide Powder SKU: M2490 CAS Number: 1313-13-9 Options: • Quantity: 100 g, 250 g
	Lithium Nickel Manganese Cobalt Oxide Powder (NMC622) SKU: M2546 CAS Number: 346417-97-8 / 182442-95-1 Options: • Quantity: 50 g, 100 g, 250 g
	Lithium Nickel Manganese Cobalt Oxide Powder (NCM111) SKU: M2547 CAS Number: 346417-97-8 / 182442-95-1 Options: • Quantity: 50 g, 100 g, 250 g
	Prussian Blue Analogue Powder SKU: M2551 CAS Number: 14038-43-8 Options: • Quantity: 10 g, 25 g
	Prussian White Powder SKU: M2552 CAS Number: 14038-43-8 Options: • Quantity: 10 g, 25 g, 50 g

Equipment

Our full range of equipment includes:

Thin Film Coating

- Automatic Film Applicator
- Bar Coater
- Doctor Blade Coater
- Dip Coater
- Slot-Die Coater
- Spin Coater
- Spin Coater Advanced
- Spin Coating Kit
- Syringe Pump
- Dual Syringe Pump

Sample Preparation

- Analytical Balance
- Microbalance
- UV Ozone Cleaner

Cleanroom Equipment

- Glove Box
- Laminar Flow Hood

Probe Station

- Manual Probe Station
- Micromanipulator

OEM Solutions

- OEM Source Measure Unit

Electrical Characterization

- Electrical Test Boards
- Four-Point Probe
- Four-Point Probe Plus
- LED Measurement System
- Manual Probe Station
- Micromanipulator
- Potentiostat
- Solar Cell I-V Test System
- Source Measure Unit

Solar Simulator

- LED Solar Simulator
- Solar Cell Testing Kit
- Indoor Light Simulator

Surface Characterization

- Contact Angle Goniometer
- Four-Point Probe
- Four-Point Probe Plus

Motion Control

- Linear Stages
- High-Precision Linear Stage

Electrochemistry Supplies

- Potentiostat
- Electrochemical Cells
- Photoelectrochemical Cells
- Reference Electrodes
- Counter Electrodes
- Working Electrodes

Spectroscopy Equipment

- Monochromator
- USB Spectrometer
- Optical Spectroscopy Kit
- LED Light Sources
- Optical Filters

Life Sciences

- Glove Box
- Laminar Flow Hood
- Analytical Balance
- Microbalance
- Micropipettes
- Syringe Pump
- Microscope Slides
- Laboratory Supplies



How to Order

All of our materials are available to purchase online through our website in specified quantities.

For large batches and alternative quantities, please contact us at **info@ossila.com** to enquire.

Request a Quote

We are happy to issue quotes, and you can order via offline channels by emailing a purchase order to **info@ossila.com**. Payment is available by credit card or bank transfer.

We offer net 30-day payment terms for offline orders as standard, but you can request a proforma invoice if you would like to pay upfront.

For orders placed through offline channels, in addition to credit card payments, we also accept payment via bank transfer in the following currencies:

- British Pound, **£ GBP**
- Canadian Dollar, **£ CAD**
- Euro, **€ EUR**
- Polish Złoty, **zł PLN**
- United States Dollar, **\$ USD**
- Romanian Leu, **Lei RON**
- Swedish Krona, **kr SEK**

Order Online

The quickest and easiest way to place your order is to add the items to your cart on www.ossila.com and then proceed to the online checkout. For online purchases, credit card payments are available in over 30 currencies.

To request a quote online, just add your products to your cart, and click "request quote". All quotes are valid for 30 days, will reflect the prices as advertised on our website in the currency of your choosing, and are covered by our price drop guarantee.