


Catholyte Reservoir

USER MANUAL

Manual version: 1.0.0
Product code: C2052D1



Contents

1. Packing List.....	3
2. Specifications.....	4
3. Related Products	5

1. Packing List

The standard items included with the **Catholyte Reservoir** are:

1. PTFE screw caps with PEEK electrode secure screw
2. Glassware body
3. Circular sampling window septa
4. PTFE rod 3 mm air inlet/outlet stopper x4
5. O-rings and spare PEEK electrode secure screw
6. PTFE tubing
7. Sampling needle




2. Specifications

The **Catholyte Reservoir** specifications are shown in **Table 2.1**.

Table 2.1. Ossila Catholyte Reservoir specifications.

Cell Window (W x H)	20 mm x 20 mm
Cell (W x H)	50 mm x 65 mm



3. Related Products



Potentiostat

Low Price Potentiostat for Cyclic Voltammetry and More

Product code: **T2006A1**



Standard Electrochemical Cell

High Quality Standard Electrochemical Cell

Product code: **C2012A1**



Ag/Ag⁺ Reference Electrodes

Non-aqueous Ag/Ag⁺ reference electrode suitable for cyclic voltammetry electrochemistry cells

Product code: **C2015A1**



L-Shaped Glassy Carbon Working Electrode

Perfect for use in standard electrochemistry cells for cyclic voltammetry and more

Product code: **C2013H1, C2013J1**