


Gas Stabilizer Cell

USER MANUAL

Manual version: 1.0.0
Product code: C2052A1



Contents

1. Packing List.....	3
2. Specifications.....	4
3. Related Products	5

1. Packing List

The standard items included with the **Gas Stabilizer Cell** are:

1. Glass body
2. Secure screw cap
3. PTFE lid with PEEK gas inlet/outlet secure screws
4. PTFE tubing
5. Gas locks
6. PEEK gas inlet/outlet secure screw and O-rings




2. Specifications

The **Gas Stabilizer Cell** specifications are shown in **Table 2.1**.

Table 2.1. Ossila Gas Stabilizer Cell specifications.

Cell Window (W x H)	40 mm x 58 mm
Cell Volume	100 ml



3. Related Products



Potentiostat

Low price Potentiostat for cyclic voltammetry and more

Product code: **T2006A1**



Standard Electrochemical Cell

High-quality standard electrochemical cell

Product code: **C2012A1**



Ag/Ag⁺ Reference Electrodes

Non-aqueous Ag/Ag⁺ reference electrode suitable for cyclic voltammetry electrochemistry cells

Product code: **C2015A1**



L-Shaped Glassy Carbon Working Electrode

Perfect for use in standard electrochemistry cells for cyclic voltammetry and more

Product code: **C2013H1, C2013J1**